
# ETHANOL CONTROL SET

2 x 5 ml

REF EH S010 CH

LOT -----

componente component composant componente	metodo method méthode método	valore value valeur valor	intervallo range intervalle intervalo	1SD 1s	unità units unités unidades
<b>LEV 1</b> Etanolo Ethanol Éthanol Etanol	metodo enzimatico enzymatic method méthode enzymatique método enzimático	-----	-----	-----	mg/dl mmol/l g/l
<b>LEV 2</b> Etanolo Ethanol Éthanol Etanol	metodo enzimatico enzymatic method méthode enzymatique método enzimático	-----	-----	-----	mg/dl mmol/l g/l

Italiano

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### Uso previsto

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## CONTROL SET

### Intended use

ETHANOL CONTROL SET contains liquid solutions for quality control of clinical chemistry methods for determination of ethanol. The solutions are used to monitor accuracy and/or precision both in automated analyzers and in manual use.

### Kit components

Ethanol Controls 2 x 5 ml

### Precautions and warnings

For in vitro diagnostic use only.

Exercise the normal precautions required for handling all laboratory reagents.

### Disposal

This product is made to be used in professional laboratories.

Please consult local regulations for a correct waste disposal.

P501: Dispose of contents according to national/international regulations.

### Assigned values and ranges

The declared value is the media of all values obtained. The corresponding control range is calculated as the target value  $\pm$  3 SD. Determinations were made using the testing reagents and analysis systems available at the time. Results must be within the defined ranges. Each laboratory should establish guidelines for corrective measures to be taken if values fall outside the range.

### Storage and stability

Store the product at 2-8°C.

The solutions are ready for use.

Ethanol is volatile! Immediately after use, close carefully the Control dropper.

After vial opening, the product is stable until expiration date if stored tightly closed at 2-8°C. Do not freeze.

### Test procedure

Dispense the required volume into a sample cup and analyze in the same way as for the samples. Because of ethanol volatility, the control solutions in open cups should be assayed within 1 h.

The control should be run daily in parallel with the patient samples and after every calibration. The control intervals should be adapted to each laboratory's individual requirements.

Additional materials required: automatic micropipette, distilled/deionized water, general laboratory equipment.

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