**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN  
OLYMPUS AU 400/480/600/640/680/2700 (Test code 891)**

TEST NAME: **RF**

SAMPLE: Volume **3 µl** Dilution **0 µl**

REAGENTS: R1 Volume **240 µl** Dilution **0 µl**  
R2 Volume **60 µl** Dilution **0 µl**

WAVELENGTH: Pri. **660** Sec. -

METHOD: **END**

REACTION SLOPE: **+**

MEASURING POINT 1: First **0** Last **27**

MEASURING POINT 2: First **0** Last **10**

REAGENT OD LIMIT: First L **-0.1** First H **0.5**  
Last L **-0.1** Last H **0.5**

DYNAMIC RANGE: L **7** H **200**

CORRELATION FACTOR: A **1** B **0**

UNIT: **IU/ml**

CALIBRATION TYPE: **5AB**

FORMULA: **POLYGONAL**

 Chema Diagnostica  
Via Campania 2/4  
60030 Monsano (AN) - ITALY - EU  
phone +39 0731 605064  
fax +39 0731 605672  
e-mail: mail@chema.com  
website: http://www.chema.com

ITALIANO rev. 19/03/2020

**FATTORE REUMATOIDE FL**

RF 6U140 2 x 56 + 2 x 14 ml

**USO**

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RF 6U140 2 x 56 + 2 x 14 ml

**INTENDED USE**

Reagent for quantitative in vitro determination of rheumatoid factor in biological fluids.

**PRINCIPLE OF THE METHOD**

Rheumatoid Factor (RF) selectively reacts with IgGs coated to latex, thus producing particle agglutination. The produced turbidity is proportional to the concentration of RF in the sample, and can be measured at the wavelength of 660 nm.

**KIT COMPONENTS**

**For in vitro diagnostic use only.**  
The components of the kit are stable until expiration date on the label.  
Keep away from direct light sources at 2-8°C.

**RF R1 2 x 56 ml (liquid) white cap**

Composition: Buffer pH 8.2, stabilizers and preservatives.

**RF R2 2 x 14 ml (liquid) red cap**

Composition: suspension of latex particles coated with human IgG, stabilizers and preservatives.

**REAGENT PREPARATION**

Use separate reagents ready to use.  
Kindly shake R2 vial before use.  
Stability: up to expiration date on labels at 2-8°C.  
Stability since first opening of vials: use preferably within 60 days at 2-8°C.

**PRECAUTIONS**

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.  
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

**SPECIMEN**

Serum, plasma.  
Keep specimens away from direct light sources.  
Samples are stable 7 days when stored at 2-8°C and 1 month at -20°C.

**EXPECTED VALUES**

Adulti < 30 IU/ml

Each laboratory should establish appropriate reference intervals related to its population.

**QUALITY CONTROL AND CALIBRATION**

It is suggested to perform an internal quality control.  
For this purpose the following human based control serum is available:  
**RHEUMATOID FACTOR CONTROL SET**

For calibration, use the following product:  
**RF STANDARD**

Please contact Customer Care for further information.

**TEST PERFORMANCE**

**Measuring range**  
Measure interval depends on the concentration of the highest standard used for calibration.  
If such a limit value is exceeded, it is suggested to dilute sample 1+4 with distilled water and to repeat the test, multiplying the result by 5.

**Hook effect**  
No Hook effect is observed with concentrations lower than 800 IU/ml.

**Sensitivity/limit of detection**  
The limit of detection is 7 IU/ml.

**Interferences**  
No interference was observed by the presence of:  
hemoglobin ≤ 400 mg/dl  
bilirubin ≤ 40 mg/dl  
lipids ≤ 500 mg/dl

**Precision**

| intra-assay (n=10) | mean (IU/ml) | SD (IU/ml) | CV%  |
|--------------------|--------------|------------|------|
| sample 1           | 24.8         | 0.60       | 2.44 |
| sample 2           | 79.6         | 0.55       | 0.69 |

| inter-assay (n=10) | mean (IU/ml) | SD (IU/ml) | CV%  |
|--------------------|--------------|------------|------|
| sample 1           | 25.2         | 2.39       | 9.48 |
| sample 2           | 79.4         | 2.37       | 2.98 |

**Methods comparison**  
A comparison between CHEMA and a commercially available product gave the following results:

$$\begin{aligned} \text{RF competitor} &= x \\ \text{RF FL CHEMA} &= y \\ n &= 48 \\ y &= 0.92x + 1.85 \text{ IU/ml} \quad r^2 = 0.97 \end{aligned}$$

**WASTE DISPOSAL**

This product is made to be used in professional laboratories.  
P501: Dispose of contents according to national/international regulations.

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