**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ  
HITACHI 911/912**

TEST:	<b>CA</b>	
APP. CODE:	<b>381</b>	
WAVELENGTH (Sec/Pri):	<b>700 - 660</b>	
ASSAY:	<b>1-POINT</b>	<b>TIME: 10 POINT: 15</b>
SAMPLE VOL:	NORMAL: <b>3</b> DECREASE: <b>2</b> INCREASE: <b>5</b>	
	R1 VOLUME: <b>300</b> R2 VOLUME: <b>0</b> R3 VOLUME: <b>0</b> R4 VOLUME: <b>0</b>	
ABS LIMIT:	<b>32000 - INC</b>	
PROZONE LIMIT:	<b>0 - UPPER</b>	
CALIB METHOD:	<b>LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)</b>	
SD LIMIT:	<b>0.250</b>	
DUPLICATE LIMIT:	<b>3%</b>	
ST. 1 CONC:	<b>0.0</b>	
EXPECTED VALUE:	<b>8.0 - 11.0</b>	
UNIT:	<b>mg/dl</b>	
INSTR. FACTOR (y=ax+b):	<b>a=1</b>	<b>b=0</b>

**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ  
OLYMPUS AU 400/480/600/640/680/2700 (Test code 862)**

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**INTENDED USE**  
Reagent for quantitative in vitro determination of calcium in biological fluids.

**PRINCIPLE OF THE METHOD**  
Arsenazo(III) combines with calcium at slight acidic pH to form a blue complex, the absorbance of which is measured at 660 nm. The reaction has high specificity and interference from magnesium is avoided, due to pH.  
For bichromatic analyzers, the reference wavelength must be set at 700 nm.

**KIT COMPONENTS**  
**For in vitro diagnostic use only.**  
The components of the kit are stable until expiration date on the label.  
Keep away from direct light sources.

**CA ASX R1 2H500: 10 x 50 ml (liquid) white cap  
6U448: 8 x 56 ml (liquid) white cap**

Composizione: arsenazo(III) 0.2 mM, Good's buffer 50 mM pH 6.8, stabilizers.

Store all components at 15-25°C.

**REAGENT PREPARATION**  
Use reagent ready to use.  
Stability: up to expiration date on labels at 15-25°C.  
Stability since first opening of vials: preferably within 60 days at 15-25°C.

**PRECAUTIONS**  
Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.  
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

**SPECIMEN**  
Serum (preferred), plasma heparinate. Do not use citrate, oxalate and EDTA as anticoagulant.  
Total calcium is stable 7 days at 2-8°C and for several months when frozen at -20°C.  
Urine specimens should be collected in 20 to 30 ml of HCl 6M per 24/h specimen (1-2 ml for random urine) in order to prevent calcium salt precipitation.  
Dilute sample urine 1:2 with redistilled water and multiply results by two.

**EXPECTED VALUES**

serum/plasma:	8.6 - 10.3 mg/dl (2.15 - 2.57 mmol/l)
urine (men):	up to 300 mg/24h (7.49 mmol/24h)
urine (women):	up to 250 mg/24h (6.24 mmol/24h)

Each laboratory should establish appropriate reference intervals related to its population.

**QUALITY CONTROL AND CALIBRATION**

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:  
**QUANTINORM CHEMA**  
with normal or close to normal control values  
**QUANTIPATH CHEMA**  
with pathological control values.  
If required, a multiparametric, human based calibrator is available:  
**AUTOCAL H**

Please contact Customer Care for further information.

**TEST PERFORMANCE**

**Linearity**  
the method is linear up to 20 mg/dl.  
If the limit value is exceeded, it is suggested to dilute sample 1+9 with distilled water and to repeat the test, multiplying the result by 10.

**Sensitivity/limit of detection (LOD)**  
the limit of detection is 0.2 mg/dl.

**Interferences**  
no interference was observed by the presence of:  
hemoglobin ≤ 450 mg/dl  
bilirubin ≤ 50 mg/dl  
Lipids interferences are possible performing readings at single wavelength of 660 nm. To avoid interferences, perform a bichromatic reading at 660 / 700 nm.

**Precision**

intra-assay (n=10)	mean (mg/dl)	SD (mg/dl)	CV%
sample 1	8.89	0.10	1.10
sample 2	13.74	0.16	1.20

inter-assay (n=20)	mean (mg/dl)	SD (mg/dl)	CV%
sample 1	9.22	0.19	2.10
sample 2	14.04	0.23	1.70

**Methods comparison**  
a comparison between Chema and a commercially available product gave the following results:

$$\text{Calcium ASX Chema} = x$$

$$\text{Calcium competitor} = y$$

$$n = 97$$

$$y = 0.98x + 0.17 \text{ mg/dl} \quad r^2 = 0.94$$

**WASTE DISPOSAL**

This product is made to be used in professional laboratories.  
P501: Dispose of contents according to national/international regulations.