**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ  
HITACHI 911/912**

TEST:	<b>EPS</b>	
APP. CODE:	<b>333</b>	
WAVELENGTH (Sec/Pri):	<b>700 - 415</b>	
ASSAY:	<b>RATE-A</b>	<i>TIME: 10 POINT: 14 - 25 DILUENT: water</i>
SAMPLE VOL:	NORMAL: <b>10</b> DECREASE: <b>8</b> INCREASE: <b>12</b>	
	R1 VOLUME: <b>240</b> R2 VOLUME: <b>60</b> R3 VOLUME: <b>0</b> R4 VOLUME: <b>0</b>	<i>DILUENT: 5 DILUENT: 5</i>
ABS LIMIT:	<b>32000 - INC</b>	
PROZONE LIMIT:	<b>0 - UPPER</b>	
CALIB METHOD:	<b>LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)</b>	
SD LIMIT:	<b>0.250</b>	
DUPLICATE LIMIT:	<b>3%</b>	
ST. 1 CONC:	<b>0.00</b>	
EXPECTED VALUE:	<b>0 - 100</b>	
UNIT:	<b>U/l</b>	
INSTR. FACTOR (y=ax+b):	a=1 b=0	

**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ  
OLYMPUS AU 400/480/600/640/680/2700 (Test code 858)**

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**INTENDED USE**  
Reagent for quantitative in vitro determination of amylase in biological fluids.

**PRINCIPLE OF THE METHOD**  
The enzyme α-amylase (EC 3.2.1.1, 1.4 α-D-glucose glucanohydrolase) hydrolyzes the EPS to release several different fragments. The fragments so formed are completely hydrolyzed to 4-nitrophenol and glucose by α-glucosidase. The 4-nitrophenol formed is detected spectrophotometrically at 405 nm to give a measurement of α-amylase activity in the sample.  
The present method has been made according to IFCC.

**KIT COMPONENTS**

**For in vitro diagnostic use only.**  
The components of the kit are stable until expiration date on the label.  
Keep away from direct light sources.

**AMY EPS R1 2H100: 4 x 20 ml (liquid) white cap  
6U280: 4 x 56 ml (liquid) white cap**

**AMY EPS R2 EA 2H100: 2 x 10 ml (liquid) red cap  
EA 6U280: 4 x 14 ml (liquid) red cap**

Composition in the test: Hepes buffer 50 mM pH 7.10, NaCl 70 mM, calcium acetate 1.0 mM, EPS-G7 5.0 mM, α-glucosidase 6 kU/l.

Store all components at 2-8°C.

**REAGENT PREPARATION**

Use reagent ready to use.  
Stability: up to expiration date on labels at 2-8°C.  
Stability since first opening of vials: preferably within 60 days at 2-8°C.

**PRECAUTIONS**

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.  
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

**SPECIMEN**

Serum, plasma (heparinate only), Urine.  
Amylase is stable in serum and plasma sample up to 2 months at 2-8°C.

**EXPECTED VALUES**

Serum - plasma: 28 - 100 U/l (0.47 - 1.67 µkat/l)  
Random urine: ≤ 460 U/l (≤ 7.68 µkat/l)

Each laboratory should establish appropriate reference intervals related to its population.

**QUALITY CONTROL AND CALIBRATION**

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

**QUANTINORM CHEMA**  
with normal or close to normal control values  
**QUANTIPATH CHEMA**  
with pathological control values.  
If required, a multiparametric, human based calibrator is available:  
**AUTOCAL H**

Please contact Customer Care for further information.

**TEST PERFORMANCE**

**Linearity**  
the method is linear up to 1500 U/l.  
If a ΔA/min of 0.500 is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

**Sensitivity/limit of detection (LOD)**  
the limit of detection is 6 U/l.

**Interferences**  
no interference was observed by the presence of:  
hemoglobin ≤ 200 mg/dl  
bilirubin ≤ 48 mg/dl  
lipids interfere in low values

**Precision**

intra-assay (n=10)	mean (U/l)	SD (U/l)	CV%
sample 1	77.90	0.74	0.90
sample 2	194.80	1.99	1.00

**inter-assay (n=20)**

mean (U/l)	SD (U/l)	CV%
sample 1 75.77	1.90	2.50
sample 2 194.15	2.39	1.20

**Methods comparison**  
a comparison between Chema Amylase FL and a commercially available product gave the following results:

$$y = 1.067x + 5.21 \quad r^2 = 0.99$$

**WASTE DISPOSAL**

This product is made to be used in professional laboratories.  
P501: Dispose of contents according to national/international regulations.

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