APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912	
TEST:	HDL
APP. CODE:	035
WAVELENGTH (Sec/Pri):	700 - 600
ASSAY:	2-POINT END <i>TIME: 10</i> <i>POINT: 3 - 31</i>
SAMPLE VOL:	NORMAL: 3 DECREASE: 2 INCREASE: 5
	R1 VOLUME: 270 R2 VOLUME: 90 R3 VOLUME: 0 R4 VOLUME: 0
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)
SD LIMIT:	0.250
DUPLICATE LIMIT:	3%
ST. 1 CONC:	0.0
EXPECTED VALUE:	30.0 - 90.0
UNIT:	mg/dl
INSTR. FACTOR (y=ax+b):	a=1 b= 0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 866)	
TEST NAME:	HDL
SAMPLE:	Volume 3 µl Dilution 0 µl
REAGENTS:	R1 Volume 270 µl Dilution 0 µl R2 Volume 90 µl Dilution 0 µl
WAVELENGTH:	Pri. 600 Sec. 700
METHOD:	END
REACTION SLOPE:	+
MEASURING POINT 1:	First 0 Last 27
MEASURING POINT 2:	First 0 Last 10
REAGENT OD LIMIT:	First L -0.1 First H 0.5 Last L -0.1 Last H 0.5
DYNAMIC RANGE:	L 1.0 H 220
CORRELATION FACTOR:	A 1 B 0
UNIT:	mg/dl
CALIBRATION TYPE:	AB
FORMULA:	Y = AX + B

 Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

ITALIANO rev. 09/05/2019

HDL-direct FL	
HD 2H160	6 x 20 + 2 x 20 ml
HD 6U240	4 x 45 + 4 x 15 ml

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INTENDED USE
Reagent for quantitative in vitro determination of HDL-cholesterol in biological fluids.

PRINCIPLE OF THE METHOD
Anti human β-lipoprotein antibody in reagent R1 binds to lipoproteins (LDL, VLDL, and chylomicrons) other than HDL. The antigen-antibody complexes formed block enzyme reactions when reagent R2 is added. Cholesterol esterase and cholesterol oxidase in reagent R2 react only with HDL-C. Hydrogen peroxide produced by the enzyme reactions with HDL-C yields a blue color complex upon oxidative condensation of F-DAOS [N-ethyl-N-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxy-4-fluoroaniline, sodium salt] and 4-aminoantipyrine in the presence of peroxidase. By measuring the absorbance of the blue color complex produced, at the optimum wavelength of 593 nm, the HDL-C concentration in the sample can be calculated when compared with the absorbance of the HDL-C Calibrator.

KIT COMPONENTS
For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label.
Keep away from direct light sources.

HDL-C R1 **2H160: 6 x 20 ml (liquid) white cap**
6U240: 4 x 45 ml (liquid) white cap

Composizione: Good's buffer pH 7.0 30 mmol/l, 4-aminoantipirina 0.9 mmol/l, POD 2400 U/l, ascorbate ossidasi 2700 U/l, and anti human-lipoprotein antibody, blend of 5-chloro-2-methyl-2-H-isothiazol-3-one [EC No 247-500-7] and 2-methyl-2-H-isothiazol-3-one [EC No 220-239-6] (3:1) in concentration 0.0015-0.06%.

HDL-C R2 **2H160: 2 x 20 ml (liquid) red cap**
6U240: 4 x 15 ml (liquid) red cap

Composizione: Good's buffer pH 7.0 30 mmol/l, colesterolo esterasi 4000 U/l, colesterolo ossidasi 20000 U/l, and F-DAOS 0.8 mmol/l.

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PRECAUZIONI
HDL-C R1: Warning. May cause an allergic skin reaction (H317). Avoid breathing vapours (P261). IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water (P303+P361+P353). If skin irritation or rash occurs: Get medical advice (P333+P313).

HDL-C R2: It is not classified as hazardous.

N-acetylcysteine (NAC), metamizole and acetaminophen may cause interference in the Trinder reaction.^(1,2) To avoid interference, the blood withdrawal should be performed before drug administration.

Store all components at 2-8°C.

REAGENT PREPARATION

Use separate reagent ready to use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-.
Caution: keep well refrigerated.

PRECAUTIONS

HDL-C R1: Warning. May cause an allergic skin reaction (H317). Avoid breathing vapours (P261). IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water (P303+P361+P353). If skin irritation or rash occurs: Get medical advice (P333+P313).

HDL-C R2: It is not classified as hazardous.

N-acetylcysteine (NAC), metamizole and acetaminophen may cause interference in the Trinder reaction.^(1,2) To avoid interference, the blood withdrawal should be performed before drug administration.

SPECIMEN
Use serum as a specimen. It is recommended to measure HDL-C immediately after collection. Ascorbic acid, bilirubin, and hemoglobin do not have a significant effect on the measurement.

EXPECTED VALUES
Adult male: 35.3 - 79.5 mg/dl
Adult female: 42.0 - 88.0 mg/dl

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION
It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:
QUANTINORM CHEMA
with normal or close to normal control values
QUANTIPATH CHEMA
with pathological control values.
If required, a multiparametric, human based calibrator is available:
AUTOCAL H

Please contact Customer Care for further information.

TEST PERFORMANCE

Linearity
the method is linear up to 220 mg/dl.
If the limit is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)
the limit of detection is 1 mg/dl.

Interferences
no interference was observed by the presence of:
hemoglobin ≤ 500 mg/dl
bilirubin (free) ≤ 50 mg/dl
bilirubin (conjugated) ≤ 40 mg/dl
ascorbic acid ≤ 50 mg/dl

Precision
intra assay (n=10) mean (mg/dl) SD (mg/dl) CV%
sample 1 32.1 0.18 0.55
sample 2 88.9 0.61 0.68

Comparison
a comparison between HDL-direct FL and other methods has shown the following results:

HDL-C Chema = x
HDL-C competitor = y
n = 50

y = 0.96x + 2.5 mg/dl r² = 0.998

WASTE DISPOSAL
This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.