APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912		
TEST:	GL	
APP. CODE:	249	
WAVELENGTH (Sec/Pri):	700 - 505	
ASSAY:	1-POINT	TIME: 10 POINT: 31
SAMPLE VOL:	NORMAL: 3 DECREASE: 2 INCREASE: 5	
	R1 VOLUME: 300 R2 VOLUME: 0 R3 VOLUME: 0 R4 VOLUME: 0	
ABS LIMIT:	32000 - INC	
PROZONE LIMIT:	0 - UPPER	
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)	
SD LIMIT:	0.250	
DUPLICATE LIMIT:	3%	
ST. 1 CONC:	0.0	
EXPECTED VALUE:	70 - 110	
UNIT:	mg/dl	
INSTR. FACTOR (y=ax+b):	a=1	b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 876)		
TEST NAME:	GL	
SAMPLE:	Volume 3 µl	Dilution 0 µl
REAGENTS:	R1 Volume 300 µl R2 Volume 0 µl	Dilution 0 µl Dilution 0 µl
WAVELENGTH:	Pri. 520	Sec. 700
METHOD:	END	
REACTION SLOPE:	+	
MEASURING POINT 1:	First 0	Last 27
MEASURING POINT 2:	First	Last
REAGENT OD LIMIT:	First L -0.1 Last L -0.1	First H 1.5 Last H 1.5
DYNAMIC RANGE:	L 1	H 500
CORRELATION FACTOR:	A 1	B 0
UNIT:	mg/dl	
CALIBRATION TYPE:	AB	
FORMULA:	Y = AX + B	

 Chema Diagnostica  
Via Campania 2/4  
60030 Monsano (AN) - ITALY - EU  
phone +39 0731 605064  
fax +39 0731 605672  
e-mail: mail@chema.com  
website: http://www.chema.com

ITALIANO rev. 02/09/2019

GLUCOSIO FL	
GL 2H500	10 x 50 ml
GL 6U448	8 x 56 ml

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**INTENDED USE**  
Reagent for quantitative in vitro determination of glucose in biological fluids.

**PRINCIPLE OF THE METHOD**  
The enzyme glucose oxidase catalyzes the oxidation of glucose to gluconic acid and H<sub>2</sub>O<sub>2</sub>. The H<sub>2</sub>O<sub>2</sub> reacts with phenol and 4-aminoantipyrine in the presence of peroxidase to form a quinoneimine dye. The intensity of color formed is proportional to the glucose concentration and can be measured photometrically between 480 and 520 nm.

**KIT COMPONENTS**  
**For in vitro diagnostic use only.**  
The components of the kit are stable until expiration date on the label.  
Keep away from direct light sources.

**GLU R1 2H500: 10 x 50 ml (liquido) capsula bianca**  
**6U448: 8 x 56 ml (liquido) capsula bianca**

Composition: phosphate buffer pH 6.50 220 mM, GOD ≥ 15000 U/l, POD ≥ 500 U/l, 4-AAP 1 mM, phenol 10 mM, surfactant.

Store all components at 2-8°C.

**REAGENT PREPARATION**  
Use separate reagent ready to use.  
Stability: up to expiration date on labels at 2-8°C.  
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-.  
Caution: keep well refrigerated.

**PRECAUTIONS**  
**GLU R1: Warning.** Causes serious eye irritation (H319). Causes skin irritation (H315). Wear protective gloves. Eye protection (P280). IF ON SKIN: Wash with plenty of water (P302+P352). IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). If eye irritation persists: get medical advice (P337+P313).

N-acetylcysteine (NAC), metamilzolo and acetaminophen may cause interference in the Trinder reaction.<sup>(1,2)</sup>  
To avoid interference, the blood withdrawal should be performed before drug administration.

**SPECIMEN**  
Serum, plasma, urine, CSF (cerebrospinal fluid).  
Immediately separate the samples from the cells.

**EXPECTED VALUES**  
Plasma/serum (fasting patient)  
adults: 70 - 105 mg/dl  
children: 70 - 105 mg/dl  
premature neonates: 25 - 80 mg/dl  
term neonates: 30 - 90 mg/dl  
CSF: 40 - 75 mg/dl  
(60% of plasma value)

Urine (fasting patient)  
random urine: < 30 mg/dl  
24h urine: < 500 mg/24h

Each laboratory should establish appropriate reference intervals related to its population.

**QUALITY CONTROL AND CALIBRATION**

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

**QUANTINORM CHEMA**  
with normal or close to normal control values

**QUANTIPATH CHEMA**  
with pathological control values.

If required, a multiparametric, human based calibrator is available:

**AUTOCAL H**  
Please contact Customer Care for further information.

**TEST PERFORMANCE**

**Linearity**  
the method is linear up to 500 mg/dl.  
If the limit value is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

**Sensitivity/limit of detection (LOD)**  
the limit of detection is 1 mg/dl.

**Interferences**  
no interference was observed by the presence of:  
hemoglobin ≤ 400 mg/dl  
bilirubin ≤ 20 mg/dl  
lipids ≤ 400 mg/dl

**Precision**  
intra-assay (n=10) mean (mg/dl) SD (mg/dl) CV%  
sample 1 91.8 0.65 0.70  
sample 2 241.1 3.34 1.39

inter-assay (n=20) mean (mg/dl) SD (mg/dl) CV%  
sample 1 92.2 2.37 2.60  
sample 2 240.6 8.11 3.40

**Methods comparison**  
a comparison between Chema and a commercially available product gave the following results:

Glucose FL Chema = x  
Glucose competitor = y  
n = 111  
y = 0.960x + 0.39 mg/dl r<sup>2</sup> = 0.984

**WASTE DISPOSAL**

This product is made to be used in professional laboratories.  
P501: Dispose of contents according to national/international regulations.