APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN HITACHI 911/912	
TEST:	HbA1c
APP. CODE:	361
WAVELENGTH (Sec/Pri):	800 - 660
ASSAY:	2 POINT END <i>TIME: 10</i> <i>POINT: 19 - 31</i> <i>DILUENT: water</i>
SAMPLE VOL:	NORMAL: 9 DECREASE: 6 INCREASE: 11
	R1 VOLUME: 225 R2 VOLUME: 0 R3 VOLUME: 75 R4 VOLUME: 0
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	Spline (POINT: 5 - SPAN: 5 - WEIGHT: 0)
SD LIMIT:	0.250
DUPLICATE LIMIT:	15%
ST. 1-2-3-4-5 CONC:	0.00 - user def.
EXPECTED VALUE:	0 - 6
UNIT:	%
INSTR. FACTOR (y=ax+b):	a=1 b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN OLYMPUS AU 400/480/600/640/680/2700 (Test code 890)	
TEST NAME:	HbA1c
SAMPLE:	Volume 6 µl Dilution 0 µl
REAGENTS:	R1 Volume 150 µl Dilution 0 µl R2 Volume 50 µl Dilution 0 µl
WAVELENGTH:	Pri. 660 Sec. 800
METHOD:	FIXED1
REACTION SLOPE:	+
MEASURING POINT 1:	First 14 Last 27
MEASURING POINT 2:	First 0 Last 0
REAGENT OD LIMIT:	First L -0.1 First H 0.5 Last L -0.1 Last H 0.5
DYNAMIC RANGE:	L 0.01 H 30.0
CORRELATION FACTOR:	A 1 B 0
UNIT:	%
CALIBRATION TYPE:	5AB
FORMULA:	SPLINE

 Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

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INTENDED USE
Reagent for quantitative in vitro determination of glycosylated hemoglobin (HbA1c) in biological fluids.

PRINCIPLE OF THE METHOD
Glycosylated hemoglobin (HbA1c) is an important biochemical marker for assessing long-term glycemic control in patients with diabetes. HbA1c is being used not only to guide diabetes treatment but also to assess quality of care, and to predict risk for development and progression of diabetes complications.

KIT COMPONENTS
For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label.
Keep away from direct light sources at 2-8°C.

HC R1 2H080: 3 x 20 ml (liquid) white cap
6U060: 1 x 45 ml (liquid) white cap

Composition: Latex suspension, stabilizers and preservatives.

HC R2 2H080: 1 x 20 ml (liquid) red cap
6U060: 1 x 15 ml (liquid) red cap

Composition: Anti-HbA1c antibody, stabilizers and preservatives.

REAGENT PREPARATION
Use separate reagents ready to use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: use preferably within 60 days at 2-8°C.

PRECAUTIONS
Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

SPECIMEN
Whole blood.
Keep specimens away from direct light sources.
Specimens are stable 2 weeks when stored at 2-8°C. Do not freeze.

Preparation of test sample:
10 µl of whole blood have to be diluted with 500 µl of deionized water, in order to hemolyze the sample.

The above procedure does not apply to reconstituted calibrators and controls.

EXPECTED VALUES
Normal values (NGSP) 4.6-6.2%

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose the following human based control set is available:

HbA1c CONTROL SET
with normal and abnormal control values.

For calibration, use the following product:
HbA1c CALIBRATOR

Please contact Customer Care for further information.

TEST PERFORMANCE
Measuring range
Measure interval depends on the concentration of the highest standard used for calibration.
If such a limit value is exceeded, sample should not be diluted and tested. The result should be reported as "high".

Interferences
No interference was observed by the presence of:
free bilirubin ≤ 20 mg/dl
conjugated bilirubin ≤ 21 mg/dl
ascorbic acid ≤ 50 mg/dl
rheumatoid factor ≤ 1100 IU/ml

Precision
intra-assay (n=10) mean (%) SD (%) CV%
sample 1 5.0 0.01 0.25
sample 2 12.8 0.02 0.17

inter-assay (n=5) mean (%) SD (%) CV%
sample 1 5.0 0.02 0.32
sample 2 12.9 0.06 0.45

Methods comparison
A comparison between CHEMA and HPLC method gave the following results:

HPLC method = x
HbA1c CHEMA = y
n = 142

y = 0.99x + 0.02% r² = 0.989

WASTE DISPOSAL
This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.

HbA1c	
<p>HC 2H080 3 x 20 + 1 x 20 ml</p> <p>HC 6U060 1 x 45 + 1 x 15 ml</p>	

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