



Cobas Mira instructions

rev 401.0.2 - 2010-07-20

Chema
D I A G N O S T I C A

Chema Diagnostica
Via Campania 2/4
I-60030 Monsano | An Italy
T +39 0731 605064 | F +39 0731 605672
www.chema.com | mail@chema.com

Cobas Mira instructions

rev 401.0.2 - 2010-07-20

Warning!

These programs should be used only as guideline. To check the reagent correct setting on the instrument, use only good quality control sera and work accordingly to with good laboratory practice. Instrumental factor must be checked. Chema Diagnostica will not be responsible for instrumental bad programming.

Analyzer: Cobas Mira

Application: **ACID PHOSPHATASE (TOTAL)** - Codes AC 0120 TC
 Preparation: POWDER SINGLE REAGENT. FOLLOW PREPARATION PROGRAM INDICATED IN INSERT SHEET.
 Storage: REFRIGERATE AT 2-8°C
 Stability: 10 DAYS

PROGRAM**GENERAL**

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: FACTOR
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 405 nm
 DECIMAL POSITION: 2
 UNIT: U/L

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 25.0 µl DIL: 25.0 µl

REAGENT CYCLE: 1
 VOL: 225 µl

CALCULATION

SAMPLE LIMIT: NO

REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 U/L
 HIGH: 60 U/L
 NORM. RANGE LOW: NO
 HIGH: NO

NUMBER OF STEPS: 1

CALC. STEP A: KINETIC
 READINGS FIRST: 10 LAST: 20
 REACTION LIMIT: 0.9000 ΔA
 POINT: T1

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: NO
 HIGH: NO
 BLANK RANGE LOW: NO
 HIGH: NO

FACTOR: 1421

CONTROL

CS1 POS: UTIL. DEF.
 CS2 POS: UTIL. DEF.
 CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **ACID PHOSPHATASE (NON PROSTATIC)** - Codes AC 0120 TC
 Preparation: POWDER SINLE REAGENT. FOLLOW PREPARATION PROGRAM INDICATED IN INSERT SHEET. ADD THE RIGHT QUANTITY OF L-TARTRATE.
 Storage: REFRIGERATE AT 2-8°C
 Stability: 10 DAYS

PROGRAM

GENERAL

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: FACTOR
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 405 nm
 DECIMAL POSITION: 2
 UNIT: U/L

CALIBRATION

CALIB. INTERVAL: ON REQUEST
 REAGENT BLANK
 REAG. RANGE LOW: NO
 HIGH: NO
 BLANK RANGE LOW: NO
 HIGH: NO
 FACTOR: 1421

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO
 SAMPLE CYCLE: 1
 VOL: 25.0 µl DIL: 25.0 µl
 REAGENT CYCLE: 1
 VOL: 225 µl

CONTROL

CS1 POS: UTIL. DEF.
 CS2 POS: UTIL. DEF.
 CS3 POS: UTIL. DEF.

CALCULATION

SAMPLE LIMIT: NO
 REAC. DIRECTION: INCREASE
 CHECK: ON
 CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000
 TEST RANGE LOW: 0 U/L
 HIGH: 60 U/L
 NORM. RANGE LOW: NO
 HIGH: NO
 NUMBER OF STEPS: 1
 CALC. STEP A: KINETIC
 READINGS FIRST: 10 LAST: 20
 REACTION LIMIT: 0.9000 ΔA
 POINT: T1

Analyzer: Cobas Mira

Application: **ALBUMIN** - Codes BC BC 0100 / 0500 / 1000 / 1500 CH
 Preparation: LIQUID READY TO USE SINGLE REAGENT
 Storage: ROOM TEMPERATURE (2-30°C)
 Stability: UP TO DECLARED EXPIRATION ON REAGENT LABEL

PROGRAM

GENERAL

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 600 nm
 DECIMAL POSITION: 1
 UNIT: g/dl

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 3.0 µl DIL: 20.0 µl

REAGENT CYCLE: 1
 VOL: 250 µl

CALCULATION

SAMPLE LIMIT: 0.1000 ΔA
 POINT: T1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0.0 g/dl
 HIGH: 6.0 g/dl
 NORM. RANGE LOW: 4.2 g/dl
 HIGH: 5.5 g/dl

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
 READINGS FIRST: CB LAST: 6

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: -0.0500 A
 HIGH: 0.3200 A
 BLANK RANGE LOW: -0.0500 A
 HIGH: 0.1000 A

CALIBRATOR CUP POS: 1
 CAL-1: UTIL. DEF

REPLICATE: TRIPL.
 DEVIATION: 6%

CONTROL

CS1 POS: UTIL. DEF.

CS2 POS: UTIL. DEF.

CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **ALKALINE PHOSPHATASE DGKC FL** - Codes AL F080 / F245 / F400 / F600 CH
 Preparation: AS INDICATED IN INSERT SHEET
 Storage: REFRIGERATE AT 2-8°C
 Stability: 30 DAYS

PROGRAM**GENERAL**

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: FACTOR
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 405 nm
 DECIMAL POSITION: 0
 UNIT: U/L

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 5.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 5.0 µl DIL: 10.0 µl

REAGENT CYCLE: 1
 VOL: 240 µl

CALCULATION

SAMPLE LIMIT: 0.2000 ΔA
 POINT: T1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 U/L
 HIGH: 2800 U/L
 NORM. RANGE LOW: 0 U/L
 HIGH: 270 U/L

NUMBER OF STEPS: 1

CALC. STEP A: KINETIC
 READINGS FIRST: 4 LAST: 8
 REACTION LIMIT: 0.4300 ΔA
 POINT: T1

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: -0.0500 A
 HIGH: 1.0000 A
 BLANK RANGE LOW: -0.0100 A
 HIGH: 0.0100 A

FACTOR: 4595

CONTROL

CS1 POS: UTIL. DEF.
 CS2 POS: UTIL. DEF.
 CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **ALKALINE PHOSPHATASE IFCC FL** - Codes AF F080 / F245 / F400 / F600 CH
 Preparation: AS INDICATED IN INSERT SHEET
 Storage: REFRIGERATE AT 2-8°C
 Stability: 30 DAYS

PROGRAM**GENERAL**

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: FACTOR
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 405 nm
 DECIMAL POSITION: 0
 UNIT: U/L

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 5.0 µl DIL: 10.0 µl

REAGENT CYCLE: 1
 VOL: 240 µl

CALCULATION

SAMPLE LIMIT: 0.2000 ΔA
 POINT: T1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 U/L
 HIGH: 3000 U/L
 NORM. RANGE LOW: 35 U/L
 HIGH: 104 U/L

NUMBER OF STEPS: 1

CALC. STEP A: KINETIC
 READINGS FIRST: 3 LAST: 9
 REACTION LIMIT: 0.4300 ΔA
 POINT: T1

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: -0.0500 A
 HIGH: 1.0000 A
 BLANK RANGE LOW: -0.0100 A
 HIGH: 0.0100 A

FACTOR: 4595

CONTROL

CS1 POS: UTIL. DEF.
 CS2 POS: UTIL. DEF.
 CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **AMYLASE FL** - Codes AM F060 / F120 / F245 CH
 Preparation: LIQUID READY TO USE SINGLE REAGENT
 Storage: REFRIGERATE AT 2-8°C
 Stability: UP TO DECLARED EXPIRATION ON REAGENT LABEL

PROGRAM**GENERAL**

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: FACTOR
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 405 nm
 DECIMAL POSITION: 0
 UNIT: U/L

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 7.0 µl DIL: 20 µl

REAGENT CYCLE: 1
 VOL: 260 µl

CALCULATION

SAMPLE LIMIT: 0.1500 ΔA
 POINT: T1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 U/L
 HIGH: 2000 U/L
 NORM. RANGE LOW: 0 U/L
 HIGH: 96 U/L

NUMBER OF STEPS: 1

CALC. STEP A: KINETIC
 READINGS FIRST: 3 LAST: 7
 REACTION LIMIT: 0.4000 ΔA
 POINT: T1

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: -0.0500 A
 HIGH: 0.2000 A
 BLANK RANGE LOW: -0.0100 A
 HIGH: 0.0100 A

FACTOR: 5296

CONTROL

CS1 POS: UTIL. DEF.
 CS2 POS: UTIL. DEF.
 CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **AMYLASE EPS FL** - Codes EA F080 / F245 CH
 Preparation: AS INDICATED IN INSERT SHEET
 Storage: REFRIGERATE AT 2-8°C
 Stability: 30 DAYS

PROGRAM

GENERAL

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: FACTOR
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 405 nm
 DECIMAL POSITION: 0
 UNIT: U/L

CALIBRATION

CALIB. INTERVAL: ON REQUEST
 REAGENT BLANK
 REAG. RANGE LOW: -0.0700 A
 HIGH: 0.8000 A
 BLANK RANGE LOW: -0.0100 A
 HIGH: 0.0200 A
 FACTOR: 5800

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 5.00
 CONC. FACTOR: NO
 SAMPLE CYCLE: 1
 VOL: 8.0 µl DIL: 20.0 µl
 REAGENT CYCLE: 1
 VOL: 220 µl

CONTROL

CS1 POS: UTIL. DEF.
 CS2 POS: UTIL. DEF.
 CS3 POS: UTIL. DEF.

CALCULATION

SAMPLE LIMIT: NO
 REAC. DIRECTION: INCREASE
 CHECK: ON
 CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000
 TEST RANGE LOW: 0 U/L
 HIGH: 1500 U/L
 NORM. RANGE LOW: 28 U/L
 HIGH: 100 U/L

NUMBER OF STEPS: 1

CALC. STEP A: KINETIC
 READINGS FIRST: 8 LAST: 16
 REACTION LIMIT: NO

Analyzer: Cobas Mira

Application: **PANCREATIC ISOAMYLASE FL** - Code PA F080 / F245 CH
 Preparation: REAGENT A - INSTALL AS "REAGENT"
 REAGENT B - INSTALL AS "START REAGENT 1"
 Storage: REFRIGERATE AT 2-8°C
 Stability: UP TO DECLARED EXPIRATION ON REAGENT LABEL

PROGRAM**GENERAL**

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S-SR1
 CALIBRATION MODE: FACTOR
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 405 nm
 DECIMAL POSITION: 0
 UNIT: U/L

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 6 µl DIL: 10.0 µl

REAGENT CYCLE: 1
 VOL: 230 µl

START REAGENT 1: CYCLE: 10
 VOL: 50 µl DIL: 10.0 µl

CALCULATION

SAMPLE LIMIT: 0.1500 ΔA
 POINT: T1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 U/L
 HIGH: 2500 U/L
 NORM. RANGE LOW: 13 U/L
 HIGH: 53 U/L

NUMBER OF STEPS: 1

CALC. STEP A: KINETIC
 READINGS FIRST: 12 LAST: 18
 REACTION LIMIT: 0.4000 ΔA
 POINT: 12

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: -0.0500 A
 HIGH: 0.2000 A
 BLANK RANGE LOW: -0.0100 A
 HIGH: 0.0100 A

FACTOR: 10460

CONTROL

CS1 POS: UTIL. DEF.
 CS2 POS: UTIL. DEF.
 CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **BICARBONATE FL** - Codes BR F060 / F245 / F400 CH
 Preparation: AS INDICATED IN INSERT SHEET
 Storage: REFRIGERATE AT 2-8°C
 Stability: 30 DAYS

PROGRAM**GENERAL**

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 405 nm
 DECIMAL POSITION: 1
 UNIT: mmol/l

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 3.0 µl DIL: 10.0 µl

REAGENT CYCLE: 1
 VOL: 300 µl

CALCULATION

SAMPLE LIMIT: 0.1000 ΔA
 POINT: T1
 REAC. DIRECTION: DECREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0.00 mmol/l
 HIGH: 50.00 mmol/l
 NORM. RANGE LOW: 8.60 mmol/l
 HIGH: 10.30 mmol/l

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
 READINGS FIRST: CB LAST: 6

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: 0.4000 A
 HIGH: 1.6000 A
 BLANK RANGE LOW: -0.2000 A
 HIGH: 0.2000 A

CALIBRATOR CUP POS: 1
 CAL-1: UTIL. DEF

REPLICATE: TRIPL.
 DEVIATION: 6%

CONTROL

CS1 POS: UTIL. DEF.

CS2 POS: UTIL. DEF.

CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **BILIRUBIN TOTAL FL** - Code DT F125 / F500 CH
 Preparation: REAGENT A - INSTALL AS "REAGENT"
 REAGENT B - INSTALL AS "START REAGENT 1"
 Storage: REFRIGERATE AT 2-8°C
 Stability: UP TO DECLARED EXPIRATION ON REAGENT LABEL

PROGRAM

GENERAL

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S-SR1
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 500 nm
 DECIMAL POSITION: 2
 UNIT: mg/dl

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 10.0 µl DIL: 10.0 µl

REAGENT CYCLE: 1
 VOL: 200 µl

START REAGENT 1: CYCLE: 13
 VOL: 50 µl DIL: 10.0 µl

CALCULATION

SAMPLE LIMIT: 0.5000 ΔA
 POINT: 1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 mg/dl
 HIGH: 20 mg/dl
 NORM. RANGE LOW: 0.2 mg/dl
 HIGH: 1.0 mg/dl

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
 READINGS FIRST: 5 LAST: 14
 REACTION LIMIT: NO

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: NO
 HIGH: NO
 BLANK RANGE LOW: NO
 HIGH: NO

CALIBRATOR CUP POS: DEF
 CAL-1: UTIL. DEF

REPLICATE: TRIPL.
 DEVIATION: 3%

CONTROL

CS1 POS: UTIL. DEF.
 CS2 POS: UTIL. DEF.
 CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **BILIRUBIN DIRECT FL** - Code DD F125 / F500 CH
 Preparation: REAGENT A - INSTALL AS "REAGENT"
 REAGENT B - INSTALL AS "START REAGENT 1"
 Storage: REFRIGERATE AT 2-8°C
 Stability: UP TO DECLARED EXPIRATION ON REAGENT LABEL

PROGRAM

GENERAL

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S-SR1
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 550 nm
 DECIMAL POSITION: 2
 UNIT: mg/dl

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 10.0 µl DIL: 10.0 µl

REAGENT CYCLE: 1
 VOL: 200 µl

START REAGENT 1: CYCLE: 12
 VOL: 50 µl DIL: 10.0 µl

CALCULATION

SAMPLE LIMIT: 0.5000 ΔA
 POINT: 1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 mg/dl
 HIGH: 13 mg/dl
 NORM. RANGE LOW: 0.0 mg/dl
 HIGH: 0.2 mg/dl

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
 READINGS FIRST: 11 LAST: 25
 REACTION LIMIT: NO

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: NO
 HIGH: NO
 BLANK RANGE LOW: NO
 HIGH: NO

CALIBRATOR CUP POS: DEF
 CAL-1: UTIL. DEF

REPLICATE: TRIPL.
 DEVIATION: 3%

CONTROL

CS1 POS: UTIL. DEF.
 CS2 POS: UTIL. DEF.
 CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **CALCIUM** - Code CA 0305 / 0505 CH
 Preparation: AS INDICATED IN INSERT SHEET
 Storage: REFRIGERATE AT 2-8°C
 Stability: 10 DAYS

PROGRAM**GENERAL**

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: SELECT
 AFTER TESTS: AMYLASE
 WAVELENGTH: 550 nm
 DECIMAL POSITION: 2
 UNIT: mg/dl

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 4.0 µl DIL: 20.0 µl

REAGENT CYCLE: 1
 VOL: 250 µl

CALCULATION

SAMPLE LIMIT: 0.1000 ΔA
 POINT: CB
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0.00 mg/dl
 HIGH: 20.00 mg/dl
 NORM. RANGE LOW: 8.60 mg/dl
 HIGH: 10.30 mg/dl

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
 READINGS FIRST: CB LAST: 4

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: -0.0500 A
 HIGH: 0.5000 A
 BLANK RANGE LOW: -0.0100 A
 HIGH: 0.5000 A

CALIBRATOR CUP POS: 1
 CAL-1: UTIL. DEF

REPLICATE: TRIPL.
 DEVIATION: 6%

CONTROL

CS1 POS: UTIL. DEF.
 CS2 POS: UTIL. DEF.
 CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **CHLORIDE** - Code CL 0100 / 0500 CH
 Preparation: LIQUID READY TO USE SINGLE REAGENT
 Storage: ROOM TEMPERATURE (2-30°C)
 Stability: UP TO DECLARED EXPIRATION ON REAGENT LABEL

PROGRAM

GENERAL

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: BEFORE
 WAVELENGTH: 500 nm
 DECIMAL POSITION: 0
 UNIT: mEq/L

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE VOL: 3.0 µl CYCLE: 1
 DIL: 20.0 µl

REAGENT VOL: 400 µl CYCLE: 1

CALCULATION

SAMPLE LIMIT: 0.1000 ΔA
 POINT: T1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 mEq/L
 HIGH: 200 mEq/L
 NORM. RANGE LOW: 98 mEq/L
 HIGH: 110 mEq/L

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
 READINGS FIRST: CB LAST: 10

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: NO
 HIGH: NO
 BLANK RANGE LOW: NO
 HIGH: NO

CALIBRATOR CUP POS: 1
 CAL-1: UTIL. DEF

REPLICATE: TRIPL.
 DEVIATION: 6%

CONTROL

CS1 POS: UTIL. DEF.

CS2 POS: UTIL. DEF.

CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **CHOLESTEROL FL** - Codes CT F100 / F400 / 100F / 150F CH
 Preparation: LIQUID READY TO USE SINGLE REAGENT
 Storage: REFRIGERATE AT 2-8°C
 Stability: UP TO DECLARED EXPIRATION ON REAGENT LABEL

PROGRAM**GENERAL**

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 500 nm
 DECIMAL POSITION: 0
 UNIT: mg/dl

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 3.0 µl DIL: 20.0 µl

REAGENT CYCLE: 1
 VOL: 250 µl

CALCULATION

SAMPLE LIMIT: 0.1000 ΔA
 POINT: T1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 mg/dl
 HIGH: 700 mg/dl
 NORM. RANGE LOW: 140 mg/dl
 HIGH: 200 mg/dl

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
 READINGS FIRST: T1 LAST: 12

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: -0.0500 A
 HIGH: 0.1600 A
 BLANK RANGE LOW: -0.0500 A
 HIGH: 0.1000 A

CALIBRATOR CUP POS: 1
 CAL-1: UTIL. DEF

REPLICATE: TRIPL.
 DEVIATION: 6%

CONTROL

CS1 POS: UTIL. DEF.

CS2 POS: UTIL. DEF.

CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **HDL-DIRECT FL** - Code HD F080 / F245 / F400 CH
 Preparation: REAGENT A - INSTALL AS "REAGENT"
 REAGENT B - INSTALL AS "START REAGENT 1"
 Storage: REFRIGERATE AT 2-8°C
 Stability: 60 DAYS ON BOARD

PROGRAM

GENERAL

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S-SR1
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 600nm
 DECIMAL POSITION: 0
 UNIT: mg/dl

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 3.0 µl DIL: 20.0 µl

REAGENT CYCLE: 1
 VOL: 180 µl

START REAGENT 1 CYCLE: 12
 VOL: 60 µl DIL: 10.0 µl

CALCULATION

SAMPLE LIMIT: 0.3000 ΔA
 POINT: 2
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 mg/dl
 HIGH: 220 mg/dl
 NORM. RANGE LOW: 35.3 mg/dl
 HIGH: 79.5 mg/dl

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
 READINGS FIRST: 10 LAST: 24

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: -0.0500 A
 HIGH: 0.1000 A
 BLANK RANGE LOW: -0.0500 A
 HIGH: 0.0500 A

CALIBRATOR CUP POS: 1
 CAL-1: UTIL. DEF

REPLICATE: TRIPL.
 DEVIATION: 6%

CONTROL

CS1 POS: UTIL. DEF.
 CS2 POS: UTIL. DEF.
 CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **LDL-DIRECT FL** - Code DL F080 CH
 Preparation: REAGENT A - INSTALL AS "REAGENT"
 REAGENT B - INSTALL AS "START REAGENT 1"
 Storage: REFRIGERATE AT 2-8°C
 Stability: 30 DAYS ON BOARD

PROGRAM**GENERAL**

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S-SR1
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 600nm
 DECIMAL POSITION: 0
 UNIT: mg/dl

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE VOL: 3.0 µl CYCLE: 1
 DIL: 20.0 µl
 REAGENT VOL: 180 µl CYCLE: 1
 START REAGENT 1 VOL: 60 µl CYCLE: 12
 DIL: 10.0 µl

CALCULATION

SAMPLE LIMIT: 0.3000 ΔA
 POINT: 2
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 mg/dl
 HIGH: 400 mg/dl
 NORM. RANGE LOW: 76 mg/dl
 HIGH: 218 mg/dl

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
 READINGS FIRST: 10 LAST: 24

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: -0.0500 A
 HIGH: 0.1000 A
 BLANK RANGE LOW: -0.0500 A
 HIGH: 0.0500 A

CALIBRATOR CUP POS: 1
 CAL-1: UTIL. DEF

REPLICATE: TRIPL.
 DEVIATION: 6%

CONTROL

CS1 POS: UTIL. DEF.
 CS2 POS: UTIL. DEF.
 CS3 POS: UTIL. DEF.

Analyzer: Cobas MiraApplication: **CHOLINESTERASE FL (DGKC)** - Codes CH F096 / F245 CHPreparation: REAGENT A - INSTALL AS "REAGENT"
REAGENT B - INSTALL AS "START REAGENT 1"

Storage: REFRIGERATE AT 2-8°C

Stability: UP TO DECLARED EXPIRATION ON REAGENT LABEL

PROGRAM**GENERAL**MEASUREMENT MODE: ABSORB
REACTION MODE: R-S-SR1
CALIBRATION MODE: FACTOR
REAGENT BLANK: REAG/DIL
CLEANER: NO
WAVELENGTH: 405 nm
DECIMAL POSITION: 0
UNIT: U/L**ANALYSIS**SAMPLE DIL. NAME: H2O
POST. DIL. FACTOR: 2.00
CONC. FACTOR: NO

SAMPLE CYCLE: 1
VOL: 4.8 µl DIL: 10.0 µl

REAGENT CYCLE: 1
VOL: 230 µl

START REAGENT 1: CYCLE: 4
VOL: 38 µl DIL: 10.0 µl**CALCULATION**SAMPLE LIMIT: 0.1500 ΔA
POINT: T1
REAC. DIRECTION: DECREASE
CHECK: ONCONVERS. FACTOR: 1.1000
OFFSET: 0.0000TEST RANGE LOW: 0 U/L
HIGH: 25000 U/L
NORM. RANGE LOW: 5600 U/L
HIGH: 11200 U/L

NUMBER OF STEPS: 1

CALC. STEP A: KINETIC
READINGS FIRST: 7 LAST: 10
REACTION LIMIT: NO**CALIBRATION**CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
REAG. RANGE LOW: 0.8000 A
HIGH: 2.0000 A
BLANK RANGE LOW: -0.0100 A
HIGH: 0.0100 A

FACTOR: 99000**CONTROL**CS1 POS: UTIL. DEF.

CS2 POS: UTIL. DEF.

CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **CK-NAC FL** - Codes CK F060 / F120 / F245 CH
 Preparation: AS INDICATED IN INSERT SHEET
 Storage: REFRIGERATE AT 2-8°C
 Stability: 30 DAYS

PROGRAM

GENERAL

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: FACTOR
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 340 nm
 DECIMAL POSITION: 0
 UNIT: U/L

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

 SAMPLE CYCLE: 1
 VOL: 10.0 µl DIL: 10.0 µl

 REAGENT CYCLE: 1
 VOL: 240 µl

CALCULATION

SAMPLE LIMIT: 0.3000 ΔA
 POINT: T1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 U/L
 HIGH: 2000 U/L
 NORM. RANGE LOW: 24 U/L
 HIGH: 204 U/L

NUMBER OF STEPS: 1

CALC. STEP A: KINETIC
 READINGS FIRST: 6 LAST: 12
 REACTION LIMIT: 0.6000 ΔA
 POINT: T1

CALIBRATION

CALIB. INTERVAL: ON REQUEST

 REAGENT BLANK
 REAG. RANGE LOW: -0.0100 A
 HIGH: 0.6000 A
 BLANK RANGE LOW: -0.0100 A
 HIGH: 0.0100 A

 FACTOR: 6878

CONTROL

CS1 POS: UTIL. DEF.

 CS2 POS: UTIL. DEF.

 CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **CK-MB FL** - Code MB F060 / F120 CH
 Preparation: AS INDICATED IN INSERT SHEET
 Storage: REFRIGERATE AT 2-8°C
 Stability: 14 DAYS

PROGRAM**GENERAL**

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: FACTOR
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 340 nm
 DECIMAL POSITION: 0
 UNIT: U/L

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 10.0 µl DIL: 10.0 µl

REAGENT CYCLE: 1
 VOL: 240 µl

CALCULATION

SAMPLE LIMIT: 0.3000 ΔA
 POINT: T1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0.0 U/L
 HIGH: 1000.0 U/L
 NORM. RANGE LOW: 0.0 U/L
 HIGH: 24.0 U/L

NUMBER OF STEPS: 1

CALC. STEP A: KINETIC
 READINGS FIRST: 10 LAST: 20
 REACTION LIMIT: 0.6000 ΔA
 POINT: T1

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: -0.0100 A
 HIGH: 0.6000 A
 BLANK RANGE LOW: -0.0100 A
 HIGH: 0.0100 A

FACTOR: 13750

CONTROL

CS1 POS: UTIL. DEF.
 CS2 POS: UTIL. DEF.
 CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **COPPER** - Code CU 0100 CH
 Preparation: AS INDICATED IN INSERT SHEET
 Storage: DO NOT REFRIGERATE!
 Stability: 14 DAYS

PROGRAM

GENERAL

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 600 nm
 DECIMAL POSITION: 0
 UNIT: microg/dl

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 20.0 µl DIL: 20.0 µl

REAGENT CYCLE: 1
 VOL: 280 µl

CALCULATION

SAMPLE LIMIT: 0.1000 ΔA
 POINT: T1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 microg/dl
 HIGH: 500 microg/dl
 NORM. RANGE LOW: 70 microg/dl
 HIGH: 140 microg/dl

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
 READINGS FIRST: CB LAST: 10

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: -0.0500 A
 HIGH: 0.5000 A
 BLANK RANGE LOW: -0.0500 A
 HIGH: 0.5000 A

CALIBRATOR CUP POS: 1
 CAL-1: UTIL. DEF

REPLICATE: TRIPL.
 DEVIATION: 6%

CONTROL

CS1 POS: UTIL. DEF.

CS2 POS: UTIL. DEF.

CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **CREATININE** - Code CR 0500 / 1000 / 1500 CH
 Preparation: AS INDICATED IN INSERT SHEET
 Storage: ROOM TEMPERATURE (2-30°C)
 Stability: 14 DAYS

PROGRAM

GENERAL

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 500 nm
 DECIMAL POSITION: 2
 UNIT: mg/dl

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 15.0 µl DIL: 20.0 µl

REAGENT CYCLE: 1
 VOL: 250 µl

CALCULATION

SAMPLE LIMIT: 0.3500 ΔA
 POINT: T1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0.00 mg/dl
 HIGH: 20.0 mg/dl
 NORM. RANGE LOW: 0.70 mg/dl
 HIGH: 1.20 mg/dl

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
 READINGS FIRST: 2 LAST: 3

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: -0.1000 A
 HIGH: 0.3500 A
 BLANK RANGE LOW: -0.0300 A
 HIGH: 0.0300 A

CALIBRATOR CUP POS: 1
 CAL-1: UTIL. DEF

REPLICATE: TRIPL.
 DEVIATION: 6%

CONTROL

CS1 POS: UTIL. DEF.

CS2 POS: UTIL. DEF.

CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **GAMMA-GT FL** - Codes GT F080 / F245 / F400 / F600 CH
 Preparation: AS INDICATED IN INSERT SHEET
 Storage: REFRIGERATE AT 2-8°C
 Stability: 60 DAYS

PROGRAM**GENERAL**

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: FACTOR
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 405 nm
 DECIMAL POSITION: 0
 UNIT: U/L

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 5.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 25.0 µl DIL: 10.0 µl

REAGENT CYCLE: 1
 VOL: 240 µl

CALCULATION

SAMPLE LIMIT: NO

REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 U/L
 HIGH: 800 U/L
 NORM. RANGE LOW: 0 U/L
 HIGH: 50 U/L

NUMBER OF STEPS: 1

CALC. STEP A: KINETIC
 READINGS FIRST: 4 LAST: 8
 REACTION LIMIT: 0.4300 ΔA
 POINT: T1

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: -0.0700 A
 HIGH: 0.8000 A
 BLANK RANGE LOW: -0.0100 A
 HIGH: 0.0200 A

FACTOR: 2130

CONTROL

CS1 POS: UTIL. DEF.
 CS2 POS: UTIL. DEF.
 CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **GLUCOSE FL** - Codes GL F400 / 100F / 150F CH
 Preparation: LIQUID READY TO USE SINGLE REAGENT
 Storage: REFRIGERATE AT 2-8°C
 Stability: UP TO DECLARED EXPIRATION ON REAGENT LABEL

PROGRAM**GENERAL**

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 500 nm
 DECIMAL POSITION: 0
 UNIT: mg/dl

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 3.0 µl DIL: 20.0 µl

REAGENT CYCLE: 1
 VOL: 250 µl

CALCULATION

SAMPLE LIMIT: 0.1000 ΔA
 POINT: T1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 mg/dl
 HIGH: 500 mg/dl
 NORM. RANGE LOW: 70 mg/dl
 HIGH: 105 mg/dl

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
 READINGS FIRST: CB LAST: 15

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: -0.0500 A
 HIGH: 0.2000 A
 BLANK RANGE LOW: -0.0500 A
 HIGH: 0.2000 A

CALIBRATOR CUP POS: 1
 CAL-1: UTIL. DEF

REPLICATE: TRIPL.
 DEVIATION: 6%

CONTROL

CS1 POS: UTIL. DEF.

CS2 POS: UTIL. DEF.

CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **GLUCOSE UV FL** - Code GL F601 CH
 Preparation: AS INDICATED IN INSERT SHEET
 Storage: REFRIGERATE AT 2-8°C
 Stability: 60 DAYS

PROGRAM**GENERAL**

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 340 nm
 DECIMAL POSITION: 0
 UNIT: mg/dl

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 3.0 µl DIL: 20.0 µl

REAGENT CYCLE: 1
 VOL: 250 µl

CALCULATION

SAMPLE LIMIT: 0.1000 ΔA
 POINT: T1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 mg/dl
 HIGH: 700 mg/dl
 NORM. RANGE LOW: 70 mg/dl
 HIGH: 105 mg/dl

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
 READINGS FIRST: CB LAST: 15

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: -0.0500 A
 HIGH: 0.2000 A
 BLANK RANGE LOW: -0.0500 A
 HIGH: 0.2000 A

CALIBRATOR CUP POS: 1
 CAL-1: UTIL. DEF

REPLICATE: TRIPL.
 DEVIATION: 6%

CONTROL

CS1 POS: UTIL. DEF.

CS2 POS: UTIL. DEF.

CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **GOT/AST FL** - Codes GO F080 / F245 / F400 / F600 CH
 Preparation: AS INDICATED IN INSERT SHEET
 Storage: REFRIGERATE AT 2-8°C
 Stability: 30 DAYS

PROGRAM**GENERAL**

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: FACTOR
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 340 nm
 DECIMAL POSITION: 0
 UNIT: U/L

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 5.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 25.0 µl DIL: 10.0 µl

REAGENT CYCLE: 1
 VOL: 240 µl

CALCULATION

SAMPLE LIMIT: NO

REAC. DIRECTION: DECREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 U/L
 HIGH: 440 U/L
 NORM. RANGE LOW: 0 U/L
 HIGH: 35 U/L

NUMBER OF STEPS: 1

CALC. STEP A: KINETIC
 READINGS FIRST: 4 LAST: 8
 REACTION LIMIT: 0.2000 ΔA
 POINT: T1

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: 0.4000 A
 HIGH: 1.6000 A
 BLANK RANGE LOW: -0.2000 A
 HIGH: 0.2000 A

FACTOR: 2910

CONTROL

CS1 POS: UTIL. DEF.
 CS2 POS: UTIL. DEF.
 CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **GPT/ALT FL** - Codes GP F080 / F245 / F400 / F600 CH
 Preparation: AS INDICATED IN INSERT SHEET
 Storage: REFRIGERATE AT 2-8°C
 Stability: 30 DAYS

PROGRAM**GENERAL**

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: FACTOR
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 340 nm
 DECIMAL POSITION: 0
 UNIT: U/L

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 5.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 25.0 µl DIL: 10.0 µl

REAGENT CYCLE: 1
 VOL: 240 µl

CALCULATION

SAMPLE LIMIT: NO

REAC. DIRECTION: DECREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 U/L
 HIGH: 440 U/L
 NORM. RANGE LOW: 0 U/L
 HIGH: 45 U/L

NUMBER OF STEPS: 1

CALC. STEP A: KINETIC
 READINGS FIRST: 4 LAST: 8
 REACTION LIMIT: 0.2000 ΔA
 POINT: T1

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: 0.4000 A
 HIGH: 1.6000 A
 BLANK RANGE LOW: -0.2000 A
 HIGH: 0.2000 A

FACTOR: 2910

CONTROL

CS1 POS: UTIL. DEF.
 CS2 POS: UTIL. DEF.
 CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **IRON FZ** - Code FE F245 / F400 CH
 Preparation: REAGENT A - INSTALL AS "REAGENT"
 REAGENT B - MIX AS INDICATED IN INSERT SHEET.
 INSTALL AS "START REAGENT 1"
 Storage: REFRIGERATE AT 2-8°C
 Stability: UP TO DECLARED EXPIRATION ON REAGENT LABEL

PROGRAM

GENERAL

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S-SR1
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 550nm
 DECIMAL POSITION: 0
 UNIT: µg/dl

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

 SAMPLE CYCLE: 1
 VOL: 50.0 µl DIL: 20.0 µl

 REAGENT CYCLE: 1
 VOL: 200 µl

 START REAGENT 1 CYCLE: 12
 VOL: 50 µl DIL: 20.0 µl

CALCULATION

SAMPLE LIMIT: 0.3000 ΔA
 POINT: 11
 REAC. DIRECTION: INCREASE
 CHECK: ON

 CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

 TEST RANGE LOW: 0 mg/dl
 HIGH: 1000 mg/dl
 NORM. RANGE LOW: 59 mg/dl
 HIGH: 158 mg/dl

NUMBER OF STEPS: 1

 CALC. STEP A: ENDPOINT
 READINGS FIRST: 11 LAST: 15

CALIBRATION

CALIB. INTERVAL: ON REQUEST

 REAGENT BLANK
 REAG. RANGE LOW: -0.0500 A
 HIGH: 0.1000 A
 BLANK RANGE LOW: -0.0500 A
 HIGH: 0.0500 A

 CALIBRATOR CUP POS: 1
 CAL-1: UTIL. DEF

 REPLICATE: TRIPL.
 DEVIATION: 6%

CONTROL

CS1 POS: UTIL. DEF.
 CS2 POS: UTIL. DEF.
 CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **LDH FL** - Codes LD F060 / F120 / F245 CH
 Preparation: AS INDICATED IN INSERT SHEET
 Storage: REFRIGERATE AT 2-8°C
 Stability: 30 DAYS

PROGRAM

GENERAL

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: FACTOR
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 340 nm
 DECIMAL POSITION: 0
 UNIT: U/L

CALIBRATION

CALIB. INTERVAL: ON REQUEST
 REAGENT BLANK
 REAG. RANGE LOW: 0.3000 A
 HIGH: 1.0000 A
 BLANK RANGE LOW: -0.2000 A
 HIGH: 0.2000 A
 FACTOR: 26715

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO
 SAMPLE CYCLE: 1
 VOL: 3.0 µl DIL: 20.0 µl
 REAGENT CYCLE: 1
 VOL: 280 µl

CONTROL

CS1 POS: UTIL. DEF.
 CS2 POS: UTIL. DEF.
 CS3 POS: UTIL. DEF.

CALCULATION

SAMPLE LIMIT: 0.2000 ΔA
 POINT: T1

REAC. DIRECTION: DECREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 U/L
 HIGH: 4000 U/L
 NORM. RANGE LOW: 225 U/L
 HIGH: 450 U/L

NUMBER OF STEPS: 1

CALC. STEP A: KINETIC
 READINGS FIRST: 3 LAST: 7
 REACTION LIMIT: 0.2000 ΔA
 POINT: T1

Analyzer: Cobas Mira

Application: **LIPASE FL** - Code LP F060 CH
 Preparation: REAGENT A - INSTALL AS "REAGENT"
 REAGENT B - INSTALL AS "START REAGENT 1"
 Storage: REFRIGERATE AT 2-8°C

PROGRAM

GENERAL

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S-SR1
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: BEFORE
 WAVELENGTH: 550 nm
 DECIMAL POSITION: 1
 UNIT: U/L

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 2.5 µl DIL: 15.0 µl

REAGENT CYCLE: 1
 VOL: 240 µl

START REAGENT 1: CYCLE: 4
 VOL: 60 µl DIL: 10.0 µl

CALCULATION

SAMPLE LIMIT: 0.5000 ΔA
 POINT: 1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 U/L
 HIGH: 250 U/L
 NORM. RANGE LOW: 0 U/L
 HIGH: 63 U/L

NUMBER OF STEPS: 1

CALC. STEP A: KINETIC
 READINGS FIRST: 6 LAST: 10
 REACTION LIMIT: NO

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: NO
 HIGH: NO
 BLANK RANGE LOW: NO
 HIGH: NO

CALIBRATOR CUP POS: DEF
 CAL-1: UTIL. DEF

REPLICATE: TRIPL.
 DEVIATION: 10%

CONTROL

CS1 POS: UTIL. DEF.
 CS2 POS: UTIL. DEF.
 CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **MAGNESIUM** - Code MG 0200 / 0500 CH
 Preparation: AS INDICATED IN INSERT SHEET
 Storage: REFRIGERATE AT 2-8°C (STABLE 90 days)
 ROOM TEMPERATURE 15-30°C (STABLE 30 days)

PROGRAM

GENERAL

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: SELECT
 AFTER TESTS: CHOL TRIG ALP
 WAVELENGTH: 500nm
 DECIMAL POSITION: 2
 UNIT: mEq/L

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: NO
 CONC. FACTOR: NO

 SAMPLE CYCLE: 1
 VOL: 3.0 µl DIL: 20.0 µl

 REAGENT CYCLE: 1
 VOL: 280 µl

CALCULATION

SAMPLE LIMIT: 0.1000 ΔA
 POINT: T1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0.00 mEq/L
 HIGH: 8.00 mEq/L
 NORM. RANGE LOW: 1.30 mEq/L
 HIGH: 2.10 mEq/L

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
 READINGS FIRST: CB LAST: 4

CALIBRATION

CALIB. INTERVAL: ON REQUEST

 REAGENT BLANK
 REAG. RANGE LOW: -0.0500 A
 HIGH: 0.5000 A
 BLANK RANGE LOW: -0.0500 A
 HIGH: 0.5000 A

 CALIBRATOR CUP POS: 1
 CAL-1: UTIL. DEF

 REPLICATE: TRIPL.
 DEVIATION: 5%

CONTROL

CS1 POS: UTIL. DEF.
 CS2 POS: UTIL. DEF.
 CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **MAGNESIUM XL** - Code MX 0300 / 0500 CH
 Preparation: LIQUID READY TO USE SINGLE REAGENT
 Storage: REFRIGERATE AT 2-8°C
 Stability: UP TO DECLARED EXPIRATION ON REAGENT LABEL

PROGRAM**GENERAL**

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: SELECT
 AFTER TESTS: CHOL TRIG ALP
 WAVELENGTH: 550 nm
 DECIMAL POSITION: 2
 UNIT: mEq/L

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: NO
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 3.0 µl DIL: 20.0 µl

REAGENT CYCLE: 1
 VOL: 300 µl

CALCULATION

SAMPLE LIMIT: 0.1000 ΔA
 POINT: T1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0.00 mEq/L
 HIGH: 6.00 mEq/L
 NORM. RANGE LOW: 1.30 mEq/L
 HIGH: 2.10 mEq/L

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
 READINGS FIRST: CB LAST: 12

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: -0.0500 A
 HIGH: 0.5000 A
 BLANK RANGE LOW: -0.0500 A
 HIGH: 0.5000 A

CALIBRATOR CUP POS: 1
 CAL-1: UTIL. DEF

REPLICATE: TRIPL.
 DEVIATION: 5%

CONTROL

CS1 POS: UTIL. DEF.

CS2 POS: UTIL. DEF.

CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **PHOSPHORUS** - Code PH 0100 / 0500 CH
 Preparation: LIQUID READY TO USE SINGLE REAGENT
 Storage: REFRIGERATE AT 2-8°C
 Stability: UP TO DECLARED EXPIRATION ON REAGENT LABEL

PROGRAM

GENERAL

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: SELECT
 AFTER TESTS: LDH - URIC - GLUC
 WAVELENGTH: 340 nm
 DECIMAL POSITION: 2
 UNIT: mg/dl

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

 SAMPLE CYCLE: 1
 VOL: 3.0 µl DIL: 20.0 µl

 REAGENT CYCLE: 1
 VOL: 250 µl

CALCULATION

SAMPLE LIMIT: 0.2000 ΔA
 POINT: T1
 REAC. DIRECTION: INCREASE
 CHECK: ON

 CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0.0 mg/dl
 HIGH: 20.0 mg/dl
 NORM. RANGE LOW: 2.5 mg/dl
 HIGH: 4.5 mg/dl

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
 READINGS FIRST: CB LAST: 12

CALIBRATION

CALIB. INTERVAL: ON REQUEST

 REAGENT BLANK
 REAG. RANGE LOW: -0.0500 A
 HIGH: 0.5000 A
 BLANK RANGE LOW: -0.0500 A
 HIGH: 0.5000 A

 CALIBRATOR CUP POS: 1
 CAL-1: UTIL. DEF

 REPLICATE: TRIPL.
 DEVIATION: 6%

CONTROL

CS1 POS: UTIL. DEF.
 CS2 POS: UTIL. DEF.
 CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **PROTEINS HS** - Code HS 0100 / 0500 CH
 Preparation: LIQUID READY TO USE SINGLE REAGENT
 Storage: ROOM TEMPERATURE (2-30°C)
 Stability: UP TO DECLARED EXPIRATION ON REAGENT LABEL

PROGRAM**GENERAL**

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 600 nm
 DECIMAL POSITION: 0
 UNIT: mg/dl

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 2.5 µl DIL: 20.0 µl

REAGENT CYCLE: 1
 VOL: 300 µl

CALCULATION

SAMPLE LIMIT: 0.3500 ΔA
 POINT: T1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 mg/dl
 HIGH: 500 mg/dl
 NORM. RANGE LOW: 28 mg/dl
 HIGH: 141 mg/dl

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
 READINGS FIRST: T1 LAST: 12

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: NO
 HIGH: NO
 BLANK RANGE LOW: NO
 HIGH: NO

CALIBRATOR CUP POS: 1
 CAL-1: UTIL. DEF

REPLICATE: TRIPL.
 DEVIATION: 10%

CONTROL

CS1 POS: UTIL. DEF.

CS2 POS: UTIL. DEF.

CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **PROTEINS (TOTAL)** - Codes TP 0100 / 0500 / 1000 / 1500 CH
 Preparation: LIQUID READY TO USE SINGLE REAGENT
 Storage: ROOM TEMPERATURE (2-30°C)
 Stability: UP TO DECLARED EXPIRATION ON REAGENT LABEL

PROGRAM

GENERAL

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 550 nm
 DECIMAL POSITION: 1
 UNIT: g/dl

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 3.0 µl DIL: 20.0 µl

REAGENT CYCLE: 1
 VOL: 280 µl

CALCULATION

SAMPLE LIMIT: 0.3500 ΔA
 POINT: T1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0.0 g/dl
 HIGH: 12.0 g/dl
 NORM. RANGE LOW: 6.30 g/dl
 HIGH: 8.30 g/dl

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
 READINGS FIRST: CB LAST: 12

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: -0.0100 A
 HIGH: 0.5000 A
 BLANK RANGE LOW: -0.0100 A
 HIGH: 0.5000 A

CALIBRATOR CUP POS: 1
 CAL-1: UTIL. DEF

REPLICATE: TRIPL.
 DEVIATION: 6%

CONTROL

CS1 POS: UTIL. DEF.

CS2 POS: UTIL. DEF.

CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **TRIGLYCERIDES FL** - Codes TR F100 / F400 / 100F / 150F CH
 Preparation: LIQUID READY TO USE SINGLE REAGENT
 Storage: REFRIGERATE AT 2-8°C
 Stability: UP TO DECLARED EXPIRATION ON REAGENT LABEL

PROGRAM

GENERAL

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 500 nm
 DECIMAL POSITION: 0
 UNIT: mg/dl

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 3.0 µl DIL: 20.0 µl

REAGENT CYCLE: 1
 VOL: 280 µl

CALCULATION

SAMPLE LIMIT: 0.1000 ΔA
 POINT: T1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 mg/dl
 HIGH: 1000 mg/dl
 NORM. RANGE LOW: 0 mg/dl
 HIGH: 200 mg/dl

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
 READINGS FIRST: CB LAST: 13

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: -0.0500 A
 HIGH: 0.1800 A
 BLANK RANGE LOW: -0.0500 A
 HIGH: 0.1800 A

CALIBRATOR CUP POS: 1
 CAL-1: UTIL. DEF

REPLICATE: TRIPL.
 DEVIATION: 6%

CONTROL

CS1 POS: UTIL. DEF.

CS2 POS: UTIL. DEF.

CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **UREA UV FL** - Codes AZ F080 / F245 / F400 / F600 CH
 Preparation: AS INDICATED IN INSERT SHEET
 Storage: REFRIGERATE AT 2-8°C
 Stability: 60 DAYS

PROGRAM

GENERAL

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: SELECT
 AFTER TESTS: IRON
 WAVELENGTH: 340 nm
 DECIMAL POSITION: 1
 UNIT: mg/dl

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

 SAMPLE CYCLE: 1
 VOL: 3.0 µl DIL: 20.0 µl

 REAGENT CYCLE: 1
 VOL: 280 µl

CALCULATION

SAMPLE LIMIT: 0.1000 ΔA
 POINT: T1
 REAC. DIRECTION: DECREASE
 CHECK: ON

 CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0.0 mg/dl
 HIGH: 300.0 mg/dl
 NORM. RANGE LOW: 10.0 mg/dl
 HIGH: 50.0 mg/dl

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
 READINGS FIRST: 3 LAST: 5

CALIBRATION

CALIB. INTERVAL: ON REQUEST

 REAGENT BLANK
 REAG. RANGE LOW: 0.4000 A
 HIGH: 1.4000 A
 BLANK RANGE LOW: -0.0300 A
 HIGH: 0.1000 A

 CALIBRATOR CUP POS: 1
 CAL-1: UTIL. DEF

 REPLICATE: TRIPL.
 DEVIATION: 6%

CONTROL

CS1 POS: UTIL. DEF.
 CS2 POS: UTIL. DEF.
 CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **URIC ACID T FL** - Code AU F402 CH
 Preparation: AS INDICATED IN INSERT SHEET
 Storage: REFRIGERATE AT 2-8°C
 Stability: 90 DAYS

PROGRAM

GENERAL

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 550 nm
 DECIMAL POSITION: 2
 UNIT: mg/dl

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 8.0 µl DIL: 20.0 µl

REAGENT CYCLE: 1
 VOL: 250 µl

CALCULATION

SAMPLE LIMIT: 0.1000 ΔA
 POINT: T1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0.00 mg/dl
 HIGH: 25.00 mg/dl
 NORM. RANGE LOW: 3.50 mg/dl
 HIGH: 7.20 mg/dl

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
 READINGS FIRST: CB LAST: 16

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: -0.0500 A
 HIGH: 0.1600 A
 BLANK RANGE LOW: -0.0500 A
 HIGH: 0.1000 A

CALIBRATOR CUP POS: 1
 CAL-1: UTIL. DEF

REPLICATE: TRIPL.
 DEVIATION: 6%

CONTROL

CS1 POS: UTIL. DEF.

CS2 POS: UTIL. DEF.

CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **URIC ACID AOX FL** - Code AX F100 / F250 / F060 CH

Preparation: REAGENT A - INSTALL AS "REAGENT"
REAGENT B - INSTALL AS "START REAGENT 1"

Storage: REFRIGERATE AT 2-8°C

PROGRAM

GENERAL

MEASUREMENT MODE: ABSORB
REACTION MODE: R-S-SR1
CALIBRATION MODE: SLOPE AVG
REAGENT BLANK: REAG/DIL
CLEANER: NO
WAVELENGTH: 550 nm
DECIMAL POSITION: 2
UNIT: mg/dl

ANALYSIS

SAMPLE DIL. NAME: H2O
POST. DIL. FACTOR: 2.00
CONC. FACTOR: NO

SAMPLE CYCLE: 1
VOL: 10.0 µl DIL: 10.0 µl

REAGENT CYCLE: 1
VOL: 190 µl

START REAGENT 1: CYCLE: 13
VOL: 50 µl DIL: 5.0 µl

CALCULATION

SAMPLE LIMIT: 0.5000 ΔA
POINT: 1
REAC. DIRECTION: INCREASE
CHECK: ON

CONVERS. FACTOR: 1.0000
OFFSET: 0.0000

TEST RANGE LOW: 0 mg/dl
HIGH: 35 mg/dl
NORM. RANGE LOW: 3.5 mg/dl
HIGH: 7.2 mg/dl

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
READINGS FIRST: 10 LAST: 25
REACTION LIMIT: NO

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
REAG. RANGE LOW: NO
HIGH: NO
BLANK RANGE LOW: NO
HIGH: NO

CALIBRATOR CUP POS: DEF
CAL-1: UTIL. DEF

REPLICATE: TRIPL.
DEVIATION: 3%

CONTROL

CS1 POS: UTIL. DEF.
CS2 POS: UTIL. DEF.
CS3 POS: UTIL. DEF.

Analyzer: Cobas Mira

Application: **ZINC** - Code ZN 0125 CH
 Preparation: AS INDICATED IN INSERT SHEET
 Storage: REFRIGERATE AT 2-8°C
 Stability: 14 DAYS

PROGRAM**GENERAL**

MEASUREMENT MODE: ABSORB
 REACTION MODE: R-S
 CALIBRATION MODE: SLOPE AVG
 REAGENT BLANK: REAG/DIL
 CLEANER: NO
 WAVELENGTH: 600 nm
 DECIMAL POSITION: 0
 UNIT: microg/dl

ANALYSIS

SAMPLE DIL. NAME: H2O
 POST. DIL. FACTOR: 2.00
 CONC. FACTOR: NO

SAMPLE CYCLE: 1
 VOL: 15.0 µl DIL: 20.0 µl

REAGENT CYCLE: 1
 VOL: 280 µl

CALCULATION

SAMPLE LIMIT: 0.1000 ΔA
 POINT: T1
 REAC. DIRECTION: INCREASE
 CHECK: ON

CONVERS. FACTOR: 1.0000
 OFFSET: 0.0000

TEST RANGE LOW: 0 microg/dl
 HIGH: 500 microg/dl
 NORM. RANGE LOW: 70 microg/dl
 HIGH: 150 microg/dl

NUMBER OF STEPS: 1

CALC. STEP A: ENDPOINT
 READINGS FIRST: CB LAST: 10

CALIBRATION

CALIB. INTERVAL: ON REQUEST

REAGENT BLANK
 REAG. RANGE LOW: -0.0500 A
 HIGH: 0.5000 A
 BLANK RANGE LOW: -0.0500 A
 HIGH: 0.5000 A

CALIBRATOR CUP POS: 1
 CAL-1: UTIL. DEF

REPLICATE: TRIPL.
 DEVIATION: 6%

CONTROL

CS1 POS: UTIL. DEF.

CS2 POS: UTIL. DEF.

CS3 POS: UTIL. DEF.