AUTOCAL H

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1) Occupational Safety and Health Standards: bloodborne pathogens. (29 CFR Part 1910.1030). Fed. Register.
2) Council Directive (2000/54/EC). Official Journal of the European Communities. No. L262 from Oct. 17, 2000.

AUTOCAL H

HUMAN BASED CALIBRATION SERUM

Intended use / summary

Autocal H is a human based multiparametric calibration serum. The concentrations and activities have been selected to ensure optimum calibration of automatic analyzers.

Package components

Calibrator (lyophilized) 5 x 3 ml / 10 x 3 ml

Reactive components in calibrator

Human serum with chemical additives and extracts of tissue of human and animal origin. The origin of biological additives is as follow:

| | |
|----------------------|------------------------------|
| Acid phosphatase | human prostate / potato |
| Alkaline phosphatase | human placenta (recombinant) |
| ALT/GPT | porcine heart |
| AST/GOT | porcine heart |
| Albumin | bovine plasma |
| Aldolase | rabbit muscle |
| Amylase, total | porcine pancreas |
| Amylase, pancreatic | porcine pancreas |
| Cholesterol | bovine plasma |
| Cholinesterase | human serum |
| Creatine kinase | rabbit muscle |
| Gamma-GT | porcine kidney |
| GLDH | bacterial recombinant |
| LDH | porcine heart |
| Lipase | human pancreas (recombinant) |
| Triglycerides | chicken egg yolk |

Calibration values

The calibration values were determined using the method mentioned in the enclosed value sheet. Determinations were performed under strictly standardized conditions on Roche analyzers using Roche and Chema reagents and Roche master calibrators. The calibration values were obtained via single assays performed in different laboratories in several independent series. The calibrator value specified is the mean of the values obtained.

Precautions and warnings

For in vitro diagnostic use only. Exercise the normal precautions required for handling all laboratory reagents. This product has been prepared exclusively from the blood of donors tested individually and shown by FDA approved methods to be free from HBsAg and antibodies to HCV and HIV. However as no test method can rule out the potential risk of infection with absolute certainty, the material should be treated just as carefully as a patient sample. In the event of exposure the directive of the responsible health authorities should be followed (1,2).

Disposal

This product is made to be used in professional laboratories. Please consult local regulations for a correct waste disposal. P501: Dispose of contents according to national/international regulations.

Handling

Carefully open one bottle of calibrator, avoiding the loss of lyophilizate, and pipette exactly 3.0 ml of distilled water. Carefully close the bottle and dissolve the content completely by occasional gentle swirling within 30 minutes. Avoid the formation of foam.

5 x 3 ml
10 x 3 ml

REF AT 0015 CH

LOT 611124

REF AT 0030 CH

2019-11

Storage and stability

Store at 2-8°C.

Criterion of the stability data specified by the manufacturer: recovery within ± 5% of initial value.

Stability of lyophilized calibrator serum at 2-8°C: up to the stated expiration date.

Stability of the components in the reconstituted calibrator:

| | |
|--------------------|----------------------------|
| at 15-25°C | 8 hours |
| at 2-8°C | 2 days |
| at (-15) - (-25)°C | 4 weeks (when frozen once) |

Exceptions:

Stability of acid phosphatase and prostatic acid phosphatase in the reconstituted calibrator (criterion: ± 10% of initial value):

| | |
|--------------------|----------------------------|
| at 15-25°C | 4 hours |
| at 2-8°C | 1 day |
| at (-15) - (-25)°C | 2 weeks (when frozen once) |

Stability of total bilirubin in the reconstituted calibrator, when stored protected from light:

| | |
|--------------------|----------------------------|
| at 15-25°C | 6 hours |
| at 2-8°C | 1 day |
| at (-15) - (-25)°C | 2 weeks (when frozen once) |

Stability of direct bilirubin in the reconstituted calibrator, when stored protected from light:

| | |
|--------------------|----------------------------|
| at 15-25°C | 3 hours |
| at 2-8°C | 8 hours |
| at (-15) - (-25)°C | 2 weeks (when frozen once) |

Test procedure

Use supplied materials as provided. Additional materials required: automatic micropipette, distilled/deionized water, general laboratory equipment.

Assay

Dispense the required volume into a sample cup and use as advised.

References

- Occupational Safety and Health Standards: bloodborne pathogens. (29 CFR Part 1910.1030). Fed. Register.
- Council Directive (2000/54/EC). Official Journal of the European Communities. No. L262 from Oct. 17, 2000.

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