

# PROTEINS HS

|            |            |
|------------|------------|
| HS 0100 CH | 2 x 50 ml  |
| HS 0500 CH | 4 x 125 ml |

## INTENDED USE

Reagent for quantitative in vitro determination of proteins in biological fluids.

## SUMMARY OF TEST

The renal excretion of moderate quantities of proteins could be an index of a beginner renal pathology, especially in diabetic subjects. The presence of light chain immunoglobulins could be an index of kidney myeloma.

The present method, based on the pigment pyrogallol red according to Watanabe, permits a reliable quantification also of proteins classes different than albumin.

## PRINCIPLE OF THE METHOD

Proteins combine with pyrogallol red to form a color complex, the absorbance of which is measured at 600 nm. Sodium dodecylsulphate is added to increase accuracy in measuring proteins other than albumine (Watanabe).

## KIT COMPONENTS

### For in vitro diagnostic use only.

The components of the kit are stable until expiration date on the label.

Keep away from direct light sources.

**HS R1**      **0100: 2 x 50 ml (liquid) blue cap**  
**0500: 4 x 125 ml (liquid) blue cap**

Composition: succinate buffer 0.05 M pH 2.50, pyrogallol red 0.04 mM, sodium molybdate 0.13 mM, sodium oxalate 1 mM, sodium benzoate 0.35 mM, SDS 0.1 mM.

**Standard:**      **proteins solution 100 mg/dl - 5 ml**

Store all components at 2-8°C.

## MATERIALS REQUIRED BUT NOT SUPPLIED

Current laboratory instrumentation. Spectrophotometer UV/VIS with thermostatic cuvette holder. Automatic micropipettes. Glass or high quality polystyrene cuvettes. Saline solution.

## REAGENT PREPARATION

Use reagent ready to use.

Stability: up to expiration date on labels at 2-8°C.

Stability since first opening of vials: preferably within 60 days at 2-8°C.

## PRECAUTIONS

**HS R1: Warning.** May cause damage to organs (H371). Do not breath vapours (P260). Wash hands thoroughly after handling (P264). Do not eat, drink or smoke when using this product (P270).



**Standard:** It is not classified as hazardous.

## SPECIMEN

Urine, cerebrospinal fluid. Stable 3 days at 2-8°C.

## TEST PROCEDURE

Wavelength: 600 nm (allowed 580 ÷ 620 nm)  
Lightpath: 1 cm  
Temperature: 25, 30 or 37°C

| dispense: | blank | standard | sample |
|-----------|-------|----------|--------|
| reagent   | 1 ml  | 1 ml     | 1 ml   |
| water     | 5 µl  | -        | -      |
| standard  | -     | 5 µl     | -      |
| sample    | -     | -        | 5 µl   |

Mix, incubate at 25, 30 or 37°C for 5 minutes.  
Read absorbances of standard (As) and samples (Ax) against reagent blank.

## RESULTS CALCULATION

Urine, cerebrospinal fluid:  
proteins mg/dl = Ax/As x 100 (standard value)

## EXPECTED VALUES

Cerebrospinal fluid: 14 - 45 mg/dl  
Urine: 28 - 141 mg/24 h

Each laboratory should establish appropriate reference intervals related to its population.

## QUALITY CONTROL

It is suggested to perform an internal quality control. For this purpose a reliable and compatible urine based control material is needed.

Please contact Customer Care for further information.

## TEST PERFORMANCE

### Linearity

the method is linear up to 500 mg/dl.  
If the limit value is exceeded, it is suggested to dilute sample 1+9 with distilled water and to repeat the test, multiplying the result by 10.

### Sensitivity/limit of detection (LOD)

the limit of detection is 0.3 mg/dl.

### Interferences

no interference was observed by the presence of:  
ascorbic acid      ≤ 500 mg/dl

### Precision

| intra-assay (n=10) | mean (mg/dl) | SD (mg/dl) | CV%  |
|--------------------|--------------|------------|------|
| sample 1           | 37.10        | 0.74       | 2.00 |
| sample 2           | 103.40       | 1.27       | 1.20 |

| inter-assay (n=20) | mean (mg/dl) | SD (mg/dl) | CV%  |
|--------------------|--------------|------------|------|
| sample 1           | 38.01        | 0.79       | 2.00 |
| sample 2           | 100.09       | 2.46       | 2.00 |

### Methods comparison

a comparison between Chema and a commercially available product gave the following results:

Proteins HS Chema = x  
Proteins competitor = y  
n = 88

$$y = 0.97x - 0.54 \text{ mg/dl} \quad r^2 = 0.978$$

## WASTE DISPOSAL

This product is made to be used in professional laboratories.

P501: Dispose of contents according to national/international regulations.

## REFERENCES

Watanabe et al. - Clin.Chem. 32/8, 1551-1554 (1986).

## MANUFACTURER

Chema Diagnostica  
Via Campania 2/4  
60030 Monsano (AN) - ITALY - UE  
phone +39 0731 605064  
fax +39 0731 605672  
e-mail: mail@chema.com  
website: http://www.chema.com

## SYMBOLS

|  |                                    |
|--|------------------------------------|
|  | in vitro diagnostic medical device |
|  | batch code                         |
|  | catalogue number                   |
|  | temperature limit                  |
|  | use by date                        |
|  | caution                            |
|  | consult instructions for use       |