CRP CALIBRATOR

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Bibliografia

"The certification of the mass concentration of C-Reactive Protein in human serum - Certified Reference Material ERM-DA474/IFCC", European Community report, 2011

Intended use

CRP Calibrator is composed of a set of sera, aimed at the preparation of calibration curves for quantitative determination of C-Reactive Protein (CRP) in human serum, through immunoturbidimetric technique.

Principle of the method

C-Reactive Protein (CRP) selectively reacts with an anti-CRP antibody and forms an immunocomplex. The produced turbidity is proportional to the concentration of CRP in the sample, and can be measured at the wavelength of 340 nm.

Package components

Calibrators (liquid) 5 x 1 ml

Composition

CRP Calibrator is a liquid product. It is a synthetic biological fluid containing stabilizers, sodium azide (< 0,1%) as preservative and human CRP.

Assigned values

CRP values have been assigned to ERM/IFCC/474.

Storage and stability

Store at 2-8°C.

The calibrators are ready for use; after vial opening, they are stable until expiration date if stored tightly closed at 2-8°C. Do not freeze.

Test procedure

Use supplied materials as provided.

Additional materials required: automatic micropipette, distilled/deionized water, general laboratory equipment.

Precautions and warnings

For in vitro diagnostic use only.

Exercise the normal precautions required for handling all laboratory reagents. This product has been tested and found free from HBsAg and antibodies to HCV and HIV. However, because no test method can offer complete assurance that products derived from human sources will not transmit infectious agents, it is recommended that this product be handled with the same precautions used for patient specimens.

Disposal

This product is made to be used in professional laboratories. Please consult local regulations for a correct waste disposal.

P501: Dispose of contents according to national/international regulations.

References

"The certification of the mass concentration of C-Reactive Protein in human serum - Certified Reference Material ERM-DA474/IFCC", European Community report, 2011

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