APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912	
TEST:	DBIL
APP. CODE:	294
WAVELENGTH (Sec/Pri):	700 - 546
ASSAY:	2 POINT END TIME: 10 POINT: 16 - 31 DILUENT: water
SAMPLE VOL:	NORMAL: 10 DECREASE: 8 INCREASE: 12
	R1 VOLUME: 200 DILUENT: 5 R2 VOLUME: 0 R3 VOLUME: 50 DILUENT: 5 R4 VOLUME: 0
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)
SD LIMIT:	0.250
DUPLICATE LIMIT:	3%
ST. 1 CONC:	0.00
EXPECTED VALUE:	0.00 - 0.20
UNIT:	mg/dl
INSTR. FACTOR (y=ax+b):	a=1 b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 851)	
TEST NAME:	DBIL
SAMPLE:	Volume 10 µl Dilution 0 µl
REAGENTS:	R1 Volume 200 µl Dilution 0 µl R2 Volume 50 µl Dilution 0 µl
WAVELENGTH:	Pri. 540 Sec. 700
METHOD:	END
REACTION SLOPE:	+
MEASURING POINT 1:	First 0 Last 27
MEASURING POINT 2:	First 0 Last 10
REAGENT OD LIMIT:	First L -0.1 First H 0.5 Last L -0.1 Last H 0.5
DYNAMIC RANGE:	L 0.05 H 13
CORRELATION FACTOR:	A 1 B 0
UNIT:	mg/dl
CALIBRATION TYPE:	AB
FORMULA:	Y = AX + B

 Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

ITALIANO rev. 26/09/2016

BILIRUBINA DIRETTA FL	
DD 2H100	4 x 20 + 2 x 10 ml
DD 6U420	6 x 56 + 6 x 14 ml

USO

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emoglobina	≤ 50 mg/dl
lipidi	≤ 500 mg/dl
acido ascorbico	≤ 30 mg/dl

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INTENDED USE

Reagent for quantitative in vitro determination of direct bilirubin in biological fluids.

PRINCIPLE OF THE METHOD

Conjugated (direct) bilirubin reacts with diazotized 2,4-dichloroaniline in acidic solution to produce an intensely coloured red diazo compound (520-560 nm). The intensity of color of this dye in solution is proportional to the concentration of direct bilirubin.

KIT COMPONENTS

For in vitro diagnostic use only.

The components of the kit are stable until expiration date on the label.

Keep away from direct light sources.

BIL D R1 2H100 4 x 20 ml (liquid) white cap
6U420 6 x 56 ml (liquid) white cap

Composition: sodium chloride 0.26 M, EDTA 0.1 mM.

BIL D R2 2H100 2 x 10 ml (liquid) red cap
6U420 6 x 14 ml (liquid) red cap

Composition: EDTA 0.1 mM, diazotized 2,4-dichloroaniline 0.1 mM, hydrochloric acid 0.18 M.

Store all components at 2-8°C.

REAGENT PREPARATION

Use separate reagent ready to use.

Stability: up to expiration date on labels at 2-8°C.

Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-

Caution: keep well refrigerated.

PRECAUTIONS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

SPECIMEN

Serum, plasma.

Specimens should be protected from direct exposure to light. Samples stored at 2-8°C in the dark are stable up to 3 days and 1 month at -20°C.

EXPECTED VALUES

adults: ≤ 0.20 mg/dl (≤ 3.4 µmol/l)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

QUANTINORM CHEMA

with normal or close to normal control values

QUANTIPATH CHEMA

with pathological control values.

If required, a multiparametric, human based calibrator is available:

AUTOCAL H

Please contact Customer Care for further information.

TEST PERFORMANCE

Linearity

the method is linear up to 13 mg/dl.

If the limit value is exceeded, it is suggested to dilute sample 1+9 with distilled water and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)

the limit of detection is 0.039 mg/dl.

Interferences

No interference was observed by the presence of:

hemoglobin	≤ 50 mg/dl
lipids	≤ 500 mg/dl
ascorbic acid	≤ 30 mg/dl

Precision

intra-assay (n=10)	mean (mg/dl)	SD (mg/dl)	CV%
sample 1	0.719	0.003	0.44
sample 2	2.430	0.019	0.78

inter-assay (n=20)	mean (mg/dl)	SD (mg/dl)	CV%
sample 1	0.735	0.039	5.31
sample 2	2.456	0.105	4.26

Methods comparison

a comparison between Chema and a commercially available product gave the following results:

Bilirubin direct Chema = x
Bilirubin direct competitor = y
n = 110

$$y = 0.911x - 0.049 \text{ mg/dl} \quad r^2 = 0.995$$

WASTE DISPOSAL

This product is made to be used in professional laboratories.

P501: Dispose of contents according to national/international regulations.