APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΨΑΜΜΑ HITACHI 911/912	
TEST:	<b>CKMB</b>
APP. CODE:	<b>060</b>
WAVELENGTH (Sec/Pri):	<b>546 - 340</b>
ASSAY:	<b>RATE-A</b> <span style="float: right;">TIME: 10 POINT: 23 - 30 DILUENT: water</span>
SAMPLE VOL:	NORMAL: 10 DECREASE: 7 INCREASE: 12
R1 VOLUME:	<b>200</b> <span style="float: right;">DILUENT: 5</span>
R2 VOLUME:	<b>0</b>
R3 VOLUME:	<b>50</b> <span style="float: right;">DILUENT: 5</span>
R4 VOLUME:	<b>0</b>
ABS LIMIT:	<b>32000 - INC</b>
PROZONE LIMIT:	<b>0 - UPPER</b>
CALIB METHOD:	<b>LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)</b>
SD LIMIT:	<b>0.250</b>
DUPLICATE LIMIT:	<b>3%</b>
ST. 1 CONC:	<b>0.00</b>
EXPECTED VALUE:	<b>0 - 24</b>
UNIT:	<b>U/I</b>
INSTR. FACTOR (y=ax+b):	a=1 b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΨΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 863)	
TEST NAME:	<b>CKMB</b>
SAMPLE:	Volume 10 µl <span style="float: right;">Dilution 0 µl</span>
REAGENTS:	R1 Volume 200 µl <span style="float: right;">Dilution 0 µl</span> R2 Volume 50 µl <span style="float: right;">Dilution 0 µl</span>
WAVELENGTH:	Pri. 340 Sec. 540
METHOD:	<b>RATE</b>
REACTION SLOPE:	<b>+</b>
MEASURING POINT 1:	First 18 <span style="float: right;">Last 26</span>
MEASURING POINT 2:	First <span style="float: right;">Last</span>
REAGENT OD LIMIT:	First L -0.1 <span style="float: right;">First H 1.0</span> Last L -0.1 <span style="float: right;">Last H 1.0</span>
DYNAMIC RANGE:	L 4 <span style="float: right;">H 2000</span>
CORRELATION FACTOR:	A 1 <span style="float: right;">B 0</span>
LINEARITY LIMIT:	<b>15%</b>
UNIT:	<b>U/I</b>
CALIBRATION TYPE:	<b>AB</b>
FORMULA:	<b>Y = AX + B</b>

 Chema Diagnostica  
 Via Campania 2/4  
 60030 Monsano (AN) - ITALY - EU  
 phone +39 0731 605064  
 fax +39 0731 605672  
 e-mail: mail@chema.com  
 website: http://www.chema.com

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**INTENDED USE**  
Reagent for quantitative in vitro determination of creatine kinase MB in biological fluids.

**PRINCIPLE OF THE METHOD**  
CK-MB consists of the subunits CK-M and CK-B. Specific antibodies against CK-M inhibits completely CK-MM activity (main part of the total CK activity) and CK-M subunit of CK-MB. Therefore only CK-B activity is measured, which is half of the CK-MB activity.

**KIT COMPONENTS**  
For in vitro diagnostic use only.  
The components of the kit are stable until expiration date on the label.  
Keep away from direct light sources.

**CK-MB R1** 2H100: 4 x 20 ml (liquid) white cap  
6U140: 2 x 56 ml (liquid) white cap

**CK-MB R2** 2H100: 2 x 10 ml (liquid) red cap  
6U140: 2 x 14 ml (liquid) red cap

Composition in the test: Buffer 100 mM pH 6.70, creatine phosphate 35 mM, glucose 20 mM, N-acetyl-L-cysteine 20 mM, magnesium acetate 10 mM, EDTA 2 mM, ADP 2 mM, NADP 2 mM, AMP 5 mM, Di(adenosine-5')pentaphosphate 10 µM, glucose-6-phosphate-dehydrogenase ≥ 1.5 kU/I, hexokinase ≥ 2.5 kU/I, Anti-CK-M monoclonal antibodies - inhibiting capacity > 2000 U/I.

Store all components at 2-8°C.

**REAGENT PREPARATION**  
Use separate reagent ready to use.  
Stability: up to expiration date on labels at 2-8°C.  
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources.  
Caution: keep well refrigerated.

**PRECAUTIONS**  
Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.  
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

SPECIMEN

Serum is the preferred specimen. Plasma containing heparin, EDTA, citrate, or fluoride may produce unpredictable reaction rates. CK activity in serum is unstable and is rapidly lost during storage. CK is inactivated both by bright daylight and by increasing specimen pH owing to loss of carbon dioxide; accordingly, specimens should be stored in the dark in tightly closed tubes. CK is susceptible to thermal denaturation; the degree of inactivation corresponds to the degree of temperature increase. Therefore, the serum specimen should be chilled to 4°C as rapidly as possible after collection. A slight degree of hemolysis can be tolerated because erythrocytes contain no CK activity. However, moderately or severely hemolyzed specimens are unsatisfactory because enzymes and intermediates liberated from the erythrocytes may affect the lag phase and the side reactions may occur in the assay system.

**EXPECTED VALUES**  
Serum: < 24 U/I ( $< 0.40 \mu\text{kat/l}$ )

Each laboratory should establish appropriate reference intervals related to its population.

**QUALITY CONTROL AND CALIBRATION**  
It is suggested to perform an internal quality control. For this purpose the following control serum is available:

**QUANTINORM CHEMA**  
If required, a multiparametric, human based calibrator is available:  
**AUTOCAL H**

Please contact Customer Care for further information.

**TEST PERFORMANCE**

**Linearity**  
The method is linear up to 2000 U/I.  
If a ΔA/min of 0.250 is exceeded, it is suggested to dilute sample 1+9 with saline solution and to repeat the test, multiplying the result by 10.

**Sensitivity/limit of detection (LOD)**  
The limit of detection is 4 U/I.

**Interferences**  
no interference was observed by the presence of:  
lipids ≤ 1700 mg/dl  
bilirubin ≤ 46 mg/dl  
hemoglobin ≤ 40 mg/dl  
ascorbic acid ≤ 47 mg/dl

**Precision**

intra-assay (n=10)	mean (U/I)	SD (U/I)	CV%
sample 1	46.21	1.01	2.18
sample 2	101.46	1.80	1.77

inter-assay (n=20)	mean (U/I)	SD (U/I)	CV%
sample 1	46.35	1.31	2.82
sample 2	101.64	1.03	1.01

**Methods comparison**  
A comparison between Chema CK-MB FL and a commercially available product gave the following results:

CK MB Chema = y  
CK MB competitor = x  
n = 82

$y = 1.00 x + 0.46 \text{ U/I}$   $r^2 = 0.999$

**WASTE DISPOSAL**  
This product is made to be used in professional laboratories.  
P501: Dispose of contents according to national/international regulations.

CK-MB FL IFCC/DGKC	
MB 2H100	4 x 20 + 2 x 10 ml
MB 6U140	2 x 56 + 2 x 14 ml

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