APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912	
TEST:	CKNAC
APP. CODE:	336
WAVELENGTH (Sec/Pri):	700 - 340
ASSAY:	RATE-A <i>TIME: 10</i> <i>POINT: 24 - 31</i> <i>DILUENT: water</i>
SAMPLE VOL:	NORMAL: 10 DECREASE: 7 INCREASE: 12
R1 VOLUME:	200 <i>DILUENT: 5</i>
R2 VOLUME:	0
R3 VOLUME:	50 <i>DILUENT: 5</i>
R4 VOLUME:	0
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)
SD LIMIT:	0.250
DUPLICATE LIMIT:	3%
ST. 1 CONC:	0.00
EXPECTED VALUE:	0 - 200
UNIT:	U/I
INSTR. FACTOR (y=ax+b):	a=1 b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 864)	
TEST NAME:	CKNAC
SAMPLE:	Volume 10 µl Dilution 0 µl
REAGENTS:	R1 Volume 200 µl Dilution 0 µl R2 Volume 50 µl Dilution 0 µl
WAVELENGTH:	Pri. 340 Sec. 700
METHOD:	RATE
REACTION SLOPE:	+
MEASURING POINT 1:	First 18 Last 26
MEASURING POINT 2:	First Last
REAGENT OD LIMIT:	First L -0.1 First H 1.0 Last L -0.1 Last H 1.0
DYNAMIC RANGE:	L 1 H 2000
CORRELATION FACTOR:	A 1 B 0
LINEARITY LIMIT:	15%
UNIT:	U/I
CALIBRATION TYPE:	AB
FORMULA:	Y = AX + B

Chema Diagnostica
 Via Campania 2/4
 60030 Monsano (AN) - ITALY - EU
 phone +39 0731 605064
 fax +39 0731 605672
 e-mail: mail@chema.com
 website: http://www.chema.com

CK-NAC FL IFCC/DGKC	
CK 2H100	4 x 20 + 2 x 10 ml
CK 6U280	4 x 56 + 4 x 14 ml

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INTENDEDO USE
 Reagent for quantitative in vitro determination of creatine kinase in biological fluids.

PRINCIPLE OF THE METHOD
 Creatine kinase (EC 2.7.3.2; adenosine triphosphate: creatine N-phosphotransferase; CK) catalyzes the conversion of creatine phosphate and ADP to creatine and ATP. ATP and glucose are converted to ADP and glucose-6-phosphate by hexokinase. Glucose-6-phosphate dehydrogenase oxidizes glucose-6-phosphate to 6-phosphogluconate, reducing NADP to NADPH. The rate of conversion of NADP/NADPH, monitored at 340 nm, is proportional to CK activity. N-acetyl cysteine (NAC) is added as an activator of CK.

KIT COMPONENTS
 For in vitro diagnostic use only.
 The components of the kit are stable until expiration date on the label.
 Keep away from direct light sources.

CK-NAC R1 2H100: 4 x 20 ml (liquid) white cap
 6U280: 4 x 56 ml (liquid) white cap

CK-NAC R2 2H100: 2 x 10 ml (liquid) red cap
 6U280: 4 x 14 ml (liquid) red cap

Composition in the test: imidazole buffer 29 mM pH 6.50, creatine phosphate 30 mM, glucose 20 mM, N-acetyl-L-cysteine 20 mM, magnesium acetate 10 mM, EDTA 2 mM, ADP 2 mM, NADP 2 mM, AMP 5 mM, D(adenosine-5') pentaphosphate 12 µM, glucose-6-phosphate-dehydrogenase ≥ 3 kU/l, hexokinase ≥ 3 kU/l.

Store all components at 2-8°C.

REAGENT PREPARATION
 Use separate reagent ready to use.
 Stability: up to expiration date on labels at 2-8°C.
 Stability since first opening of vials: preferably within 60 days at 2-8°C - away from light sources.
 Caution: keep well refrigerated.

PRECAUTIONS
 Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.
 Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

SPECIMEN
 Serum is the preferred specimen. Plasma containing heparin, EDTA, citrate, or fluoride may produce unpredictable reaction rates. CK activity in serum is unstable and is rapidly lost during storage. CK is inactivated both by bright daylight and by increasing specimen pH owing to loss of carbon dioxide; accordingly, specimens should be stored in the dark in tightly closed tubes. CK is susceptible to thermal denaturation; the degree of inactivation corresponds to the degree of temperature increase. Therefore, the serum specimen should be chilled to 4°C as rapidly as possible after collection. A slight degree of hemolysis can be tolerated because erythrocytes contain no CK activity. However, moderately or severely hemolyzed specimens are unsatisfactory because enzymes and intermediates liberated from the erythrocytes may affect the lag phase and the side reactions occurring in the assay system.

EXPECTED VALUES		
Men	24 - 204 U/I	(0.39 - 3.40 µkat/l)
Women	24 - 173 U/I	(0.39 - 2.90 µkat/l)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION
 It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:
QUANTINORM CHEMA
 with normal or close to normal control values
QUANTIPATH CHEMA
 with pathological control values.
 If required, a multiparametric, human based calibrator is available:
AUTOCAL H

Please contact Customer Care for further informations.

TEST PERFORMANCE
Linearity
 the method is linear up to 2000 U/I.
 If a ΔA/min of 0.250 is exceeded, it is suggested to dilute sample 1+9 with saline solution and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)
 the limit of detection is 1.6 U/I.

Interferences
 no interference was observed by the presence of:
 hemoglobin ≤ 400 mg/dl
 bilirubin ≤ 40 mg/dl
 lipids ≤ 660 mg/dl

Precision	mean (U/I)	SD (U/I)	CV%
intra-assay (n=10) sample 1	148.21	0.94	0.64
sample 2	464.75	3.98	0.86
inter-assay (n=20) sample 1	148.35	1.33	0.90
sample 2	461.34	4.62	1.00

Methods comparison
 a comparison between Chema CK-NAC FL and a commercially available product gave the following results:

CK NAC Chema = x
 CK-NAC competitor = y
 n = 100
 y = 1.04x - 3.10 U/I r² = 0.9985

WASTE DISPOSAL
 This product is made to be used in professional laboratories.
 P501: Dispose of contents according to national/international regulations.

CK-NAC FL IFCC/DGKC	
CK 2H100	4 x 20 + 2 x 10 ml
CK 6U280	4 x 56 + 4 x 14 ml

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