APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ <b>HITACHI 911/912</b>	
TEST:	<b>CREA</b>
APP. CODE:	<b>339</b>
WAVELENGTH (Sec/Pri):	<b>570 - 505</b>
ASSAY:	<b>2 POINT RATE</b> TIME: <b>10</b> POINT: <b>19-23</b>
SAMPLE VOL:	NORMAL: <b>25</b> DECREASE: <b>20</b> INCREASE: <b>30</b>
	R1 VOLUME: <b>125</b> R2 VOLUME: <b>0</b> R3 VOLUME: <b>125</b> R4 VOLUME: <b>0</b>
ABS LIMIT:	<b>32000 - INC</b>
PROZONE LIMIT:	<b>0 - UPPER</b>
CALIB METHOD:	<b>LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)</b>
SD LIMIT:	<b>0.250</b>
DUPLICATE LIMIT:	<b>3%</b>
ST. 1 CONC:	<b>0.00</b>
EXPECTED VALUE:	<b>0.6 - 1.2</b>
UNIT:	<b>mg/dl</b>
INSTR. FACTOR (y=ax+b):	a=1 b= <b>-0.3</b>

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ <b>OLYMPUS AU 400/480/600/640/680/2700 (Test code 869)</b>	
TEST NAME:	<b>CREA</b>
SAMPLE:	Volume <b>25</b> µl Dilution <b>0</b> µl
REAGENTS:	R1 Volume <b>125</b> µl Dilution <b>0</b> µl R2 Volume <b>125</b> µl Dilution <b>0</b> µl
WAVELENGTH:	Pri. <b>520</b> Sec. <b>570</b>
METHOD:	<b>FIXED</b>
REACTION SLOPE:	<b>+</b>
MEASURING POINT 1:	First <b>14</b> Last <b>18</b>
MEASURING POINT 2:	First Last
REAGENT OD LIMIT:	First L <b>-0.1</b> First H <b>0.5</b> Last L <b>-0.1</b> Last H <b>0.5</b>
DYNAMIC RANGE:	L <b>0.2</b> H <b>20</b>
CORRELATION FACTOR:	A <b>1</b> B <b>-0.3</b>
UNIT:	<b>mg/dl</b>
CALIBRATION TYPE:	<b>AB</b>
FORMULA:	<b>Y = AX + B</b>

 Chema Diagnostica  
 Via Campania 2/4  
 60030 Monsano (AN) - ITALY - EU  
 phone +39 0731 605064  
 fax +39 0731 605672  
 e-mail: mail@chema.com  
 website: http://www.chema.com

ITALIANO rev. 03/1/2017

CREATININA	
CR 2H400	4 x 50 + 4 x 50 ml
CR 6U448	4 x 56 + 4 x 56 ml

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**INTENDED USE**  
Reagent for quantitative in vitro determination of creatinine in biological fluids.

**PRINCIPLE OF THE METHOD**  
Creatinine reacts with picric acid in alkaline environment to form a color complex. Developing of this red color may be followed photometrically at 500-520 nm. The association on surfactant and sodium tetraborate keeps interferences at minimum.

**KIT COMPONENTS**

**For in vitro diagnostic use only.**  
The components of the kit are stable until expiration date on the label.  
Keep away from direct light sources.

**CREA R1** 2H400 4 x 50 ml (liquid) white cap  
6U448 4 x 56 ml (liquid) white cap

**CREA R2** 2H400 4 x 50 ml (liquid) red cap  
6U448 4 x 56 ml (liquid) red cap

Composition in the test: picric acid 14 mM, NaOH 0.18 M, sodium tetraborate 10 mM, surfactant.  
Store all components at 15-25°C.

**REAGENT PREPARATION**

Use separate reagent ready to use.  
Stability of working reagent: preferably within 30 days at 15-25°C, well capped and away from light sources.  
Stability of unmixed reagents: up to expiration date on labels at 15-25°C;  
Stability since first opening of vials of unmixed reagents: preferably within 60 days at 15-25°C.

PRECAUTIONS

**CREA R1: Warning.** Causes serious eye irritation (H319). Causes skin irritation (H315). Wear protective gloves. Eye protection (P280). IF ON SKIN: Wash with plenty of water (P302+P352). IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). If eye irritation persists: get medical advice (P337+P313).

**CREA R2:** It is not classified as hazardous.

**SPECIMEN**

Serum, plasma, Urine.  
Creatinine is stable 24 hours at 2-8°C. Freeze samples for prolonged storage.  
Dilute urine sample 1:100 with deionized water. It could be convenient a slight acidification of urine with HCl.

**EXPECTED VALUES**

Serum/plasma samples:  
Men: 0.7 - 1.2 mg/dl (62 - 105 µmol/l)  
Women: 0.6 - 1.1 mg/dl (53 - 97 µmol/l)

24h urine:  
Men: 1000 - 2000 mg/24h (8.85 - 17.70 mmol/24h)  
Women: 800 - 1800 mg/24h (7.08 - 15.93 mmol/24h)

Each laboratory should establish appropriate reference intervals related to its population.

**QUALITY CONTROL AND CALIBRATION**

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

**QUANTINORM CHEMA**  
with normal or close to normal control values  
**QUANTIPATH CHEMA**  
with pathological control values.  
If required, a multiparametric, human based calibrator is available:  
**AUTOCAL H**

Please contact Customer Care for further information.

**TEST PERFORMANCE**

**Linearity**  
the method is linear up to 20 mg/dl.  
If the value is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

**Sensitivity/limit of detection (LOD)**  
the limit of detection is 0.2 mg/dl.

**Interferences**  
no interference was observed by the presence of:  
hemoglobin ≤ 500 mg/dl  
lipids ≤ 1250 mg/dl  
Bilirubins give interference at low levels.

Precision	mean (mg/dl)	SD (mg/dl)	CV%
intra-assay (n=10)			
sample 1	1.25	0.03	2.60
sample 2	3.87	0.07	1.90
inter-assay (n=20)			
sample 1	1.31	0.04	2.90
sample 2	3.80	0.14	3.80

**Methods comparison**  
a comparison between Chema and a commercially available product gave the following results:

Creatinine Chema = x  
Creatinine competitor = y  
n = 104  
 $y = 0.982x - 0.081 \text{ mg/dl} \quad r^2 = 0.94$

**WASTE DISPOSAL**

This product is made to be used in professional laboratories.  
P501: Dispose of contents according to national/international regulations.