**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ПРОГРАММА**  
**HITACHI 911/912**

TEST: ALP DGKC

APP. CODE: 342

WAVELENGTH (Sec/Pri): 700 - 415

ASSAY: RATE-A

TIME: 10

POINT: 20 - 30

DILUENT: water

SAMPLE VOL: NORMAL: 7

DECREASE: 5

INCREASE: 10

R1 VOLUME: 280 DILUENT: 5

R2 VOLUME: 0

R3 VOLUME: 70 DILUENT: 5

R4 VOLUME: 0

ABS LIMIT: 32000 - INC

PROZONE LIMIT: 0 - UPPER

CALIB METHOD: LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)

SD LIMIT: 0.250

DUPLICATE LIMIT: 3%

ST. 1 CONC: 0.0

EXPECTED VALUE: 0 - 270

UNIT: U/I

INSTR. FACTOR (y=ax+b): a=1 b= 0

**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ПРОГРАММА**  
**OLYMPUS AU 400/480/600/640/680/2700 (Test code 872)**

TEST NAME: ALP DGKC

SAMPLE: Volume 7  $\mu$ l Dilution 0  $\mu$ l

REAGENTS: R1 Volume 280  $\mu$ l Dilution 0  $\mu$ l  
R2 Volume 70  $\mu$ l Dilution 0  $\mu$ l

WAVELENGHT: Pri. 410 Sec. 700

METHOD: RATE

REACTION SLOPE: +

MEASURING POINT 1: First 14 Last 27

MEASURING POINT 2: First Last

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**INTENDED USE**

Reagent for quantitative in vitro determination of alkaline phosphatase in biological fluids.

**PRINCIPLE OF THE METHOD**

The enzyme alkaline phosphatase (EC 3.1.3.1, orthophosphoric-monoester phosphohydrolase) hydrolyzes the 4-NPP to release 4-nitrophenol, under alkaline conditions. The 4-nitrophenol formed is detected spectrophotometrically at 405 nm to give a measurement of alkaline phosphatase activity in the sample.

The present method has been made according to DGKC.

**KIT COMPONENTS**

**For in vitro diagnostic use only.**  
The components of the kit are stable until expiration date on the label.

Keep away from direct light sources.

**ALP DGKC R1** 2H100: 4 x 20 ml (liquid) white cap  
2H400: 8 x 40 ml (liquid) white cap  
6U420: 6 x 56 ml (liquid) white cap

**ALP DGKC R2** 2H100: 2 x 10 ml (liquid) red cap  
2H400: 4 x 20 ml (liquid) red cap  
6U420: 6 x 14 ml (liquid) red cap

Composition in the test: DEA buffer pH 9.8 1 M, MgCl<sub>2</sub> 0.5 mM, 4-Nitrophenylphosphate 10 mM.

Store all components at 2-8°C.

**REAGENT PREPARATION**

Use separate reagent ready to use.

Stability: up to expiration date on labels at 2-8°C.

Stability since first opening of vials: preferably within 60 days at 2-8°C - away from light sources.

Caution: keep well refrigerated.

**PRECAUTIONS**

**ALP DGKC R1: Danger.** May cause damage to organs through prolonged or repeated exposure (H373). Causes serious eye damage (H318). Causes skin irritation (H315). Wear protective gloves. Eye protection (P280). Do not breath vapours (P260). IF ON SKIN: Wash with plenty of water (P305+P352). IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). Get medical advice if you feel unwell (P314).

**ALP DGKC R2:** It is not classified as hazardous.

**SPECIMEN**

Serum, plasma (heparinate only).  
Sera kept at room temperatures usually show a slight but real increase in activity, which varies from 1% over a 6-h period to 3 to 6% over a 1 to 4 days period. Even in sera stored at refrigerator temperature, activity increases slowly. In frozen sera, activity decreases but slowly recovers after thawing the serum.  
A similar enhancement of activity, but of greater magnitude, occurs with reconstituted lyophilized preparations, such as those available as control sera or calibrators. In reconstituted material the increases with storage at 4 and 20°C are about 10 and 30%, respectively. Enhancement of activity continues for several days, but at a decreasing rate. The cause of this phenomenon is not known but may be attributed to renaturation of partially denatured enzyme or to dissociation, on warming, of a phosphate-lipoprotein complex or a multimer of the enzyme that was formed in the freeze-drying process.

**EXPECTED VALUES**

Men:	< 270 U/l	(< 4.50 $\mu$ kat/l)
Women:	< 240 U/l	(< 4.00 $\mu$ kat/l)

**QUALITY CONTROL AND CALIBRATION**

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

**QUANTINORM CHEMA**

with normal or close to normal control values.

**QUANTIPATH CHEMA**

with pathological control values.  
If required, a multiparametric, human based calibrator is available:

**AUTOCAL H**

Please contact Customer Care for further informations.

**TEST PERFORMANCE**

**Linearity**

the method is linear up to 3000 U/l.  
If a AA/min of 0.500 is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

**Sensitivity/limit of detection (LOD)**

the limit of detection is 7 U/l.

**Interferences**

no interference was observed by the presence of:  
hemoglobin ≤ 400 mg/dl  
bilirubin ≤ 27 mg/dl  
lipids ≤ 1000 mg/dl

**Precision**

intra-assay (n=10)	mean (U/l)	SD (U/l)	CV%
sample 1	175.70	0.95	0.50
sample 2	426.70	2.41	0.60

inter-assay (n=20)	mean (U/l)	SD (U/l)	CV%
sample 1	167.26	3.99	2.40
sample 2	408.28	8.61	2.10

**Methods comparison**

a comparison between Chema and a commercially available product gave the following results:

ALP Chema = x

ALP competitor = y

n = 112

y = 0.96x - 2.17 U/l  $r^2 = 0.999$

**WASTE DISPOSAL**

This product is made to be used in professional laboratories.

P501: Dispose of contents according to national/international regulations.

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