APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912	
TEST:	ALP IFCC
APP. CODE:	343
WAVELENGTH (Sec/Pri):	700 - 415
ASSAY:	RATE-A <i>TIME: 10</i> <i>POINT: 20 - 30</i> <i>DILUENT: water</i>
SAMPLE VOL:	NORMAL: 7 DECREASE: 5 INCREASE: 10
R1 VOLUME: 280	<i>DILUENT: 5</i>
R2 VOLUME: 0	
R3 VOLUME: 70	<i>DILUENT: 5</i>
R4 VOLUME: 0	
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)
SD LIMIT:	0.250
DUPLICATE LIMIT:	3%
ST. 1 CONC:	0.0
EXPECTED VALUE:	30 - 130
UNIT:	U/l
INSTR. FACTOR (y=ax+b):	a=1 b= 0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 873)	
TEST NAME:	ALP IFCC
SAMPLE:	Volume 7 µl Dilution 0 µl
REAGENTS:	R1 Volume 280 µl Dilution 0 µl R2 Volume 70 µl Dilution 0 µl
WAVELENGTH:	Pri. 410 Sec. 700
METHOD:	RATE
REACTION SLOPE:	+
MEASURING POINT 1:	First 14 Last 27
MEASURING POINT 2:	First Last
REAGENT OD LIMIT:	First L -0.1 First H 1.0 Last L -0.1 Last H 1.0
DYNAMIC RANGE:	L 5 H 3000
CORRELATION FACTOR:	A 1 B 0
LINEARITY LIMIT:	15%
UNIT:	U/l
CALIBRATION TYPE:	AB
FORMULA:	Y = AX + B

 Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

FOSFATASI ALCALINA FL IFCC	
AF 2H100	4 x 20 + 2 x 10 ml
AF 2H400	8 x 40 + 4 x 20 ml
AF 6U420	6 x 56 + 6 x 14 ml

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INTENDED USE
Reagent for quantitative in vitro determination of alkaline phosphatase in biological fluids.

PRINCIPLE OF THE METHOD
The enzyme alkaline phosphatase (EC 3.1.3.1., orthophosphoric-monoester phosphohydrolase) hydrolyzes the 4-NPP to release 4-nitrophenol, under alkaline conditions. The 4-nitrophenol formed is detected spectrophotometrically at 405 nm to give a measurement of alkaline phosphatase activity in the sample.
The present method has been made according to IFCC.

KIT COMPONENTS
For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label.
Keep away from direct light sources.

ALP IFCC R1 2H100: 4 x 20 ml (liquid) white cap
2H400: 8 x 40 ml (liquid) white cap
6U420: 6 x 56 ml (liquid) white cap

ALP IFCC R2 2H100: 2 x 10 ml (liquid) red cap
2H400: 4 x 20 ml (liquid) red cap
6U420: 6 x 14 ml (liquid) red cap

Composition in the test: 2-amino-2-methyl-1-propanol buffer 0.35 M pH 10.40 (30°C), magnesium acetate 2 mM, zinc sulfate 1 mM, HEDTA 2 mM, 4-NPP 16 mM.

Store all components at 2-8°C.

REAGENT PREPARATION
Use separate reagent ready to use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-
Caution: keep well refrigerated.

PRECAUTIONS
ALP IFCC R1: Warning. Causes serious eye irritation (H319). Causes skin irritation (H315). Wear protective gloves. Eye protection (P280). IF ON SKIN: wash with plenty of water (P302+P352). IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). If eye irritation persists: get medical advice (P337+P313).

ALP IFCC R2: It is not classified as hazardous.

SPECIMEN
Serum, plasma (heparinized only).
Sera kept at room temperatures usually show a slight but real increase in activity, which varies from 1% over a 6-h period to 3 to 6% over a 1 to 4 days period. Even in sera stored at refrigerator temperature, activity increases slowly. In frozen sera, activity decreases but slowly recovers after thawing the serum.
A similar enhancement of activity, but of greater magnitude, occurs with reconstituted lyophilized preparations, such as those available as control sera or calibrators. In reconstituted material the increases with storage at 4 and 20°C are about 10 and 30%, respectively. Enhancement of activity continues for several days, but at a decreasing rate. The cause of this phenomenon is not known but may be attributed to renaturation of partially denatured enzyme or to dissociation, on warming, of a phosphate-lipoprotein complex or a multimer of the enzyme that was formed in the freeze-drying process.

EXPECTED VALUES

Men:	40 - 129 U/l	(0.67 - 2.15 µkat/l)
Women:	35 - 104 U/l	(0.58 - 1.74 µkat/l)

QUALITY CONTROL AND CALIBRATION
It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:
QUANTINORM CHEMA with normal or close to normal control values
QUANTIPATH CHEMA with pathological control values.
If required, a multiparametric, human based calibrator is available:
AUTOCAL H

Please contact Customer Care for further informations.

TEST PERFORMANCE
Linearity
the method is linear up to 3000 U/l.
If a ΔA/min of 0.500 is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)
the limit of detection is 5.2 U/l.

Interferences
no interference was observed by the presence of:
hemoglobin ≤ 400 mg/dl
bilirubin ≤ 40 mg/dl
lipids ≤ 900 mg/dl

Precision

intra-assay (n=10)	mean (U/l)	SD (U/l)	CV%
sample 1	84.40	2.41	2.86
sample 2	222.40	5.74	2.58

inter-assay (n=20)	mean (U/l)	SD (U/l)	CV%
sample 1	86.66	2.66	3.07
sample 2	210.39	6.08	2.89

Methods comparison
a comparison between Chema and a commercially available product gave the following results:

ALP Chema = x	ALP competitor = y	n = 150
y = 1.03x - 2.57 U/l	r ² = 0.998	

WASTE DISPOSAL
This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.

REFERENCES
Clin. Chim. Acta, (1983) 339F - 367F
Tietz Textbook of Clinical Chemistry, Second Edition, Burtis-Ashwood (1994).

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