APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ <b>HITACHI 911/912</b>	
TEST:	<b>GOT</b>
APP. CODE:	<b>327</b>
WAVELENGTH (Sec/Pri):	<b>700 - 340</b>
ASSAY:	<b>RATE-A</b> <i>TIME: 10 POINT: 20-30 DILUENT: water</i>
SAMPLE VOL:	NORMAL: <b>25</b> DECREASE: <b>20</b> INCREASE: <b>30</b>
R1 VOLUME:	<b>200</b> <i>DILUENT: 5</i>
R2 VOLUME:	<b>0</b>
R3 VOLUME:	<b>50</b> <i>DILUENT: 5</i>
R4 VOLUME:	<b>0</b>
ABS LIMIT:	<b>6500 - DEC</b>
PROZONE LIMIT:	<b>0 - UPPER</b>
CALIB METHOD:	<b>LINEAR</b> (POINT: 2 - SPAN: 2 - WEIGHT: 0)
SD LIMIT:	<b>0.250</b>
DUPLICATE LIMIT:	<b>3%</b>
ST. 1 CONC:	<b>0.0</b>
EXPECTED VALUE:	<b>0.0 - 35.0</b>
UNIT:	<b>U/I</b>
INSTR. FACTOR (y=ax+b):	a=1 b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ  
**OLYMPUS AU 400/480/600/640/680/2700 (Test code 853)**

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**INTENDED USE**  
Reagent for quantitative in vitro determination of GOT in biological fluids.

**PRINCIPLE OF THE METHOD**  
The enzyme aspartate aminotransferase (EC 2.6.1.1; L-Aspartate:2-Oxoglutarate Aminotransferase, AST or AspAT; Glutamate Oxaloacetate Transaminase, GOT) catalyzes the transaminase reaction between L-Aspartate and 2-Oxoglutarate. The 2-Oxalacetate formed, is reduced to malate in the presence of MDH. As the reactions proceed, NADH is oxidized to NAD. The disappearance of NADH per unit time is followed by measuring the decrease in absorbance at 340 nm. The present method has been made according to IFCC (2002).

**KIT COMPONENTS**  
**For in vitro diagnostic use only.**  
The components of the kit are stable until expiration date on the label.  
Keep away from direct light sources.

**GOT R1** 2H400: 8 x 40 ml (liquid) white cap  
6U420: 6 x 56 ml (liquid) white cap

**GOT R2** 2H400: 4 x 20 ml (liquid) red cap  
6U420: 6 x 14 ml (liquid) red cap

Composition in the test: Tris buffer 80 mM pH 7.65, L-aspartate 240 mM, 2-Oxoglutarate 12 mM, NADH 0.18 mM, MDH ≥ 600 U/I, LDH ≥ 900 U/I.

Store all components at 2-8°C.

**REAGENT PREPARATION**  
Use separate reagent ready to use.  
Stability: up to expiration date on labels at 2-8°C.  
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-  
Caution: keep well refrigerated.

**PRECAUTIONS**  
**GOT R1: Warning.** Causes serious eye irritation (H319).  
 Causes skin irritation (H315). Wear protective gloves. Eye protection (P280). IF ON SKIN: Wash with plenty of water (P302+P352). IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). If eye irritation persists: get medical advice (P337+P313).

**GOT R2:** It is not classified as hazardous.

**SPECIMEN**

Serum, plasma.  
Specimens should be protected from direct exposure to light. Samples stored at 2-8°C in the dark are stable up to 3 days and 1 month at -20°C.

**EXPECTED VALUES**  
Men: < 35 U/I (< 0.58 µkat/l)  
Women: < 31 U/I (< 0.52 µkat/l)

Each laboratory should establish appropriate reference intervals related to its population.

**QUALITY CONTROL AND CALIBRATION**  
It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:  
**QUANTINORM CHEMA**  
with normal or close to normal control values  
**QUANTIPATH CHEMA**  
with pathological control values.  
If required, a multiparametric, human based calibrator is available:  
**AUTOCAL H**

Please contact Customer Care for further information.

**TEST PERFORMANCE**  
**Linearity**  
the method is linear up to 440 U/I.  
If a ΔA/min of 0.200 is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

**Sensitivity/limit of detection (LOD)**  
the limit of detection is 0.463 U/I.

**Interferences**  
no interference was observed by the presence of:  
hemoglobin ≤ 250 mg/dl  
bilirubin ≤ 45 mg/dl  
lipids ≤ 500 mg/dl

**Precision**  
intra-assay (n=10) mean (U/I) SD (U/I) CV%  
sample 1 46.19 0.31 0.67  
sample 2 137.25 0.92 0.67

inter-assay (n=20) mean (U/I) SD (U/I) CV%  
sample 1 46.18 2.04 4.41  
sample 2 137.76 6.30 4.57

**Methods comparison**  
a comparison between Chema and a commercially available product gave the following results:

GOT Chema = x  
GOT competitor = y  
n = 83  
y = 1.003x - 0.560 U/I r<sup>2</sup> = 0.990

**WASTE DISPOSAL**  
This product is made to be used in professional laboratories.  
P501: Dispose of contents according to national/international regulations.

GOT/AST FL IFCC		
GO 2H400	8 x 40 + 4 x 20 ml	
GO 6U420	6 x 56 + 6 x 14 ml	

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