APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912			
TEST:	GPT		
APP. CODE:	328		
WAVELENGTH (Sec/Pri):	700 - 340		
ASSAY:	RATE-A	<i>TIME: 10</i>	<i>POINT: 20-30</i>
		<i>DILUENT: water</i>	
SAMPLE VOL:	NORMAL: 25 DECREASE: 20 INCREASE: 30		
	R1 VOLUME: 200 R2 VOLUME: 0 R3 VOLUME: 50 R4 VOLUME: 0	<i>DILUENT: 5</i>	<i>DILUENT: 5</i>
ABS LIMIT:	6500 - DEC		
PROZONE LIMIT:	0 - UPPER		
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)		
SD LIMIT:	0.250		
DUPLICATE LIMIT:	3%		
ST. 1 CONC:	0.0		
EXPECTED VALUE:	0.0 - 45.0		
UNIT:	U/l		
INSTR. FACTOR (y=ax+b):	a=1 b=0		

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 854)			
TEST NAME:	GPT		
SAMPLE:	Volume 25 μ l	Dilution 0 μ l	
REAGENTS:	R1 Volume 200 μ l R2 Volume 50 μ l	Dilution 0 μ l Dilution 0 μ l	
WAVELENGTH:	Pri. 340 Sec. 700		
METHOD:	RATE		
REACTION SLOPE:	-		
MEASURING POINT 1:	First 14	Last 23	
MEASURING POINT 2:	First	Last	
REAGENT OD LIMIT:	First L 0.8 Last L 0.8	First H 2.0 Last H 2.0	
DYNAMIC RANGE:	L 0.2	H 440	
CORRELATION FACTOR:	A 1	B 0	
LINEARITY LIMIT:	15%		
UNIT:	U/l		
CALIBRATION TYPE:	AB		
FORMULA:	Y = AX + B		

 Chema Diagnostica
 Via Campania 2/4
 60030 Monsano (AN) - ITALY - EU
 phone +39 0731 605064
 fax +39 0731 605672
 e-mail: mail@chema.com
 website: http://www.chema.com

ITALIANO rev. 26/09/2016

GPT/ALT FL IFCC		
GP 2H400	8 x 40 + 4 x 20 ml	
GP 6U420	6 x 56 + 6 x 14 ml	

USO

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INTENDED USE

Reagent for quantitative in vitro determination of GPT in biological fluids.

PRINCIPLE OF THE METHOD

The enzyme alanine aminotransferase (EC 2.6.1.2; L-Alanine:2-Oxoglutarate Aminotransferase, ALT or A1aAT; Glutamate Pyruvate Transaminase, GPT) catalyzes the transaminase reaction between L-Alanine and 2-Oxoglutarate. The pyruvate formed, is reduced to lactate in the presence of LDH. As the reactions proceed, NADH is oxidized to NAD. The disappearance of NADH per unit time is followed by measuring the decrease in absorbance at 340 nm.

The present method has been made according to IFCC (2002).

KIT COMPONENTS

For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label.
Keep away from direct light sources.

GPT R1 2H400: 8 x 40 ml (liquid) white cap
6U420: 6 x 56 ml (liquid) white cap

GPT R2 2H400: 4 x 20 ml (liquid) red cap
6U420: 6 x 14 ml (liquid) red cap

Composition in the test: Tris buffer 100 mM pH 7.15, L-Alanine 500 mM, 2-Oxoglutarate 15 mM, NADH 0.18 mM, LDH \geq 1700 U/l.

Store all components at 2-8°C.

REAGENT PREPARATION

Use separate reagent ready to use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: preferably within 60 days at 2-8°C away from light sources.
Caution: keep well refrigerated.

PRECAUTIONS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

SPECIMEN

Serum (preferred). Plasma is not recommended.
Collect blood with a minimum of venous stasis.
GPT is stable up to 4 days at 2-8°C or 1 month at -20°C.

EXPECTED VALUES

Men: < 45 U/l (< 0.74 μ kat/l)
Women: < 34 U/l (< 0.56 μ kat/l)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

QUANTINORM CHEMA

with normal or close to normal control values

QUANTIPATH CHEMA

with pathological control values.

If required, a multiparametric, human based calibrator is available:

AUTOCAL H

Please contact Customer Care for further information.

TEST PERFORMANCE

Linearity

the method is linear up to 440 U/l.
If a Δ A/min of 0.200 is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)

the limit of detection is 0.169 U/l.

Interferences

no interference was observed by the presence of:

hemoglobin	\leq 500 mg/dl
bilirubin	\leq 45 mg/dl
lipids	\leq 450 mg/dl

Precision

intra-assay (n=10)	mean (U/l)	SD (U/l)	CV%
sample 1	49.29	0.35	0.71
sample 2	132.15	0.57	0.43

inter-assay (n=20)	mean (U/l)	SD (U/l)	CV%
sample 1	49.31	1.66	3.37
sample 2	132.85	4.28	3.22

Methods comparison

a comparison between Chema and a commercially available product gave the following results:

GPT Chema = x
GPT competitor = y
n = 126
y = 0.992x - 0.299 $r^2 = 0.999$

WASTE DISPOSAL

This product is made to be used in professional laboratories.

P501: Dispose of contents according to national/international regulations.

FRANÇAIS	rev. 26/09/2016
GPT/ALT FL IFCC	
GP 2H400	8 x 40 + 4 x 20 ml
GP 6U420	6 x 56 + 6 x 14 ml

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