APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ <b>HITACHI 911/912</b>	
TEST:	<b>GL-UV</b>
APP. CODE:	<b>346</b>
WAVELENGTH (Sec/Pri):	<b>700 - 505</b>
ASSAY:	<b>2 POINT END</b> <i>TIME: 10</i> <i>POINT: 16-31</i> <i>DILUENT: Water</i>
SAMPLE VOL:	NORMAL: <b>3</b> DECREASE: <b>2</b> INCREASE: <b>5</b>
R1 VOLUME:	<b>240</b>
R2 VOLUME:	<b>0</b>
R3 VOLUME:	<b>60</b> <i>DILUENT: 5</i>
R4 VOLUME:	<b>0</b>
ABS LIMIT:	<b>32000 - INC</b>
PROZONE LIMIT:	<b>0 - UPPER</b>
CALIB METHOD:	<b>LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)</b>
SD LIMIT:	<b>0.250</b>
DUPLICATE LIMIT:	<b>3%</b>
ST. 1 CONC:	<b>0.0</b>
EXPECTED VALUE:	<b>70 - 110</b>
UNIT:	<b>mg/dl</b>
INSTR. FACTOR (y=ax+b):	a=1    b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ <b>OLYMPUS AU 400/480/600/640/680/2700 (Test code 877)</b>	
TEST NAME:	<b>GL-UV</b>
SAMPLE:	Volume <b>3 µl</b> Dilution <b>0 µl</b>
REAGENTS:	R1 Volume <b>240 µl</b> Dilution <b>0 µl</b> R2 Volume <b>60 µl</b> Dilution <b>0 µl</b>
WAVELENGTH:	Pri. <b>340</b> Sec. <b>700</b>
METHOD:	<b>END</b>
REACTION SLOPE:	<b>+</b>
MEASURING POINT 1:	First <b>0</b> Last <b>27</b>
MEASURING POINT 2:	First <b>0</b> Last <b>10</b>
REAGENT OD LIMIT:	First L <b>-0.1</b> First H <b>0.5</b> Last L <b>-0.1</b> Last H <b>0.5</b>
DYNAMIC RANGE:	L <b>1</b> H <b>700</b>
CORRELATION FACTOR:	A <b>1</b> B <b>0</b>
UNIT:	<b>mg/dl</b>
CALIBRATION TYPE:	<b>AB</b>
FORMULA:	<b>Y = AX + B</b>

 Chema Diagnostica  
 Via Campania 2/4  
 60030 Monsano (AN) - ITALY - EU  
 phone +39 0731 605064  
 fax +39 0731 605672  
 e-mail: mail@chema.com  
 website: http://www.chema.com

GLUCOSIO UV FL	
GL 2H251	4 x 50 + 2 x 25 ml
GL 6U421	6 x 56 + 6 x 14 ml

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**INTENDED USE**  
Reagent for quantitative in vitro determination of glucose in biological fluids.

**PRINCIPLE OF THE METHOD**  
Glucose, in presence of hexokinase, reacts with ATP forming glucose-6-phosphate and ADP. The glucose-6-phosphate reacts with NAD<sup>+</sup> in presence of G-6-PDH to form D-glucono-δ-lactone-6-phosphate and NADH. The intensity of absorbance at 340 nm is proportional to the glucose concentration and can be measured photometrically.

**KIT COMPONENTS**  
**For in vitro diagnostic use only.**  
The components of the kit are stable until expiration date on the label.  
Keep away from direct light sources.

<b>GLU-UV R1</b>	<b>2H251: 4 x 50 ml (liquid) white cap</b>
	<b>6U421: 6 x 56 ml (liquid) white cap</b>
<b>GLU-UV R2</b>	<b>2H251: 2 x 25 ml (liquid) red cap</b>
	<b>6U421: 6 x 14 ml (liquid) red cap</b>

Composition in the test: TRIS pH 7.40 80 mM, MgCl<sub>2</sub> 5 mM, ATP 2mM, NAD 2 mM, hexokinase > 2 KU/l, glucose-6-phosphate dehydrogenase > 2 KU/l.  
Store all components at 2-8°C.

**REAGENT PREPARATION**  
Use separate reagent ready to use.  
Stability: up to expiration date on labels at 2-8°C.  
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-  
Caution: keep well refrigerated.

**PRECAUTIONS**  
Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully, avoiding contact with skin and swallow.  
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

Serum, plasma, urine, CSF (cerebrospinal fluid).  
Separated and nonhemolyzed samples are stable 8 hours at 25°C and 3 days at 2-8°C. Variable stability is observed with longer storage periods. Glycolysis decreases serum glucose by approximately 5 to 7% in 1 h (5 to 10 mg/dl) in normal uncentrifuged coagulated blood at room temperature. The rate of in vitro glycolysis is higher in the presence of leukocytosis or bacterial contamination. Plasma, removed from the cells after moderate centrifugation, contains leukocytes that also metabolize glucose, although cell-free sterile plasma has no glycolytic activity. Glycolysis can be inhibited and glucose stabilized for as long as 3 d at room temperature by adding sodium iodoacetate or sodium fluoride (NaF) to the specimen. Although fluoride maintains long-term blood glucose stability, the rate of decline in the first hour after sample collection is not altered. Cerebrospinal fluid (CSF) may be contaminated with bacteria or other cells and should be analyzed for glucose immediately. If a delay in measurement is unavoidable, the sample should be centrifuged and stored at 4°C or -20 °C. In 24-h collections of urine, glucose may be preserved by adding 5 ml of glacial acetic acid to the container before starting the collection. The final pH of the urine is usually between 4 and 5, which inhibits bacterial activity. Urine samples may lose as much as 40% of their glucose after 24 h at room temperature.

**EXPECTED VALUES**

Plasma/serum (fasting patient)	
adults:	70 - 105 mg/dl
children:	70 - 105 mg/dl
premature neonates:	25 - 80 mg/dl
term neonates:	30 - 90 mg/dl
CSF:	40 - 75 mg/dl
	(60% of plasma value)
Urine (fasting patient)	
random urine:	< 30 mg/dl
24h urine:	< 500 mg/24h

Each laboratory should establish appropriate reference intervals related to its population.

**QUALITY CONTROL AND CALIBRATION**  
It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:  
**QUANTINORM CHEMA**  
with normal or close to normal control values  
**QUANTIPATH CHEMA**  
with pathological control values.  
If required, a multiparametric, human based calibrator is available:  
**AUTOCAL H**

Please contact Customer Care for further information.

**TEST PERFORMANCE**  
**Linearity**  
the method is linear up to 700 mg/dl.  
If the limit value is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

**Sensitivity/limit of detection (LOD)**  
the limit of detection is 1 mg/dl.

**Interferences**  
no interference was observed by the presence of:  
hemoglobin ≤ 500 mg/dl  
bilirubin ≤ 30 mg/dl  
lipids ≤ 1000 mg/dl

**Precision**

intra-assay (n=10)	mean (mg/dl)	SD (mg/dl)	CV%
sample 1	95.20	1.32	1.40
sample 2	224.30	2.36	1.10

inter-assay (n=20)	mean (mg/dl)	SD (mg/dl)	CV%
sample 1	96.47	2.78	2.90
sample 2	252.06	9.56	3.80

**Methods comparison**  
a comparison between Chema and a commercially available product gave the following results:

Glucose UV FL Chema = x  
 Glucose competitor = y  
 n = 100  
 $y = 0.953x + 1.05 \text{ mg/dl} \quad r^2 = 0.99$

**WASTE DISPOSAL**  
This product is made to be used in professional laboratories.  
P501: Dispose of contents according to national/international regulations.