APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN HITACHI 911/912	
TEST:	IGM
APP. CODE:	355
WAVELENGTH (Sec/Pri):	340
ASSAY:	2 POINT END <span style="float: right;">TIME: 10 POINT: 16 - 31 DILUENT: water</span>
SAMPLE VOL:	NORMAL: 3 DECREASE: 2 INCREASE: 5
R1 VOLUME: 240 R2 VOLUME: 0 R3 VOLUME: 60 R4 VOLUME: 0	DILUENT: 5
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LOGIT-LOG (5P) (POINT: 6 - SPAN: 6 - WEIGHT: 0)
SD LIMIT:	0.50
DUPLICATE LIMIT:	10%
ST. 1 CONC:	0.0
EXPECTED VALUE:	0.0 - 2.3
UNIT:	g/l
INSTR. FACTOR (y=ax+b):	a=1 b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN OLYMPUS AU 400/480/600/640/680/2700 (Test code 888)	
TEST NAME:	IGM
SAMPLE:	Volume 3 µl <span style="float: right;">Dilution 0 µl</span>
REAGENTS:	R1 Volume 240 µl <span style="float: right;">Dilution 0 µl</span> R2 Volume 60 µl <span style="float: right;">Dilution 0 µl</span>
WAVELENGTH:	Pri. 340 Sec.
METHOD:	END
REACTION SLOPE:	+
MEASURING POINT 1:	First 0 <span style="float: right;">Last 27</span>
MEASURING POINT 2:	First 0 <span style="float: right;">Last 10</span>
REAGENT OD LIMIT:	First L -0.1 <span style="float: right;">First H 0.5</span> Last L -0.1 <span style="float: right;">Last H 0.5</span>
DYNAMIC RANGE:	L 0.02 <span style="float: right;">H 4.0</span>
CORRELATION FACTOR:	A 1 <span style="float: right;">B 0</span>
UNIT:	g/l
CALIBRATION TYPE:	5AB
FORMULA:	POLYGONAL

 Chema Diagnostica  
 Via Campania 2/4  
 60030 Monsano (AN) - ITALY - EU  
 phone +39 0731 605064  
 fax +39 0731 605672  
 e-mail: mail@chema.com  
 website: http://www.chema.com

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**INTENDED USE**  
 Reagent for quantitative in vitro determination of IgM in biological fluids.

**PRINCIPLE OF THE METHOD**  
 Immunoglobulins M (IgM) selectively react with an anti-IgM antibody and form an immunocomplex. The produced turbidity is proportional to the concentration of IgM in the sample and can be measured at the wavelength of 340 nm.

**KIT COMPONENTS**  
**For in vitro diagnostic use only.**  
 The components of the kit are stable until expiration date on the label.  
 Keep away from direct light sources at 2-8°C.

**IGM R1** 2H100 2 x 40 ml (liquid) white cap  
 6U140 2 x 56 ml (liquid) white cap

Composition: Buffer pH 7.50, PEG ≥ 2%, stabilizers and preservatives.

**IGM R2** 2H100 2 x 10 ml (liquid) red cap  
 6U140 2 x 14 ml (liquid) red cap

Composition: Anti-human IgM antibody ≥ 2%, stabilizers and preservatives.

**REAGENT PREPARATION**  
 Use separate reagent ready to use.  
 Stability: up to expiration date on labels at 2-8°C.  
 Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-  
 Caution: keep well refrigerated.

**PRECAUTIONS**  
**IGM R1: Danger.** Causes serious eye damage (H318).  
 Wear protective gloves. Eye protection (P280).  
 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). Immediately call a doctor (P310).

**IGM R2:** It is not classified as hazardous.

**SPECIMEN**  
 Serum, plasma.  
 Keep specimens away from direct light sources.  
 Samples are stable 7 days when stored at 2-8°C and 1 month at -20°C.

EXPECTED VALUES	
Newborns	0.05-0.3 g/l
Adults	0.4-2.3 g/l

Each laboratory should establish appropriate reference intervals related to its population

**QUALITY CONTROL AND CALIBRATION**  
 It is suggested to perform an internal quality control. For this purpose the following human based control serum is available:  
**QUANTINORM CHEMA**  
 with normal or close to normal control values.  
 If required, a multiparametric, human based calibrator is available:  
**REFERENCE P MULTICALIBRATOR**

Please contact Customer Care for further information.

**TEST PERFORMANCE**

**Measure interval**  
 Measure interval depends on the concentration of the highest standard used for calibration.  
 If such a limit value is exceeded, it is suggested to dilute sample 1+4 with distilled water and to repeat the test, multiplying the result by 5.

**Hook effect**  
 No Hook effect is observed with concentrations lower than 30 g/l.

**Sensitivity/limit of detection**  
 The limit of detection is 0.016 g/l.

**Interferences**  
 No interference was observed by the presence of:  
 hemoglobin ≤ 1000 mg/dl  
 bilirubin ≤ 45 mg/dl  
 lipids ≤ 770 mg/dl  
 rheumatoid factor ≤ 630 IU/ml

**Precision**  
 intra-assay (n=10) mean (g/l) SD (g/l) CV%  
 sample 1 0.70 0.005 0.69  
 sample 2 1.40 0.009 0.66

inter-assay (n=20) mean (g/l) SD (g/l) CV%  
 sample 1 0.70 0.024 3.38  
 sample 2 1.40 0.061 4.32

**Methods comparison**  
 A comparison between CHEMA and a commercially available product gave the following results:

$$\begin{aligned}
 \text{IgM competitor} &= x \\
 \text{IgM FL CHEMA} &= y \\
 n &= 20 \\
 y &= 1.186x - 0.058 \text{ g/l} \quad r^2 = 0.99
 \end{aligned}$$

**WASTE DISPOSAL**  
 This product is made to be used in professional laboratories.  
 P501: Dispose of contents according to national/international regulations.

IgM FL	
GM 2H100	2 x 40 + 2 x 10 ml
GM 6U140	2 x 56 + 2 x 14 ml

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