APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN HITACHI 911/912	
TEST:	<b>IGM</b>
APP. CODE:	<b>355</b>
WAVELENGTH (Sec/Pri):	<b>340</b>
ASSAY:	<b>2 POINT END</b> <i>TIME: 10 POINT: 16 - 31 DILUENT: water</i>
SAMPLE VOL:	NORMAL: <b>3</b> DECREASE: <b>2</b> INCREASE: <b>5</b>
R1 VOLUME:	<b>240</b>
R2 VOLUME:	<b>0</b>
R3 VOLUME:	<b>60</b> <i>DILUENT: 5</i>
R4 VOLUME:	<b>0</b>
ABS LIMIT:	<b>32000 - INC</b>
PROZONE LIMIT:	<b>0 - UPPER</b>
CALIB METHOD:	<b>LOGIT-LOG (5P)</b> <i>(POINT: 6 - SPAN: 6 - WEIGHT: 0)</i>
SD LIMIT:	<b>0.50</b>
DUPLICATE LIMIT:	<b>10%</b>
ST. 1 CONC:	<b>0.0</b>
EXPECTED VALUE:	<b>0.0 - 2.3</b>
UNIT:	<b>g/l</b>
INSTR. FACTOR (y=ax+b):	a=1    b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN OLYMPUS AU 400/480/600/640/680/2700 (Test code 888)	
TEST NAME:	<b>IGM</b>
SAMPLE:	Volume <b>3</b> $\mu$ l      Dilution <b>0</b> $\mu$ l
REAGENTS:	R1 Volume <b>240</b> $\mu$ l      Dilution <b>0</b> $\mu$ l R2 Volume <b>60</b> $\mu$ l      Dilution <b>0</b> $\mu$ l
WAVELENGTH:	Pri. <b>340</b> Sec.
METHOD:	<b>END</b>
REACTION SLOPE:	<b>+</b>
MEASURING POINT 1:	First <b>0</b> Last <b>27</b>
MEASURING POINT 2:	First <b>0</b> Last <b>10</b>
REAGENT OD LIMIT:	First L <b>-0.1</b> First H <b>0.5</b> Last L <b>-0.1</b> Last H <b>0.5</b>
DYNAMIC RANGE:	L <b>0.02</b> H <b>4.0</b>
CORRELATION FACTOR:	A <b>1</b> B <b>0</b>
UNIT:	<b>g/l</b>
CALIBRATION TYPE:	<b>5AB</b>
FORMULA:	<b>POLYGONAL</b>

 Chema Diagnostica  
 Via Campania 2/4  
 60030 Monsano (AN) - ITALY - EU  
 phone +39 0731 605064  
 fax +39 0731 605672  
 e-mail: mail@chema.com  
 website: http://www.chema.com

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**INTENDED USE**  
Reagent for quantitative in vitro determination of IgM in biological fluids.

**PRINCIPLE OF THE METHOD**  
Immunoglobulins M (IgM) selectively react with an anti-IgM antibody and form an immunocomplex. The produced turbidity is proportional to the concentration of IgM in the sample and can be measured at the wavelength of 340 nm.

**KIT COMPONENTS**  
**For in vitro diagnostic use only.**  
The components of the kit are stable until expiration date on the label.  
Keep away from direct light sources at 2-8°C.

**IGM R1** 2H100 2 x 40 ml (liquid) white cap  
6U140 2 x 56 ml (liquid) white cap

Composition: Buffer pH 7.50, PEG  $\geq$  2%, stabilizers and preservatives.

**IGM R2** 2H100 2 x 10 ml (liquid) red cap  
6U140 2 x 14 ml (liquid) red cap

Composition: Anti-human IgM antibody  $\geq$  2%, stabilizers and preservatives.

**REAGENT PREPARATION**  
Use separate reagent ready to use.  
Stability: up to expiration date on labels at 2-8°C.  
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-  
Caution: keep well refrigerated.

**PRECAUTIONS**  
**IGM R1: Danger.** Causes serious eye damage (H318). Wear protective gloves. Eye protection (P280).  IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). Immediately call a doctor (P310).

**IGM R2:** It is not classified as hazardous.

**SPECIMEN**  
Serum, plasma.  
Keep specimens away from direct light sources.  
Samples are stable 7 days when stored at 2-8°C and 1 month at -20°C.

EXPECTED VALUES	
Newborns	0.05-0.3 g/l
Adults	0.4-2.3 g/l

Each laboratory should establish appropriate reference intervals related to its population

**QUALITY CONTROL AND CALIBRATION**  
It is suggested to perform an internal quality control. For this purpose the following human based control serum is available:  
**QUANTINORM CHEMA**  
with normal or close to normal control values.  
If required, a multiparametric, human based calibrator is available:  
**REFERENCE P MULTICALIBRATOR**

Please contact Customer Care for further information.

**TEST PERFORMANCE**

**Measure interval**  
Measure interval depends on the concentration of the highest standard used for calibration.  
If such a limit value is exceeded, it is suggested to dilute sample 1+4 with distilled water and to repeat the test, multiplying the result by 5.

**Hook effect**  
No Hook effect is observed with concentrations lower than 30 g/l.

**Sensitivity/limit of detection**  
The limit of detection is 0.016 g/l.

**Interferences**  
No interference was observed by the presence of:  
hemoglobin  $\leq$  1000 mg/dl  
bilirubin  $\leq$  45 mg/dl  
lipids  $\leq$  770 mg/dl  
rheumatoid factor  $\leq$  630 IU/ml

**Precision**  
intra-assay (n=10) mean (g/l) SD (g/l) CV%  
sample 1 0.70 0.005 0.69  
sample 2 1.40 0.009 0.66  
inter-assay (n=20) mean (g/l) SD (g/l) CV%  
sample 1 0.70 0.024 3.38  
sample 2 1.40 0.061 4.32

**Methods comparison**  
A comparison between CHEMA and a commercially available product gave the following results:

$$y = 1.186x - 0.058 \text{ g/l} \quad r^2 = 0.99$$

**WASTE DISPOSAL**  
This product is made to be used in professional laboratories.  
P501: Dispose of contents according to national/international regulations.

<b>IgM FL</b>	
GM 2H100	2 x 40 + 2 x 10 ml
GM 6U140	2 x 56 + 2 x 14 ml

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