APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912	
TEST:	AMY-P
APP. CODE:	335
WAVELENGTH (Sec/Pri):	700 - 415
ASSAY:	RATE-A <i>TIME: 10 POINT: 24 - 30 DILUENT: water</i>
SAMPLE VOL:	NORMAL: 7 DECREASE: 5 INCREASE: 10
R1 VOLUME:	280 <i>DILUENT: 5</i>
R2 VOLUME:	0
R3 VOLUME:	70 <i>DILUENT: 5</i>
R4 VOLUME:	0
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LINEAR (<i>POINT: 2 - SPAN: 2 - WEIGHT: 0</i>)
SD LIMIT:	0.250
DUPLICATE LIMIT:	3%
ST. 1 CONC:	0.00
EXPECTED VALUE:	0 - 55
UNIT:	U/l
INSTR. FACTOR (y=ax+b):	a=1 b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 860)	
TEST NAME:	AMY-P
SAMPLE:	Volume 5 µl Dilution 0 µl
REAGENTS:	R1 Volume 200 µl Dilution 0 µl R2 Volume 50 µl Dilution 0 µl
WAVELENGTH:	Pri. 410 Sec. 700
METHOD:	RATE
REACTION SLOPE:	+
MEASURING POINT 1:	First 18 Last 26
MEASURING POINT 2:	First Last
REAGENT OD LIMIT:	First L -0.1 First H 0.5 Last L -0.1 Last H 0.5
DYNAMIC RANGE:	L 2 H 2500
CORRELATION FACTOR:	A 1 B 0
LINEARITY LIMIT:	15%
UNIT:	U/l
CALIBRATION TYPE:	AB
FORMULA:	Y = AX + B

 Chema Diagnostics
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

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INTENDED USE

Reagent for quantitative in vitro determination of pancreatic isoamylase in biological fluids.

PRINCIPLE OF THE METHOD

The enzyme α -amylase (EC 3.2.1.1, 1,4 α -D-glucose glucanohydrolase) hydrolyzes the EPS to release several different fragments. The fragments so formed are completely hydrolyzed to 4-nitrophenol and glucose by α -glucosidase. The selective inhibition of S-amylase is performed by use of two different monoclonal antibodies.

The 4-nitrophenol formed is detected spectrophotometrically at 405 nm to give a measurement of α -amylase activity in the sample.

The present method has been made according to IFCC.

KIT COMPONENTS

For in vitro diagnostic use only.

The components of the kit are stable until expiration date on the label.

Keep away from direct light sources.

AMY-P R1 2H100: 4 x 20 ml (liquid) white cap
6U140: 2 x 56 ml (liquid) white cap

AMY-P R2 2H100: 2 x 10 ml (liquid) red cap
6U140: 2 x 14 ml (liquid) red cap

Composition in the test: Hepes buffer 50 mM pH 7.10, NaCl 70 mM, calcium acetate 1.0 mM, 4,6-Ethylidene(G1)-4-nitrophenyl(G7)- α -(1->4)-D-maltoheptaoside 5.0 mM, α -glucosidase 6 kU/l, monoclonal antibodies (mouse) \geq 25 mg/l.

Store all components at 2-8°C.

REAGENT PREPARATION

Use separate reagent ready to use.

Stability: up to expiration date on labels at 2-8°C.

Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-.

Caution: keep well refrigerated.

PRECAUTIONS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.

Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

SPECIMEN

Serum, plasma (heparinate only).

Amylase is stable in serum and plasma sample up to 2 months at 2-8°C.

EXPECTED VALUES

Serum - plasma:	13 - 53 U/l (0.22 - 0.88 μ kat/l)
Random urine:	\leq 350 U/l (\leq 5.84 μ kat/l)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

QUANTINORM CHEMA

with normal or close to normal control values

QUANTIPATH CHEMA

with pathological control values.

If required, a multiparametric, human based calibrator is available:

AUTOCAL H

Please contact Customer Care for further information.

TEST PERFORMANCE

Linearity

the method is linear up to 2500 U/l.

If a $\Delta A/\text{min}$ of 0.500 is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)

the limit of detection is 2 U/l.

Interferences

no interference was observed by the presence of:

hemoglobin \leq 500 mg/dl

bilirubin \leq 25 mg/dl

lipids interfer in low values

Precision

intra-assay (n=10)	mean (U/l)	SD (U/l)	CV%
sample 1	38.00	0.67	1.80
sample 2	103.00	1.41	1.40

inter-assay (n=20)	mean (U/l)	SD (U/l)	CV%
sample 1	38.71	0.97	2.50
sample 2	102.61	1.62	1.60

Methods comparison

a comparison between Chema Isoamylase FL and a commercially available product gave the following results:

Isoamylase Chema = x
Isoamylase competitor = y
n = 108

$$y = 1.02x - 0.605 \text{ U/l} \quad r^2 = 0.997$$

WASTE DISPOSAL

This product is made to be used in professional laboratories.

P501: Dispose of contents according to national/international regulations.

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