APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ
HITACHI 911/912

TEST:	ALP DGKC	
APP. CODE:	342	
WAVELENGTH (Sec/Pri):	700 - 415	
ASSAY:	RATE-A	<i>TIME: 10 POINT: 20 - 30 DILUENT: water</i>
SAMPLE VOL:	NORMAL: 7 DECREASE: 5 INCREASE: 10	
	R1 VOLUME: 280	<i>DILUENT: 5</i>
	R2 VOLUME: 0	
	R3 VOLUME: 70	<i>DILUENT: 5</i>
	R4 VOLUME: 0	
ABS LIMIT:	32000 - INC	
PROZONE LIMIT:	0 - UPPER	
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)	
SD LIMIT:	0.250	
DUPLICATE LIMIT:	3%	
ST. 1 CONC:	0.0	
EXPECTED VALUE:	0 - 270	
UNIT:	U/l	
INSTR. FACTOR (y=ax+b):	a=1	b= 0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ
OLYMPUS AU 400/480/600/640/680/2700 (Test code 872)

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INTENDED USE

Reagent for quantitative in vitro determination of alkaline phosphatase in biological fluids.

PRINCIPLE OF THE METHOD

The enzyme alkaline phosphatase (EC 3.1.3.1, ortho-phosphoric-monoester phosphohydrolase) hydrolyzes the 4-NPP to release 4-nitrophenol, under alkaline conditions. The 4-nitrophenol formed is detected spectrophotometrically at 405 nm to give a measurement of alkaline phosphatase activity in the sample.

The present method has been made according to DGKC.

KIT COMPONENTS

For in vitro diagnostic use only.

The components of the kit are stable until expiration date on the label.
Keep away from direct light sources.

ALP DGKC R1 2H100: 4 x 20 ml (liquid) white cap
2H400: 8 x 40 ml (liquid) white cap
6U420: 6 x 56 ml (liquid) white cap

ALP DGKC R2 2H100: 2 x 10 ml (liquid) red cap
2H400: 4 x 20 ml (liquid) red cap
6U420: 6 x 14 ml (liquid) red cap

Composition in the test: DEA buffer pH 9.8 1 M, MgCl₂ 0,5 mM, 4-Nitrophenylphosphate 10 mM.

Store all components at 2-8°C.

REAGENT PREPARATION

Use separate reagent ready to use.

Stability: up to expiration date on labels at 2-8°C.

Stability since first opening of vials: preferably within 60 days at 2-8°C - away from light sources.

Caution: keep well refrigerated.

PRECAUTIONS

ALP DGKC R1: Danger. May cause damage to organs through prolonged or repeated exposure (H373). Causes serious eye damage (H318). Causes skin irritation (H315). Wear protective gloves. Eye protection (P280). Do not breathe vapours (P260). IF ON SKIN: Wash with plenty of water (P302+P352). IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). Get medical advice if you feel unwell (P314).

ALP DGKC R2: It is not classified as hazardous.

SPECIMEN

Serum, plasma (heparinate only).

Sera kept at room temperatures usually show a slight but real increase in activity, which varies from 1% over a 6-h period to 3 to 6% over a 1 to 4 days period. Even in sera stored at refrigerator temperature, activity increases slowly. In frozen sera, activity decreases but slowly recovers after thawing the serum.

A similar enhancement of activity, but of greater magnitude, occurs with reconstituted lyophilized preparations, such as those available as control sera or calibrators. In reconstituted material the increases with storage at 4 and 20°C are about 10 and 30%, respectively. Enhancement of activity continues for several days, but at a decreasing rate. The cause of this phenomenon is not known but may be attributed to renaturation of partially denatured enzyme or to dissociation, on warming, of a phosphate-lipoprotein complex or a multimer of the enzyme that was formed in the freeze-drying process.

EXPECTED VALUES

Men: < 270 U/l (< 4.50 µkat/l)
Women: < 240 U/l (< 4.00 µkat/l)

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

QUANTINORM CHEMA

with normal or close to normal control values

QUANTIPATH CHEMA

with pathological control values.

If required, a multiparametric, human based calibrator is available:

AUTOCAL H

Please contact Customer Care for further informations.

TEST PERFORMANCE

Linearity

the method is linear up to 3000 U/l.

If a ΔA/min of 0.500 is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)

the limit of detection is 7 U/l.

Interferences

no interference was observed by the presence of:

hemoglobin ≤ 400 mg/dl
bilirubin ≤ 27 mg/dl
lipids ≤ 1000 mg/dl

Precision

intra-assay (n=10)	mean (U/l)	SD (U/l)	CV%
sample 1	175.70	0.95	0.50
sample 2	426.70	2.41	0.60

inter-assay (n=20)	mean (U/l)	SD (U/l)	CV%
sample 1	167.26	3.99	2.40
sample 2	408.28	8.61	2.10

Methods comparison

a comparison between Chema and a commercially available product gave the following results:

ALP Chema = x	ALP competitor = y
n = 112	

y = 0.96x - 2.17 U/l r² = 0.999

WASTE DISPOSAL

This product is made to be used in professional laboratories.

P501: Dispose of contents according to national/international regulations.

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