APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ <b>HITACHI 911/912</b>	
TEST:	<b>ALP IFCC</b>
APP. CODE:	<b>343</b>
WAVELENGTH (Sec/Pri):	<b>700 - 415</b>
ASSAY:	<b>RATE-A</b> <i>TIME: 10</i> <i>POINT: 20 - 30</i> <i>DILUENT: water</i>
SAMPLE VOL:	NORMAL: 7 DECREASE: 5 INCREASE: 10
R1 VOLUME: 280	<i>DILUENT: 5</i>
R2 VOLUME: 0	
R3 VOLUME: 70	<i>DILUENT: 5</i>
R4 VOLUME: 0	
ABS LIMIT:	<b>32000 - INC</b>
PROZONE LIMIT:	<b>0 - UPPER</b>
CALIB METHOD:	<b>LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)</b>
SD LIMIT:	<b>0.250</b>
DUPLICATE LIMIT:	<b>3%</b>
ST. 1 CONC:	<b>0.0</b>
EXPECTED VALUE:	<b>30 - 130</b>
UNIT:	<b>U/l</b>
INSTR. FACTOR (y=ax+b):	a=1    b= 0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ <b>OLYMPUS AU 400/480/600/640/680/2700 (Test code 873)</b>	
TEST NAME:	<b>ALP IFCC</b>
SAMPLE:	Volume 7 µl      Dilution 0 µl
REAGENTS:	R1 Volume 280 µl      Dilution 0 µl R2 Volume 70 µl      Dilution 0 µl
WAVELENGTH:	Pri. 410    Sec. 700
METHOD:	<b>RATE</b>
REACTION SLOPE:	<b>+</b>
MEASURING POINT 1:	First 14      Last 27
MEASURING POINT 2:	First      Last
REAGENT OD LIMIT:	First L -0.1    First H 1.0 Last L -0.1    Last H 1.0
DYNAMIC RANGE:	L 5      H 3000
CORRELATION FACTOR:	A 1      B 0
LINEARITY LIMIT:	<b>15%</b>
UNIT:	<b>U/l</b>
CALIBRATION TYPE:	<b>AB</b>
FORMULA:	<b>Y = AX + B</b>

 Chema Diagnostica  
Via Campania 2/4  
60030 Monsano (AN) - ITALY - EU  
phone +39 0731 605064  
fax +39 0731 605672  
e-mail: mail@chema.com  
website: http://www.chema.com

FOSFATASI ALCALINA FL IFCC	
AF 2H100	4 x 20 + 2 x 10 ml
AF 2H400	8 x 40 + 4 x 20 ml
AF 6U420	6 x 56 + 6 x 14 ml

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**INTENDED USE**  
Reagent for quantitative in vitro determination of alkaline phosphatase in biological fluids.

**PRINCIPLE OF THE METHOD**  
The enzyme alkaline phosphatase (EC 3.1.3.1., orthophosphoric-monoester phosphohydrolase) hydrolyzes the 4-NPP to release 4-nitrophenol, under alkaline conditions. The 4-nitrophenol formed is detected spectrophotometrically at 405 nm to give a measurement of alkaline phosphatase activity in the sample. The present method has been made according to IFCC.

**KIT COMPONENTS**  
**For in vitro diagnostic use only.**  
The components of the kit are stable until expiration date on the label.  
Keep away from direct light sources.

**ALP IFCC R1** 2H100: 4 x 20 ml (liquid) white cap  
2H400: 8 x 40 ml (liquid) white cap  
6U420: 6 x 56 ml (liquid) white cap

**ALP IFCC R2** 2H100: 2 x 10 ml (liquid) red cap  
2H400: 4 x 20 ml (liquid) red cap  
6U420: 6 x 14 ml (liquid) red cap

Composition in the test: 2-amino-2-methyl-1-propanol buffer 0.35 M pH 10.40 (30°C), magnesium acetate 2 mM, zinc sulfate 1 mM, HEDTA 2 mM, 4-NPP 16 mM.

Store all components at 2-8°C.

**REAGENT PREPARATION**  
Use separate reagent ready to use.  
Stability: up to expiration date on labels at 2-8°C.  
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-  
Caution: keep well refrigerated.

**WASTE DISPOSAL**  
This product is made to be used in professional laboratories.  
P501: Dispose of contents according to national/international regulations.

**REFERENCES**  
Clin. Chim. Acta, (1983) 339F - 367F  
Tietz Textbook of Clinical Chemistry, Second Edition, Burtis-Ashwood (1994).

**PRECAUTIONS**  
**ALP IFCC R1: Warning.** Causes serious eye irritation (H319). Causes skin irritation (H315). Wear protective gloves. Eye protection (P280). IF ON SKIN: wash with plenty of water (P302+P352). IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). If eye irritation persists: get medical advice (P337+P313).

**ALP IFCC R2:** It is not classified as hazardous.

**SPECIMEN**  
Serum, plasma (heparinized only).  
Sera kept at room temperatures usually show a slight but real increase in activity, which varies from 1% over a 6-h period to 3 to 6% over a 1 to 4 days period. Even in sera stored at refrigerator temperature, activity increases slowly. In frozen sera, activity decreases but slowly recovers after thawing the serum.  
A similar enhancement of activity, but of greater magnitude, occurs with reconstituted lyophilized preparations, such as those available as control sera or calibrators. In reconstituted material the increases with storage at 4 and 20°C are about 10 and 30%, respectively. Enhancement of activity continues for several days, but at a decreasing rate. The cause of this phenomenon is not known but may be attributed to renaturation of partially denatured enzyme or to dissociation, on warming, of a phosphate-lipoprotein complex or a multimer of the enzyme that was formed in the freeze-drying process.

**EXPECTED VALUES**  
Men: 40 - 129 U/l (0.67 - 2.15 µkat/l)  
Women: 35 - 104 U/l (0.58 - 1.74 µkat/l)

**QUALITY CONTROL AND CALIBRATION**  
It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:  
**QUANTINORM CHEMA** with normal or close to normal control values  
**QUANTIPATH CHEMA** with pathological control values.  
If required, a multiparametric, human based calibrator is available:  
**AUTOCAL H**

Please contact Customer Care for further informations.

**TEST PERFORMANCE**  
**Linearity**  
the method is linear up to 3000 U/l.  
If a ΔA/min of 0.500 is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

**Sensitivity/limit of detection (LOD)**  
the limit of detection is 5.2 U/l.

**Interferences**  
no interference was observed by the presence of:  
hemoglobin ≤ 400 mg/dl  
bilirubin ≤ 40 mg/dl  
lipids ≤ 900 mg/dl

**Precision**  
intra-assay (n=10)    mean (U/l)    SD (U/l)    CV%  
sample 1    84.40    2.41    2.86  
sample 2    222.40    5.74    2.58

inter-assay (n=20)    mean (U/l)    SD (U/l)    CV%  
sample 1    86.66    2.66    3.07  
sample 2    210.39    6.08    2.89

**Methods comparison**  
a comparison between Chema and a commercially available product gave the following results:

ALP Chema = x  
ALP competitor = y  
n = 150  
  
y = 1.03x - 2.57 U/l      r² = 0.998