APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912	
TEST:	<b>HS</b>
APP. CODE:	<b>349</b>
WAVELENGTH (Sec/Pri):	<b>700 - 600</b>
ASSAY:	<b>1-POINT</b> <i>TIME: 10</i> <i>POINT: 31</i>
SAMPLE VOL:	NORMAL: <b>3</b> DECREASE: <b>2</b> INCREASE: <b>5</b>
	R1 VOLUME: <b>300</b> R2 VOLUME: <b>0</b> R3 VOLUME: <b>0</b> R4 VOLUME: <b>0</b>
ABS LIMIT:	<b>32000 - INC</b>
PROZONE LIMIT:	<b>0 - UPPER</b>
CALIB METHOD:	<b>LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)</b>
SD LIMIT:	<b>0.250</b>
DUPLICATE LIMIT:	<b>3%</b>
ST. 1 CONC:	<b>0.00</b>
EXPECTED VALUE:	<b>10 - 150</b>
UNIT:	<b>mg/dl</b>
INSTR. FACTOR (y=ax+b):	a=1    b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 881)	
TEST NAME:	<b>HS</b>
SAMPLE:	Volume <b>3 µl</b> Dilution <b>0 µl</b>
REAGENTS:	R1 Volume <b>300 µl</b> Dilution <b>0 µl</b> R2 Volume <b>0 µl</b> Dilution <b>0 µl</b>
WAVELENGTH:	Pri. <b>600</b> Sec. <b>700</b>
METHOD:	<b>END</b>
REACTION SLOPE:	<b>+</b>
MEASURING POINT 1:	First <b>0</b> Last <b>27</b>
MEASURING POINT 2:	First      Last
REAGENT OD LIMIT:	First L <b>0.01</b> First H <b>1.0</b> Last L <b>0.01</b> Last H <b>1.0</b>
DYNAMIC RANGE:	L <b>0.3</b> H <b>500</b>
CORRELATION FACTOR:	A <b>1</b> B <b>0</b>
UNIT:	<b>mg/dl</b>
CALIBRATION TYPE:	<b>AB</b>
FORMULA:	<b>Y = AX + B</b>

 Chema Diagnostica  
Via Campana 2/4  
60030 Monsano (AN) - ITALY - EU  
phone +39 0731 605064  
fax +39 0731 605672  
e-mail: mail@chema.com  
website: http://www.chema.com

ITALIANO rev. 26/09/2016

PROTEINE HS	
HS 2H200	10 x 20 ml
HS 6U168	3 x 56 ml

#### USO

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#### INTENDED USE

Reagent for quantitative in vitro determination of proteins in biological fluids.

#### PRINCIPLE OF THE METHOD

Proteins combine with pyrogallol red to form a color complex, the absorbance of which is measured at 600 nm. Sodium dodecylsulphate is added to increase accuracy in measuring proteins other than albumine (Watanabe).

#### KIT COMPONENTS

**For in vitro diagnostic use only.**  
The components of the kit are stable until expiration date on the label.  
Keep away from direct light sources.

**HS R1 2H200: 10 x 20 ml (liquid) white cap**  
**6U168: 3 x 56 ml (liquid) white cap**

Composizione: succinate buffer 0.05 M pH 2.50, pyrogallol red 0.04 mM, sodium molibdate 0.13 mM, sodium oxalate 1 mM, sodium benzoate 0.35 mM, SDS 0.1 mM.

Store all components at 2-8°C.

#### REAGENT PREPARATION

Use separate reagent ready to use.  
Stability: up to expiration date on labels at 2-8°C.  
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-.  
Caution: keep well refrigerated.

#### PRECAUTIONS

**HS R1: Warning.** May cause damage to organs (H371). Do not breath vapours (P260). Wash hands thoroughly after handling (P264). Do not eat, drink or smoke when using this product (P270).

#### SPECIMEN

Urine, cerebrospinal fluid. Stable 3 days at 2-8°C.

#### EXPECTED VALUES

Cerebrospinal fluid: 14 - 45 mg/dl  
Urine: 28 - 141 mg/24 h

Each laboratory should establish appropriate reference intervals related to its population.

#### QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose a reliable and compatible urine based control material is needed.

Please contact Customer Care for further information.

#### TEST PERFORMANCE

**Linearity**  
the method is linear up to 500 mg/dl.  
If the limit value is exceeded, it is suggested to dilute sample 1+9 with distilled water and to repeat the test, multiplying the result by 10.

**Sensitivity/limit of detection (LOD)**  
the limit of detection is 0.3 mg/dl.

#### Interferences

no interference was observed by the presence of:  
ascorbic acid ≤ 500 mg/dl

#### Precision

intra-assay (n=10) mean (mg/dl) SD (mg/dl) CV%  
sample 1 37.10 0.74 2.00  
sample 2 103.40 1.27 1.20

inter-assay (n=20) mean (mg/dl) SD (mg/dl) CV%  
sample 1 38.01 0.79 2.00  
sample 2 100.09 2.46 2.00

#### Methods comparison

a comparison between Chema and a commercially available product gave the following results:

Proteins HS Chema = x  
Proteins competitor = y  
n = 88

$$y = 0.97x - 0.54 \text{ mg/dl} \quad r^2 = 0.978$$

#### WASTE DISPOSAL

This product is made to be used in professional laboratories.  
P501: Dispose of contents according to national/international regulations.

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