Componente component composant componente компонент	metodo method méthode método метод	VALORE VALUE VALEUR VALOR ЗНАЧЕНИЕ СHEMA	valore value valor значение	intervallo range intervale intervallo диапазон	1SD 1s	unità units unités unidades единица
GGT γ-Glutamiltransferasi γ-Glutamyltransferase γ-Glutamyltransferase γ-Glutamyltransferasa γ-Глутамилтрансфераза	liquido std. Szasz liquid std. Szasz liquide std. Szasz liquido std. Szasz жидкий станд. Зейца liquido std. IFCC liquid std. IFCC liquide std. IFCC liquido std. IFCC жидкий станд. IFCC		37°C 205 3.42	169 - 241 2.79 - 4.05	12 0.21	U/l Ед./л μkat/l мккат/л
GLDH Glutamato deidrogenasi Glutamate dehydrogenase Glutamate déshydrogénase Glutamato deshidrogenasa Глутамат дегидрогеназа			37°C 43.3 0.723	32.8 - 53.8 0.549 - 0.897	3.5 0.058	U/l Ед./л μkat/l мккат/л
GLUC Glucosio Glucose Glucose Glucosa Глюкоза	HK G6P-PT-CCO-PAP	VALORE VALUE VALEUR VALOR ЗНАЧЕНИЕ СHEMA	256 14.2 2.56	217 - 295 12.1 - 16.3 2.17 - 2.95	13 0.7 0.13	mg/dl мг/дл mmol/l ммоль/л g/l г/л
HBDH α-Hydroxybutyrate deidrogenasi α-Hydroxybutyrate dehydrogenase α-Hydroxybutyrate déshydrogénase α-Hydroxybutyrate deshidrogenasa α-Гидроксибутират дегидрогеназы	DGKC		37°C 331 5.53	271 - 391 4.54 - 6.52	20 0.33	U/l Ед./л μkat/l мккат/л
K Potassio Potassium Potassium Potasio Калий	ISE potenziometria indiretta ISE indirect potentiometry ISE potentiométrie indirecte ISE potenziometria indirecta Непрямая потенциометрия ISE		6.59 25.8	5.99 - 7.19 23.4 - 28.2	0.20 0.8	mmol/l ммоль/л mg/dl мг/дл
LACT L-Lattato L-Lactate L-Lactate L-Lactato L-Лактат	metodo enzimatico colorimetrico enzymatic colorimetric method méthode colorimétrique enzymatique método enzimático colorimétrico энзиматический колориметрический метод		33.2 3.69 332	27.2 - 39.2 3.03 - 4.35 272 - 392	2.0 0.22 20	mg/dl мг/дл mmol/l ммоль/л mg/l мг/л
LDH-L Lattato deidrogenasi Lactate dehydrogenase Lactate déshydrogénase Lactato deshidrogenasa Лактат дегидрогеназы	IFCC liquida/liquid/liquide/liquido/жидкий		37°C 319 5.33	262 - 376 4.37 - 6.29	19 0.32	U/l Ед./л μkat/l мккат/л
LDH-P Lattato deidrogenasi Lactate dehydrogenase Lactate déshydrogénase Lactato deshidrogenasa Лактат дегидрогеназы	DGKC	VALORE VALUE VALEUR VALOR ЗНАЧЕНИЕ СHEMA	37°C 302 10.1	494 - 710 8.3 - 11.9	36 0.6	U/l Ед./л μkat/l мккат/л
Li Litio Lithium Lithium Litio Литий	ISE potenziometria diretta ISE direct potentiometry ISE potentiométrie directe ISE potenziometria directa Прямая потенциометрия ISE		2.03 1.41	1.79 - 2.27 1.23 - 1.59	0.08 0.06	mmol/l ммоль/л mg/dl мг/дл
LIP Lipasi Lipase Lipase Lipasa Липаза	metodo enzimatico colorimetrico enzymatic colorimetric method méthode colorimétrique enzymatique método enzimático colorimétrico энзиматический колориметрический метод	VALORE VALUE VALEUR VALOR ЗНАЧЕНИЕ СHEMA	37°C 108 1.80	90 - 126 1.47 - 2.13	6 0.11	U/l Ед./л μkat/l мккат/л
MG Magnesio Magnesium Magnesium Magnésium Magnesio Магний	blu di xilidide xylydyl blue bleu de xylydyle azul de xilidil ксилидил синий	VALORE VALUE VALEUR VALOR ЗНАЧЕНИЕ СHEMA	1.41 2.82 3.43	1.23 - 1.59 2.49 - 3.15 3.01 - 3.85	0.06 0.11 0.14	mmol/l ммоль/л μkat/l мккат/л mg/dl мг/дл

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LOT

Test procedure

Use supplied materials as provided.

Additional materials required: automatic micropipette, distilled/deionized water, general laboratory equipment.

Assay

Dispense the required volume into a sample cup and analyze in the same way as for the samples. The control should be run daily in parallel with the patient samples and after every calibration. The control intervals should be adapted to each laboratory's individual requirements. Each laboratory should establish QC procedures that conform with local, state and/or federal regulations or accreditation required.

Assigned values and ranges

Target value determination procedure: single assays in 7 external laboratories, with ≥ 3 independent series, under the direction of manufacturer department for external value assignment; this department also performed the evaluation.

The declared value is the median of all values obtained. The corresponding control range is calculated as the target value ± 3 SD, with the standard deviation (SD) being the value obtained from several target value determinations.

Determinations were made using the testing reagents and analysis systems available at the time. Results must be within the defined ranges. Each laboratory should establish guidelines for corrective measures to be taken if values fall outside the range.

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English

HUMAN BASED UNIVERSAL CONTROL SERUM

Intended use

Quantipath Chema is for use in the quality control of general clinical chemistry methods for the quantitative determination of substrates, electrolytes, lipids, enzymes and proteins. The control is used to monitor accuracy or precision both for manual techniques and assays on automated clinical chemistry analyzers.

Summary

Quantipath Chema is a lyophilized control serum based on human serum with concentrations / activities found in the pathological range.

Kit components

Lyophilized control serum 5 x 5 ml / 10 x 5 ml

Reactive lyophilized components

Human serum with chemical additives and tissue extracts of human and animal origin. The origin of biological additives is as follows:

Alkaline phosphatase	placenta (human, recombinant)
ALT/GPT	porcine heart
AST/GOT	human, recombinant
Albumin	human serum
Amylase, total	human saliva / porcine pancreas
Amylase, pancreatic	porcine pancreas
Cholesterol	bovine plasma
Cholinesterase	human serum
Creatine kinase	rabbit muscle
Gamma-GT	human, recombinant
LDH	porcine heart
Lipasi	pancreas (human, recombinant)
GLDH	bacterial, recombinant
Total protein	human serum

The concentrations/activities are lot specific. The exact values and ranges are listed in the enclosed values sheet.

Precautions and warnings

For *in vitro* diagnostic use only.

Exercise the normal precautions required for handling all laboratory reagents. This product has been prepared exclusively from the blood of donors tested individually and shown by FDA approved methods to be free from HBsAg and antibodies to HCV and HIV. However as no test method can rule out the potential risk of infection with absolute certainty, the material should be handled just as carefully as a patient sample. In the event of exposure the directive of the responsible health authorities should be followed. (a,b)

Disposal

This product is made to be used in professional laboratories.

Please consult local regulations for a correct waste disposal.

P501: Dispose of contents according to national/international regulations.

Handling

Carefully open one bottle of serum, avoiding the loss of lyophilizate, and pipette exactly 5.0 ml of distilled / deionized water. Carefully close the bottle and dissolve the contents completely by occasional gentle swirling within 30 minutes. Avoid the formation of foam. **Important:** with the exception of **alkaline phosphatase**, all enzymes can be measured immediately. To reactivate the alkaline phosphatase, allow the reconstituted control to stand for one hour at +25°C.

Storage and stability

Store at 2-8°C.

Criterion of the stability data specified by the manufacturer: recovery within ±10% of initial value. Stability of lyophilized calibrator serum at 2-8°C: up to the stated expiration date.

STABILITY IN RECONSTITUTED CONTROL SERUM			
	Components	Total bilirubin (protected from light)	Direct bilirubin (protected from light)
15-25°C	12 hours	8 hours	4 hours
2-8°C	5 days	1 day	8 hours
(-15)-(-25)°C	1 month (when frozen once)	2 weeks (when frozen once)	2 weeks (when frozen once)

The possible appearance of a slight green coloration has no effect on the recovery of the values. Store controls tightly capped when not in use.