APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN HITACHI 911/912	
TEST:	TRF
APP. CODE:	356
WAVELENGTH (Sec/Pri):	505
ASSAY:	2 POINT END <span style="float: right;">TIME: 10 POINT: 16 - 31 DILUENT: water</span>
SAMPLE VOL:	NORMAL: 2 DECREASE: 2 INCREASE: 5 <span style="float: right;">DILUENT: 5</span>
R1 VOLUME:	240
R2 VOLUME:	0
R3 VOLUME:	60
R4 VOLUME:	0
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LOGIT-LOG (5P) (POINT: 6 - SPAN: 6 - WEIGHT: 0)
SD LIMIT:	0.50
DUPLICATE LIMIT:	10%
ST. 1 CONC:	0.0
EXPECTED VALUE:	1.0 - 4.0
UNIT:	g/l
INSTR. FACTOR (y=ax+b):	a=1 b= 0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN OLYMPUS AU 400/480/600/640/680/2700 (Test code 889)	
TEST NAME:	TRF
SAMPLE:	Volume 2 µl <span style="float: right;">Dilution 0 µl</span>
REAGENTS:	R1 Volume 240 µl <span style="float: right;">Dilution 0 µl</span> R2 Volume 60 µl <span style="float: right;">Dilution 0 µl</span>
WAVELENGTH:	Pri. 520 Sec.
METHOD:	END
REACTION SLOPE:	+
MEASURING POINT 1:	First 0 <span style="float: right;">Last 27</span>
MEASURING POINT 2:	First 0 <span style="float: right;">Last 10</span>
REAGENT OD LIMIT:	First L -0.1 <span style="float: right;">First H 0.5</span> Last L -0.1 <span style="float: right;">Last H 0.5</span>
DYNAMIC RANGE:	L 0.01 <span style="float: right;">H 9.0</span>
CORRELATION FACTOR:	A 1 <span style="float: right;">B 0</span>
UNIT:	g/l
CALIBRATION TYPE:	5AB
FORMULA:	POLYGONAL

 Chema Diagnostica  
Via Campania 2/4  
60030 Monsano (AN) - ITALY - EU  
phone +39 0731 605064  
fax +39 0731 605672  
e-mail: mail@chema.com  
website: http://www.chema.com

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**INTENDED USE**  
Reagent for quantitative in vitro determination of transferrin in biological fluids.

**PRINCIPLE OF THE METHOD**  
Transferrin selectively reacts with an anti-transferrin antibody and forms an immunocomplex. The produced turbidity is proportional to the concentration of transferrin in the sample, and can be measured at the wavelength of 505 nm.

**KIT COMPONENTS**  
**For in vitro diagnostic use only.**  
The components of the kit are stable until expiration date on the label.  
Keep away from direct light sources.

**TRF R1** 2H100: 2 x 40 ml (liquid) white cap  
6U140: 2 x 56 ml (liquid) white cap

Composizione: Buffer pH 7.00, PEG ≥ 2%, stabilizers and preservatives.

**TRF R2** 2H100: 2 x 10 ml (liquid) red cap  
6U140: 2 x 14 ml (liquid) red cap

Composizione: Anti-human transferrin antibody ≥ 2%, stabilizers and preservatives.

Store all components at 2-8°C.

**REAGENT PREPARATION**  
Use separate reagent ready to use.  
Stability: up to expiration date on labels at 2-8°C.  
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-  
Caution: keep well refrigerated.

**PRECAUTIONS**  
Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.  
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

**SPECIMEN**  
Serum, plasma.  
Keep specimens away from direct light sources.  
Samples are stable 7 days when stored at 2-8°C and 1 month at -20°C.

**EXPECTED VALUES**  
Newborns 1.2-2.5 g/l  
Adults 2.0-3.6 g/l  
Each laboratory should establish appropriate reference intervals related to its population

**QUALITY CONTROL AND CALIBRATION**  
It is suggested to perform an internal quality control. For this purpose the following human based control serum is available:  
**QUANTINORM CHEMA**  
with normal or close to normal control values.  
If required, a multiparametric, human based calibrator is available:  
**REFERENCE P MULTICALIBRATOR**

Please contact Customer Care for further information.

**TEST PERFORMANCE**  
**Measure interval**  
Measure interval depends on the concentration of the highest standard used for calibration.  
If such a limit value is exceeded, it is suggested to dilute sample 1+4 with distilled water and to repeat the test, multiplying the result by 5.

**Hook effect**  
No Hook effect is observed with concentrations lower than 40 g/l.

**Sensitivity/limit of detection**  
The limit of detection is 0.01 g/l.

**Interferences**  
No interference was observed by the presence of:  
hemoglobin ≤ 1000 mg/dl  
bilirubin ≤ 32 mg/dl  
lipids ≤ 1200 mg/dl  
rheumatoid factor ≤ 630 IU/ml

**Precision**  
intra-assay (n=10) mean (g/l) SD (g/l) CV%  
sample 1 2.20 0.02 1.15  
sample 2 4.13 0.04 0.99  
inter-assay (n=10) mean (g/l) SD (g/l) CV%  
sample 1 2.19 0.05 2.36  
sample 2 4.12 0.08 2.05

**Methods comparison**  
A comparison between CHEMA and a commercially available product gave the following results:

Transferrin competitor = x  
Transferrin FL CHEMA = y  
n = 82  
 $y = 1.004x + 0.018 \text{ g/l}$   $r^2 = 0.991$

**WASTE DISPOSAL**  
This product is made to be used in professional laboratories.  
P501: Dispose of contents according to national/international regulations.

TRANSFERRINE FL	
TF 2H100	2 x 40 + 2 x 10 ml
TF 6U140	2 x 56 + 2 x 14 ml

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