APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912	
TEST:	UREA
APP. CODE:	392
WAVELENGTH (Sec/Pri):	700 - 340
ASSAY:	2 POINT RATE <i>TIME: 10 POINT: 19 - 23 DILUENT: water</i>
SAMPLE VOL:	NORMAL: 3 DECREASE: 2 INCREASE: 5
R1 VOLUME:	240 <i>DILUENT: 5</i>
R2 VOLUME:	0
R3 VOLUME:	60 <i>DILUENT: 5</i>
R4 VOLUME:	0
ABS LIMIT:	6500 - DEC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)
SD LIMIT:	0.250
DUPLICATE LIMIT:	3%
ST. 1 CONC:	0.00
EXPECTED VALUE:	10 - 50
UNIT:	mg/dl
INSTR. FACTOR (y=ax+b):	a=1 b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 884)	
TEST NAME:	UREA
SAMPLE:	Volume 3 µl Dilution 0 µl
REAGENTS:	R1 Volume 240 µl Dilution 0 µl R2 Volume 60 µl Dilution 0 µl
WAVELENGTH:	Pri. 340 Sec. 700
METHOD:	RATE
REACTION SLOPE:	-
MEASURING POINT 1:	First 13 Last 17
MEASURING POINT 2:	First Last
REAGENT OD LIMIT:	First L 0.8 First H 2.0 Last L 0.8 Last H 2.0
DYNAMIC RANGE:	L 1 H 300
CORRELATION FACTOR:	A 1 B 0
LINEARITY LIMIT:	15%
UNIT:	mg/dl
CALIBRATION TYPE:	AB
FORMULA:	Y = AX + B

 Chema Diagnostica
 Via Campania 2/4
 60030 Monsano (AN) - ITALY - EU
 phone +39 0731 605064
 fax +39 0731 605672
 e-mail: mail@chema.com
 website: http://www.chema.com

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ENGLISH	
UREA UV FL	
AZ 2H400	8 x 40 + 4 x 20 ml
AZ 6U420	6 x 56 + 6 x 14 ml
INTENDED USE	
Reagent for quantitative in vitro determination of urea in biological fluids.	
PRINCIPLE OF THE METHOD	
The urease hydrolyzes urea in sample to release ammonium ions, which react with 2-oxoglutarate and NADH in presence of glutamate dehydrogenase to form glutamate and NAD ⁺ . The decrease of absorbance is measured at 340 nm.	
KIT COMPONENTS	
For in vitro diagnostic use only.	
The components of the kit are stable until expiration date on the label.	
Keep away from direct light sources.	
UREA R1	2H400: 8 x 40 ml (liquid) white cap 6U420: 6 x 56 ml (liquid) white cap
UREA R2	2H400 4 x 20 ml (liquid) red cap 6U420: 6 x 14 ml (liquid) red cap
Composition in the test: CAPSO buffer 8 mM pH 7.60, 2-Oxoglutarate 7.5 mM, Urease > 8 KU/I, GLDH > 800 U/I, NADH 0.25 mM, stabilizers.	
Store all components at 2-8°C.	
REAGENT PREPARATION	
Use separate reagent ready to use.	
Stability: up to expiration date on labels at 2-8°C.	
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-	
Caution: keep well refrigerated.	

PRECAUTIONS	
UREA R1: It is not classified as hazardous.	
UREA R2: Warning. Causes serious eye irritation (H319). Causes skin irritation (H315). Wear protective gloves. Eye protection (P280). IF ON SKIN: Wash with plenty of water (P302+P352). IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). If eye irritation persists: get medical advice (P337+P313).	
SPECIMEN	
Serum, plasma (avoid ammonium heparinate). Urine. Urea is stable 3 days at 2-8°C.	
Dilute urine sample 1:100 with deionized water.	
EXPECTED VALUES	
Adults:	10 - 50 mg/dl (1.7 - 8.3 mmol/l)
Urine:	20 - 35 g/24h (332 - 580 mmol/24h)
Each laboratory should establish appropriate reference intervals related to its population.	
QUALITY CONTROL AND CALIBRATION	
It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:	
QUANTINORM CHEMA with normal or close to normal control values	
QUANTIPATH CHEMA with pathological control values.	
If required, a multiparametric, human based calibrator is available: AUTOCAL H	
Please contact Customer Care for further information.	
TEST PERFORMANCE	
Linearity	
the method is linear up to 300 mg/dl.	
If the value is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.	
Sensitivity/limit of detection (LOD)	
the limit of detection is 1 mg/dl.	
Interferences	
no interference was observed by the presence of:	
hemoglobin	≤ 500 mg/dl
bilirubin	≤ 44 mg/dl
lipids	≤ 600 mg/dl
Precision	
intra-assay (n=10)	mean (mg/dl) SD (mg/dl) CV%
sample 1	46.19 0.65 1.40
sample 2	140.89 2.72 1.90
inter-assay (n=20)	mean (mg/dl) SD (mg/dl) CV%
sample 1	42.77 1.91 4.50
sample 2	144.29 6.72 4.70
Methods comparison	
a comparison between Chema and a commercially available product gave the following results:	
Urea UV FL Chema = x	
Urea competitor = y	
n = 100	
$y = 0.9746x + 3.03 \text{ mg/dl} \quad r^2 = 0.986$	
WASTE DISPOSAL	
This product is made to be used in professional laboratories.	
P501: Dispose of contents according to national/international regulations.	