APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912	
TEST:	CREA
APP. CODE:	339
WAVELENGTH (Sec/Pri):	570 - 505
ASSAY:	2 POINT RATE <i>TIME: 10</i> <i>POINT: 19-23</i>
SAMPLE VOL:	NORMAL: 25 DECREASE: 20 INCREASE: 30
	R1 VOLUME: 125 R2 VOLUME: 0 R3 VOLUME: 125 R4 VOLUME: 0
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)
SD LIMIT:	0.250
DUPLICATE LIMIT:	3%
ST. 1 CONC:	0.00
EXPECTED VALUE:	0.6 - 1.2
UNIT:	mg/dl
INSTR. FACTOR (y=ax+b):	a=1 b= -0.3

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 869)	
TEST NAME:	CREA
SAMPLE:	Volume 25 µl Dilution 0 µl
REAGENTS:	R1 Volume 125 µl Dilution 0 µl R2 Volume 125 µl Dilution 0 µl
WAVELENGTH:	Pri. 520 Sec. 570
METHOD:	FIXED
REACTION SLOPE:	+
MEASURING POINT 1:	First 14 Last 18
MEASURING POINT 2:	First Last
REAGENT OD LIMIT:	First L -0.1 First H 0.5 Last L -0.1 Last H 0.5
DYNAMIC RANGE:	L 0.2 H 20
CORRELATION FACTOR:	A 1 B -0.3
UNIT:	mg/dl
CALIBRATION TYPE:	AB
FORMULA:	Y = AX + B

 Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

ITALIANO rev. 03/1/2017

CREATININA	
CR 2H400	4 x 50 + 4 x 50 ml
CR 6U448	4 x 56 + 4 x 56 ml

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INTENDED USE
Reagent for quantitative in vitro determination of creatinine in biological fluids.

PRINCIPLE OF THE METHOD
Creatinine reacts with picric acid in alkaline environment to form a color complex. Developing of this red color may be followed photometrically at 500-520 nm. The association on surfactant and sodium tetraborate keeps interferences at minimum.

KIT COMPONENTS

For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label.
Keep away from direct light sources.

CREA R1 2H400 4 x 50 ml (liquid) white cap
6U448 4 x 56 ml (liquid) white cap

CREA R2 2H400 4 x 50 ml (liquid) red cap
6U448 4 x 56 ml (liquid) red cap

Composition in the test: picric acid 14 mM, NaOH 0.18 M, sodium tetraborate 10 mM, surfactant.
Store all components at 15-25°C.

REAGENT PREPARATION

Use separate reagent ready to use.
Stability of working reagent: preferably within 30 days at 15-25°C, well capped and away from light sources.
Stability of unmixed reagents: up to expiration date on labels at 15-25°C;
Stability since first opening of vials of unmixed reagents: preferably within 60 days at 15-25°C.

PRECAUTIONS

CREA R1: Warning. Causes serious eye irritation (H319). Causes skin irritation (H315). Wear protective gloves. Eye protection (P280). IF ON SKIN: Wash with plenty of water (P302+P352). IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). If eye irritation persists: get medical advice (P337+P313).

CREA R2: It is not classified as hazardous.

SPECIMEN

Serum, plasma, Urine.
Creatinine is stable 24 hours at 2-8°C. Freeze samples for prolonged storage.
Dilute urine sample 1:100 with deionized water. It could be convenient a slight acidification of urine with HCl.

EXPECTED VALUES

Serum/plasma samples:
Men: 0.7 - 1.2 mg/dl (62 - 105 µmol/l)
Women: 0.6 - 1.1 mg/dl (53 - 97 µmol/l)

24h urine:
Men: 1000 - 2000 mg/24h (8.85 - 17.70 mmol/24h)
Women: 800 - 1800 mg/24h (7.08 - 15.93 mmol/24h)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

QUANTINORM CHEMA
with normal or close to normal control values
QUANTIPATH CHEMA
with pathological control values.
If required, a multiparametric, human based calibrator is available:
AUTOCAL H

Please contact Customer Care for further information.

TEST PERFORMANCE

Linearity
the method is linear up to 20 mg/dl.
If the value is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)
the limit of detection is 0.2 mg/dl.

Interferences
no interference was observed by the presence of:
hemoglobin ≤ 500 mg/dl
lipids ≤ 1250 mg/dl
Bilirubins give interference at low levels.

Precision
intra-assay (n=10) mean (mg/dl) SD (mg/dl) CV%
sample 1 1.25 0.03 2.60
sample 2 3.87 0.07 1.90

inter-assay (n=20) mean (mg/dl) SD (mg/dl) CV%
sample 1 1.31 0.04 2.90
sample 2 3.80 0.14 3.80

Methods comparison
a comparison between Chema and a commercially available product gave the following results:

Cretinine Chema = x
Creatinine competitor = y
n = 104
 $y = 0.982x - 0.081 \text{ mg/dl} \quad r^2 = 0.94$

WASTE DISPOSAL

This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.