APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN HITACHI 911/912	
TEST:	ASO
APP. CODE:	357
WAVELENGTH (Sec/Pri):	600
ASSAY:	2 POINT END <i>TIME: 10</i> <i>POINT: 16 - 31</i> <i>DILUENT: water</i>
SAMPLE VOL:	NORMAL: 3 DECREASE: 2 INCREASE: 5
	R1 VOLUME: 250 R2 VOLUME: 0 R3 VOLUME: 40 <i>DILUENT: 3</i> R4 VOLUME: 0
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)
SD LIMIT:	0.50
DUPLICATE LIMIT:	5%
ST. 1 CONC:	0.0
EXPECTED VALUE:	0.0 - 200
UNIT:	IU/ml
INSTR. FACTOR (y=ax+b):	a=1 b= 0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN OLYMPUS AU 400/480/600/640/680/2700 (Test code 892)	
TEST NAME:	ASO
SAMPLE:	Volume 3 µl Dilution 0 µl
REAGENTS:	R1 Volume 250 µl Dilution 0 µl R2 Volume 40 µl Dilution 0 µl
WAVELENGTH:	Pri. 600 Sec.
METHOD:	END
REACTION SLOPE:	+
MEASURING POINT 1:	First 0 Last 27
MEASURING POINT 2:	First 0 Last 10
REAGENT OD LIMIT:	First L -0.1 First H 0.5 Last L -0.1 Last H 0.5
DYNAMIC RANGE:	L 10 H 400
CORRELATION FACTOR:	A 1 B 0
UNIT:	IU/ml
CALIBRATION TYPE:	AB
FORMULA:	y = Ax + B

 Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

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INTENDED USE
Reagent for quantitative in vitro determination of antistreptolysin O in biological fluids.

PRINCIPLE OF THE METHOD
Antibodies of the sample selectively react with Streptolysin O coated to latex, thus producing particle agglutination. The produced turbidity is proportional to the concentration of ASO in the sample, and can be measured at the wavelength of 600 nm.

KIT COMPONENTS
For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label.
Keep away from direct light sources.

ASO R1 2H100 2 x 42 ml (liquid) white cap
6U140 2 x 54 ml (liquid) white cap

Composizione: Buffer pH 7.4, stabilizers and preservatives.

ASO R2 2H100 1 x 14 ml (liquid) red cap
6U140 2 x 9 ml (liquid) red cap

Composizione: suspension of latex particles coated with Streptolysin O, stabilizers and preservatives.

Store all components at 2-8°C.

REAGENT PREPARATION
Use separate reagents ready to use.
Kindly shake R2 vial before use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: use preferably within 60 days at 2-8°C.

PRECAUTIONS
Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

SPECIMEN
Serum, plasma.
Keep specimens away from direct light sources.
Samples are stable 2 days when stored at 2-8°C and 1 month at -20°C.

EXPECTED VALUES
Adults < 200 IU/ml

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

QUANTINORM CHEMA - MULTINORM CHEMA with normal or close to normal control values,
QUANTIPATH CHEMA - MULTIPATH CHEMA with pathological or close to pathological control values.

Please contact Customer Care for further information.

TEST PERFORMANCE
Linearity
the method is linear up to 400 IU/ml. If the value is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Hook effect
No Hook effect is observed with concentrations lower than 4100 IU/ml.

Sensitivity/limit of detection
The limit of detection is 10 IU/ml.

Interferences
No interference was observed by the presence of:
hemoglobin ≤ 1000 mg/dl
bilirubin ≤ 60 mg/dl
lipids ≤ 2500 mg/dl
rheumatoid factor ≤ 250 IU/ml

Precision
intra-assay (n=10) mean (IU/ml) SD (IU/ml) CV%
sample 1 74.7 1.56 2.09
sample 2 173.1 2.84 1.64

inter-assay (n=20) mean (IU/ml) SD (IU/ml) CV%
sample 1 75.7 3.44 4.54
sample 2 172.7 4.58 2.65

Methods comparison
A comparison between CHEMA and a commercially available product gave the following results:

ASO competitor = x
ASO FL CHEMA = y
n = 43
 $y = 0.93x + 2.88 \text{ IU/ml}$ $r^2 = 0.97$

WASTE DISPOSAL
This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.

ASO FL	
AS 2H100	2 x 42 + 1 x 14 ml
AS 6U140	2 x 54 + 2 x 9 ml

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