**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN
OLYMPUS AU 400/480/600/640/680/2700 (Test code 890)**

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INTENDED USE

Reagent for quantitative in vitro determination of glycated hemoglobin (HbA1c) in biological fluids.

PRINCIPLE OF THE METHOD

Glycated hemoglobin (HbA1c) is an important biochemical marker for assessing long-term glycemic control in patients with diabetes. HbA1c is being used not only to guide diabetes treatment but also to assess quality of care, and to predict risk for development and progression of diabetes complications.

KIT COMPONENTS

For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label.
Keep away from direct light sources at 2-8°C.

HC R1 1 x 45 ml (liquid) white cap

Composition: Latex suspension, stabilizers and preservatives.

HC R2 1 x 15 ml (liquid) red cap

Composition: Anti-HbA1c antibody, stabilizers and preservatives.

REAGENT PREPARATION

Use separate reagents ready to use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: use preferably within 60 days at 2-8°C.

PRECAUTIONS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

SPECIMEN

Whole blood.
Keep specimens away from direct light sources.
Specimens are stable 2 weeks when stored at 2-8°C. Do not freeze.

Preparation of test sample:

10 µl of whole blood have to be diluted with 500 µl of deionized water, in order to hemolyze the sample.

The above procedure does not apply to reconstituted calibrators and controls.

EXPECTED VALUES

Normal values (NGSP) 4.6-6.2%

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose the following human based control set is available:

HbA1c CONTROL SET

with normal and abnormal control values.

For calibration, use the following product:

HbA1c CALIBRATOR

Please contact Customer Care for further information.

TEST PERFORMANCE

Measuring range
Measure interval depends on the concentration of the highest standard used for calibration.
If such a limit value is exceeded, sample should not be diluted and tested. The result should be reported as "high".

Interferences

No interference was observed by the presence of:

free bilirubin	≤ 20 mg/dl
conjugated bilirubin	≤ 21 mg/dl
ascorbic acid	≤ 50 mg/dl
rheumatoid factor	≤ 1100 IU/ml

Precision

intra-assay (n=10)	mean (%)	SD (%)	CV%
sample 1	5.0	0.01	0.25
sample 2	12.8	0.02	0.17

inter-assay (n=5)	mean (%)	SD (%)	CV%
sample 1	5.0	0.02	0.32
sample 2	12.9	0.06	0.45

Methods comparison

A comparison between CHEMA and HPLC method gave the following results:

$$\begin{aligned} \text{HPLC method} &= x \\ \text{HbA1c CHEMA} &= y \\ n &= 142 \end{aligned}$$

$$y = 0.99x + 0.02\% \quad r^2 = 0.989$$

WASTE DISPOSAL

This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.

HbA1c	
HC 6U060	1 x 45 + 1 x 15 ml
UTILISATION	

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