**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ
HITACHI 911/912**

TEST: **CO2**
 APP. CODE: **378**
 WAVELENGTH (Sec/Pri): **505 - 415**
 ASSAY: **2 POINT RATE** TIME: **10**
 POINT: **3-18**
 SAMPLE VOL: NORMAL: **3**
 DECREASE: **2**
 INCREASE: **5**
 R1 VOLUME: **300**
 R2 VOLUME: **0**
 R3 VOLUME: **0**
 R4 VOLUME: **0**
 ABS LIMIT: **1500 - DEC**
 PROZONE LIMIT: **0 - UPPER**
 CALIB METHOD: **LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)**
 SD LIMIT: **0.250**
 DUPLICATE LIMIT: **3%**
 ST. 1 CONC: **0.0**
 EXPECTED VALUE: **0.0 - 30**
 UNIT: **mmol/l**
 INSTR. FACTOR (y=ax+b): a=1 b=0

**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ
OLYMPUS AU 400/480/600/640/680/2700 (Test code 861)**

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
BR 2H120 6 x 20 ml
 BR 6U168 3 x 56 ml

INTENDED USE


Reagent for quantitative in vitro determination of bicarbonate in biological fluids.

PRINCIPLE OF THE METHOD

Bicarbonate reacts with phosphoenolpyruvate (PEP), in the presence of phosphoenolpyruvate carboxylase (PEP-C), to form oxaloacetate and phosphate.

The oxaloacetate is then converted to malate by the action of malate dehydrogenase (MDH) and reduced nicotinamide adenine dinucleotide analog (NADH-analog).

The decrease in absorbance at 405 or 415 nm resulting from the oxidation of NADH analog is proportional to the amount of CO₂ in the sample.

KIT COMPONENTS

For in vitro diagnostic use only.
 The components of the kit are stable until expiration date on the label.
 Keep away from direct light sources.

CO2 R1 2H120 6 x 20 ml (liquid) white cap
6U168 3 x 56 ml (liquid) white cap

Composition: Buffer 0.1 M, PEP 10 mM, PEC-C > 100 U/l, MDH > 1000 U/l, NAD-analog 0.5 mM, stabilizer and preservative.

Store all components at 2-8°C.

REAGENT PREPARATION

Use separate reagent ready to use.
 Stability: up to expiration date on labels at 2-8°C.
 Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-.
 Caution: keep well refrigerated.

PRECAUTIONS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.
 Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

SPECIMEN

Serum, heparin plasma.
 Specimens should be protected from direct exposure to light. Samples stored at 2-8°C in the dark and in closed vials are stable up to 3 days and 1 month at -20°C.

EXPECTED VALUES

Newborn: 13 - 22 mmol/l - mEq/l
 Infant, child: 20 - 28 mmol/l - mEq/l
 Adults: 22 - 29 mmol/l - mEq/l
 Adults > 60yr: 23 - 31 mmol/l - mEq/l

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose a suitable human based control sera has to be used.
 Please contact Customer Care for further information.

TEST PERFORMANCE

Linearity
 the method is linear up to 50 mmol/l.
 If the limit value is exceeded, it is suggested to dilute sample 1+4 with distilled water and to repeat the test, multiplying the result by 5.

Sensitivity/limit of detection (LOD)
 the limit of detection is 0.7 mmol/l.

Interferences

no interference was observed by the presence of:
 lipids ≤ 1700 mg/dl
 bilirubin ≤ 70 mg/dl
 hemoglobin ≤ 1000 mg/dl

Precision

intra-assay (n=10) mean (mmol/l) SD (mmol/l) CV%
 sample 1 13.78 0.12 0.87
 sample 2 35.72 0.27 0.76

inter-assay (n=20) mean (mmol/l) SD (mmol/l) CV%
 sample 1 14.43 0.86 5.97
 sample 2 35.59 1.09 3.05

Methods comparison

a comparison between Chema and a commercially available product gave the following results:

$$\begin{aligned} &\text{Bicarbonato Chema} = x \\ &\text{Bicarbonato competitor} = y \\ &n = 86 \end{aligned}$$

$$y = 1.037x - 0.191 \text{ mmol/l} \quad r^2 = 0.992$$

WASTE DISPOSAL

This product is made to be used in professional laboratories.
 P501: Dispose of contents according to national/international regulations.