**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ
HITACHI 911/912**

TEST:	DBIL	
APP. CODE:	294	
WAVELENGTH (Sec/Pri):	700 - 546	
ASSAY:	2 POINT END	<i>TIME: 10 POINT: 16 - 31 DILUENT: water</i>
SAMPLE VOL:	NORMAL: 10 DECREASE: 8 INCREASE: 12	
	R1 VOLUME: 200 R2 VOLUME: 0 R3 VOLUME: 50 R4 VOLUME: 0	<i>DILUENT: 5 DILUENT: 5</i>
ABS LIMIT:	32000 - INC	
PROZONE LIMIT:	0 - UPPER	
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)	
SD LIMIT:	0.250	
DUPLICATE LIMIT:	3%	
ST. 1 CONC:	0.00	
EXPECTED VALUE:	0.00 - 0.20	
UNIT:	mg/dl	
INSTR. FACTOR (y=ax+b):	a=1 b=0	

**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ
OLYMPUS AU 400/480/600/640/680/2700 (Test code 851)**

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INTENDED USE
Reagent for quantitative in vitro determination of direct bilirubin in biological fluids.

PRINCIPLE OF THE METHOD
Conjugated (direct) bilirubin reacts with diazotized 2,4-dichloroaniline in acidic solution to produce an intensely coloured red diazo compound (520-560 nm). The intensity of color of this dye in solution is proportional to the concentration of direct bilirubin.

KIT COMPONENTS
For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label.
Keep away from direct light sources.

**BIL D R1 2H100 4 x 20 ml (liquid) white cap
6U420 6 x 56 ml (liquid) white cap**

Composition: sodium chloride 0.26 M, EDTA 0.1 mM.

**BIL D R2 2H100 2 x 10 ml (liquid) red cap
6U420 6 x 14 ml (liquid) red cap**

Composition: EDTA 0.1 mM, diazotized 2,4-dichloroaniline 0.1 mM, hydrochloric acid 0.18 M.

Store all components at 2-8°C.

REAGENT PREPARATION
Use separate reagent ready to use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-
Caution: keep well refrigerated.

PRECAUTIONS
Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

SPECIMEN
Serum, plasma.
Specimens should be protected from direct exposure to light. Samples stored at 2-8°C in the dark are stable up to 3 days and 1 month at -20°C.

EXPECTED VALUES
adults: ≤ 0.20 mg/dl (≤ 3.4 µmol/l)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:
QUANTINORM CHEMA with normal or close to normal control values
QUANTIPATH CHEMA with pathological control values.
If required, a multiparametric, human based calibrator is available:
AUTOCAL H

Please contact Customer Care for further information.

TEST PERFORMANCE

Linearity
the method is linear up to 13 mg/dl.
If the limit value is exceeded, it is suggested to dilute sample 1+9 with distilled water and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)
the limit of detection is 0.039 mg/dl.

Interferences
No interference was observed by the presence of:
hemoglobin ≤ 50 mg/dl
lipids ≤ 500 mg/dl
ascorbic acid ≤ 30 mg/dl

Precision
intra-assay (n=10) mean (mg/dl) SD (mg/dl) CV%
sample 1 0.719 0.003 0.44
sample 2 2.430 0.019 0.78

inter-assay (n=20) mean (mg/dl) SD (mg/dl) CV%
sample 1 0.735 0.039 5.31
sample 2 2.456 0.105 4.26

Methods comparison
a comparison between Chema and a commercially available product gave the following results:

Bilirubin direct Chema = x
Bilirubin direct competitor = y
n = 110

y = 0.911x - 0.049 mg/dl r² = 0.995

WASTE DISPOSAL

This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.

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