APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912	
TEST:	PH
APP. CODE:	344
WAVELENGTH (Sec/Pri):	700 - 340
ASSAY:	1-POINT <i>TIME: 10</i> <i>POINT: 31</i>
SAMPLE VOL:	NORMAL: 3 DECREASE: 2 INCREASE: 5
	R1 VOLUME: 300 R2 VOLUME: 0 R3 VOLUME: 0 R4 VOLUME: 0
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)
SD LIMIT:	0.250
DUPLICATE LIMIT:	3%
ST. 1 CONC:	0.00
EXPECTED VALUE:	2.0 - 7.0
UNIT:	mg/dl
INSTR. FACTOR (y=ax+b):	a=1 b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 874)	
TEST NAME:	PH
SAMPLE:	Volume 3 µl Dilution 0 µl
REAGENTS:	R1 Volume 300 µl Dilution 0 µl R2 Volume 0 µl Dilution 0 µl
WAVELENGTH:	Pri. 340 Sec. 700
METHOD:	END
REACTION SLOPE:	+
MEASURING POINT 1:	First 0 Last 17
MEASURING POINT 2:	First Last
REAGENT OD LIMIT:	First L -0.1 First H 0.5 Last L -0.1 Last H 0.5
DYNAMIC RANGE:	L 0.4 H 20
CORRELATION FACTOR:	A 1 B 0
UNIT:	mg/dl
CALIBRATION TYPE:	AB
FORMULA:	Y = AX + B

 Chema Diagnostica
Via Campana 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

FOSFORO UV	
PH 2H500	10 x 50 ml
PH 6U448	8 x 56 ml

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INTENDED USE
Reagent for quantitative in vitro determination of phosphorus in biological fluids.

PRINCIPLE OF THE METHOD
The phosphate ions react with ammonium molybdate to form a phosphomolybdate complex. The colorless phosphomolybdate complex can be measured directly by ultraviolet (UV) absorption at 340 nm. An acid pH is necessary for the formation of complexes.

KIT COMPONENTS
For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label.
Keep away from direct light sources.

PHOS R1 **2H500: 10 x 50 ml (liquid) white cap**
6U448: 8 x 56 ml (liquid) white cap
Composition: ammonium molybdate 0.4 mmol/l, sulphuric acid 0.21 mol/l, surfactant.
Store all components at 2-8°C.

REAGENT PREPARATION
Use separate reagent ready to use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: preferably within 60 days at 2-8°C away from light sources.
Caution: keep well refrigerated.

PRECAUTIONS
Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

Fosforus Chema = x
Fosforus competitor = y
n = 102

$$y = 1.005x - 0.109 \text{ mg/dl} \quad r^2 = 0.975$$

Serum is the preferred specimen. Although heparinized plasma is acceptable, levels of inorganic phosphate are about 0.2 to 0.3 mg/dl lower than in serum. Anticoagulants such as citrate, oxalate, and EDTA interfere with formation of the phosphomolybdate complex and should not be used.
Inorganic phosphate in whole blood specimens may either decrease or increase with time, depending on the type of specimen, temperature, and duration of storage. Levels in plasma or serum are increased by prolonged storage with cells at room temperature or 37 °C; it is important to promptly separate serum or plasma from erythrocytes. Hemolyzed specimens are unacceptable because erythrocytes contain high concentrations of organic phosphate esters, which can be hydrolyzed to inorganic phosphate during storage. Inorganic phosphate increases by 4 to 5 mg/dl per day in hemolyzed specimens stored at 4°C. Glucose phosphate, creatine phosphate, and other organic phosphates may also be hydrolyzed by assay conditions, resulting in overestimation of inorganic phosphate levels.
Phosphate is considered to be stable in serum that has been separated from the clot for days at 4°C and months when frozen.
Urine samples should be collected in 6 mol/l HCl, 20-30 ml for a 24 hours specimen, to avoid precipitation of phosphate complexes.
Dilute urine samples 1:20 with purified water before assay.

EXPECTED VALUES
serum/plasma (adults): 2.5 - 4.5 mg/dl (0.81 - 1.45 mmol/l)
serum/plasma (children): 4.0 - 7.0 mg/dl (1.29 - 2.26 mmol/l)
urine (nonrestricted diet): 0.4 - 1.3 g/24h (12.9 - 42.2 mmol/24h)

Each laboratory should establish appropriate reference intervals related to its population.
QUALITY CONTROL AND CALIBRATION
It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:
QUANTINORM CHEMA
with normal or close to normal control values
QUANTIPATH CHEMA
with pathological control values.
If required, a multiparametric, human based calibrator is available:
AUTOCAL H

Please contact Customer Care for further information.
TEST PERFORMANCE
Linearity
the method is linear up to 20 mg/dl.
If the limit value is exceeded, it is suggested to dilute sample 1+9 with distilled water and to repeat the test, multiplying the result by 10.
Sensitivity/limit of detection (LOD)
the limit of detection is 0.4 mg/dl.

Interferences
no interference was observed by the presence of:
bilirubin ≤ 25 mg/dl
hemoglobin ≤ 100 mg/dl
Hemolysis interferes.
Both positive and negative interferences with lipemic samples has been observed.

Precision

intra-assay (n=10)	mean (mg/dl)	SD (mg/dl)	CV%
sample 1	3.50	0.05	1.50
sample 2	5.87	0.11	1.90

inter-assay (n=20)	mean (mg/dl)	SD (mg/dl)	CV%
sample 1	3.41	0.08	2.40
sample 2	5.84	0.11	1.90

Methods comparison
a comparison between Chema and a commercially available product gave the following results:

phosphorus Chema = x
phosphorus competitor = y
n = 102

$$y = 1.005x - 0.109 \text{ mg/dl} \quad r^2 = 0.975$$

WASTE DISPOSAL
This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.

PHOSPHORE UV	
PH 2H500	10 x 50 ml
PH 6U448	8 x 56 ml

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