APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912	
TEST:	LDL
APP. CODE:	059
WAVELENGTH (Sec/Pri):	700 - 600
ASSAY:	2-POINT END <i>TIME: 10</i> <i>POINT: 16 - 31</i>
SAMPLE VOL:	NORMAL: 3 DECREASE: 2 INCREASE: 5
	R1 VOLUME: 270 R2 VOLUME: 0 R3 VOLUME: 90 R4 VOLUME: 0
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)
SD LIMIT:	0.250
DUPLICATE LIMIT:	3%
ST. 1 CONC:	0.0
EXPECTED VALUE:	70.0 - 220.0
UNIT:	mg/dl
INSTR. FACTOR (y=ax+b):	a=1 b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 867)	
TEST NAME:	LDL
SAMPLE:	Volume 3 µl Dilution 0 µl
REAGENTS:	R1 Volume 270 µl Dilution 0 µl R2 Volume 90 µl Dilution 0 µl
WAVELENGTH:	Pri. 600 Sec. 700
METHOD:	END
REACTION SLOPE:	+
MEASURING POINT 1:	First 0 Last 27
MEASURING POINT 2:	First 0 Last 10
REAGENT OD LIMIT:	First L -0.1 First H 0.5 Last L -0.1 Last H 0.5
DYNAMIC RANGE:	L 1.0 H 500
CORRELATION FACTOR:	A 1 B 0
UNIT:	mg/dl
CALIBRATION TYPE:	AB
FORMULA:	Y = AX + B

 Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

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INTENDED USE

Reagent for quantitative in vitro determination of LDL-cholesterol in biological fluids.

PRINCIPLE OF THE METHOD

When a sample is mixed with Reagent R1, the protecting reagent binds to LDL and protects LDL from enzyme reactions. Cholesterol esterase (CHE) and cholesterol oxidase (CO) react with non-LDL lipoproteins [chylomicrons (CM), very low density lipoproteins (VLDL) and HDL]. Hydrogen peroxide produced by the enzyme reactions with non-LDL cholesterol is decomposed by catalase in Reagent R1. When Reagent R2 is added, the protecting reagent is removed from LDL and catalase is inactivated. In this second process, CHE and CO react only with LDL-C. Hydrogen peroxide produced by the enzyme reactions with LDL-C yields a color complex upon oxidase condensation with N-(3-sulfopropyl)-3-methoxy-5-methylaniline (HMMP5) and 4-aminoantipyrine (4AA) in the presence of peroxidase (POD). By measuring the absorbance of the blue color complex produced at approximately 600 nm, the LDL-C concentration in the sample can be calculated when compared vs. the absorbance of the LDL-C Calibrator.

KIT COMPONENTS

For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label. Keep away from direct light sources.

LDL-C R1 2H080: 3 x 20 ml (liquid) white cap
6U120: 2 x 45 ml (liquid) white cap

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Store at 2-8°C. Do not freeze.

REAGENT PREPARATION

Use separate reagents ready to use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: 30 days at 2-8°C.

PRECAUTIONS

LDL-C R1: Warning. May cause an allergic skin reaction (H317). Wear protective gloves/protective clothing/eye protection/face protection (P280). IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower] (P303+P361+P353). If skin irritation or rash occurs: Get medical advice (P333+P313).

LDL-C R2: It is not classified as hazardous.

N-acetylcysteine (NAC), metamilzolo and acetaminophen may cause interference in the Trinder reaction.^(4,5) To avoid interference, the blood withdrawal should be performed before drug administration.

SPECIMEN

Serum, plasma heparinate.
Anticoagulants such as heparin, citrate, oxalate and EDTA do not have significant influences on the assay when they are used in their usual amounts.
Samples with triglyceride concentrations exceeding 1000 mg/dl should be diluted and reanalyzed.

Use fresh specimens. Do not use specimens that repeated freeze-thaw, because lipoproteins may be denatured.

EXPECTED VALUES

normal values: 76 - 218 mg/dl
NCEP ATP's Decision cut-off points for LDL-C:
desirable: < 130 mg/dl
borderline high risk for CHD 130 - 159 mg/dl
high risk for CHD ≥ 160 mg/dl

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:
QUANTINORM CHEMA with normal or close to normal control values
QUANTIPATH CHEMA with pathological control values.
If required, a multiparametric, human based calibrator is available:
AUTOCAL H

Please contact Customer Care for further information.

TEST PERFORMANCE

Linearity

the method is linear up to 500 mg/dl.
If the limit is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)
the limit of detection is 1 mg/dl.

Interferences

no interference was observed by the presence of:
hemoglobin ≤ 500 mg/dl
bilirubin (free) ≤ 50 mg/dl
bilirubin (conjugated) ≤ 40 mg/dl
ascorbic acid ≤ 50 mg/dl
lipid ≤ 2000 mg/dl

Precision

intra assay (n=21)	mean (mg/dl)	SD (mg/dl)	CV%
sample 1	85.6	1.46	1.71
sample 2	129.6	2.28	1.76

inter-assay (n=9)

mean (mg/dl)	SD (mg/dl)	CV%
sample 3	87.8	1.69
sample 4	129.6	1.99

Methods comparison

a comparison between LDL-direct FL and a CDC Reference method (beta-quantification) gave the following results:

LDL-direct FL Chema = x	Reference method = y
n = 25	
y = 1.0015x - 0.715 mg/dl	r ² = 0.986

WASTE DISPOSAL

This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.

LDL-direct FL	
DL 2H080	3 x 20 + 1 x 20 ml
DL 6U120	2 x 45 + 2 x 15 ml
UTILISATION	

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BIBLIOGRAFIA / BIBLIOGRAPHY / BIBLIOGRAPHIE / BIBLIOGRAFÍA / БИБЛИОГРАФИЯ

1) Burtis, C. A and Ashwood, E. R., Ed. Tietz Textbook of Clinical Chemistry, 2nd Ed., Saunders, Philadelphia, 1994.
2) NIH Publication No 95-3044, Recommendations on Lipoprotein Measurement (1995).
3) Japan Atherosclerosis Society: Guidelines for Prevention of Atherosclerotic Cardiovascular Diseases, 5-7 (2002).
4) N-acetylcysteine interference of Trinder-based assays. Genzen JR, Hunsaker JJ, Nelson LS, Faine BA, Krasowski MD. Clin Biochem. 2016 Jan;49(1-2):100-4
5) Drug interference in Trinder reaction. Wiewiorka O, Čermáková Z, Dastych M. Euromedlab 2017. ISSN 1437-4431

	dispositivo medico-diagnostico <i>in vitro</i> <i>in vitro</i> diagnostic medical device dispositif médical de diagnostic <i>in vitro</i> producto sanitario para diagnóstico <i>in vitro</i> <i>in vitro</i> diagnosticheeske medische instrumente
	numero di lotto batch code numéro de lot numero de lote лот выпуска
	numero di catalogo catalogue number réfêrence catalogue numero de catalogue номер по каталогу
	limite di temperatura temperature limit limite de température limite de temperatura диапазон температуры при хранении
	usare entro la data use-by date utiliser avant la date utilizar por fecha срок годности
	attenzione caution attention atención внимание