APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN <b>HITACHI 911/912</b>	
TEST:	<b>CREA-E</b>
APP. CODE:	<b>340</b>
WAVELENGTH (Sec/Pri):	<b>700 - 546</b>
ASSAY:	<b>2 POINT END</b> <span style="float: right;"><i>TIME: 10</i> <i>POINT: 16 - 31</i> <i>DILUENT: Water</i></span>
SAMPLE VOL:	NORMAL: <b>5</b> DECREASE: <b>3</b> INCREASE: <b>8</b>
	R1 VOLUME: <b>240</b> R2 VOLUME: <b>0</b> R3 VOLUME: <b>60</b> <span style="float: right;"><i>DILUENT: 5</i></span> R4 VOLUME: <b>0</b>
ABS LIMIT:	<b>32000 - INC</b>
PROZONE LIMIT:	<b>0 - UPPER</b>
CALIB METHOD:	<b>LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)</b>
SD LIMIT:	<b>0.250</b>
DUPLICATE LIMIT:	<b>3%</b>
ST. 1 CONC:	<b>0.00</b>
EXPECTED VALUE:	<b>0.5 - 1.2</b>
UNIT:	<b>mg/dl</b>
INSTR. FACTOR (y=ax+b):	<b>a=1 b=0</b>

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN <b>OLYMPUS AU 400/480/600/640/680/2700 (Test code 870)</b>	
TEST NAME:	<b>CREA-E</b>
SAMPLE:	Volume <b>5 µl</b> <span style="float: right;">Dilution <b>0 µl</b></span>
REAGENTS:	R1 Volume <b>240 µl</b> <span style="float: right;">Dilution <b>0 µl</b></span> R2 Volume <b>60 µl</b> <span style="float: right;">Dilution <b>0 µl</b></span>
WAVELENGTH:	Pri. <b>540</b> Sec. <b>700</b>
METHOD:	<b>END</b>
REACTION SLOPE:	<b>+</b>
MEASURING POINT 1:	First <b>0</b> Last <b>27</b>
MEASURING POINT 2:	First <b>0</b> Last <b>10</b>
REAGENT OD LIMIT:	First L <b>-0.1</b> First H <b>0.5</b> Last L <b>-0.1</b> Last H <b>0.5</b>
DYNAMIC RANGE:	L <b>0.05</b> H <b>50</b>
CORRELATION FACTOR:	A <b>1</b> B <b>0</b>
UNIT:	<b>mg/dl</b>
CALIBRATION TYPE:	<b>AB</b>
FORMULA:	<b>Y = AX + B</b>

 Chema Diagnostica  
Via Campania 2/4  
60030 Monsano (AN) - ITALY - EU  
phone +39 0731 605064  
fax +39 0731 605672  
e-mail: mail@chema.com  
website: http://www.chema.com

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**INTENDED USE**  
Reagent for quantitative in vitro determination of creatinine in biological fluids.

**PRINCIPLE OF THE METHOD**  
Through a series of enzymatic reactions, creatinine is converted in glycine, whilst endogenous components such as creatine and sarcosine are eliminated in the first step of the sequence. The formed hydrogen peroxide reacts with TOPS in the presence of peroxidase, to give a quinoneimine dye. The intensity of color, measured at 546 nm, is proportional to creatinine concentration in the sample.

**KIT COMPONENTS**  
**For in vitro diagnostic use only.**  
The components of the kit are stable until expiration date on the label at 2-8°C.  
Keep away from direct light sources.

**CREA-E R1** 2H250 4 x 50 ml (liquid) white cap  
6U280 4 x 56 ml (liquid) white cap  
**CREA-E R2** 2H250 2 x 25 ml (liquid) red cap  
6U280 4 x 14 ml (liquid) red cap

Composition in the test: Creatinase ≥ 10 kU/l, Creatininase ≥ 10 kU/l, Sarcosine Oxidase ≥ 1 kU/l, Peroxidase ≥ 5 kU/l, TOPS ≥ 3 mM, 4-aminoantipyrine ≥ 20 mg/l.

Store all components at 2-8°C.

**REAGENT PREPARATION**  
Use separate reagent ready to use.  
Stability: up to expiration date on labels at 2-8°C.  
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-.  
Caution: keep well refrigerated.

**PRECAUTIONS**  
Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.  
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

N-acetylcysteine (NAC), metamizole and acetaminophen may cause interference in the Trinder reaction.<sup>(1,2)</sup>  
To avoid interference, the blood withdrawal should be performed before drug administration.

**SPECIMEN**  
Serum, plasma. Urine.  
Creatinine is stable 24 hours at 2-8°C. Freeze samples for prolonged storage.  
Dilute urine sample 1:100 with deionized water.

**EXPECTED VALUES**  
Serum/plasma:  
Men: 0.67 - 1.17 mg/dl (59 - 104 µmol/l)  
Women: 0.51 - 0.95 mg/dl (45 - 84 µmol/l)

24h urine:  
Men: 1000 - 2000 mg/24h (8.85 - 17.70 mmol/24h)  
Women: 800 - 1800 mg/24h (7.08 - 15.93 mmol/24h)

Each laboratory should establish appropriate reference intervals related to its population.

**QUALITY CONTROL AND CALIBRATION**  
It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:  
**QUANTINORM CHEMA**  
with normal or close to normal control values.  
**QUANTIPATH CHEMA**  
with pathological control values.

If required, a multiparametric, human based calibrator is available:  
**AUTOCAL H**  
Please contact Customer Care for further information.

**TEST PERFORMANCE**  
**Linearity**  
The method is linear up to 50 mg/dl.  
If the value is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

**Sensitivity/limit of detection (LOD)**  
The limit of detection is 0.04 mg/dl.

**Interferences**  
No interference was observed by the presence of:  
hemoglobin ≤ 1000 mg/dl  
bilirubin ≤ 28 mg/dl  
lipids ≤ 1400 mg/dl  
ascorbic acid ≤ 50 mg/dl

**Precision**  
intra-assay (n=10) mean (mg/dl) SD (mg/dl) CV%  
sample 1 0.98 0.01 1.10  
sample 2 3.98 0.02 0.56  
inter-assay (n=20) mean (mg/dl) SD (mg/dl) CV%  
sample 1 0.97 0.03 3.28  
sample 2 3.97 0.11 2.90

**Methods comparison**  
A comparison between Chema and a commercially available product gave the following results with 99 samples:

$$\begin{aligned} \text{Creatinine competitor} &= x \\ \text{Creatinine Chema} &= y \\ y &= 1.004 x + 0.037 \text{ mg/dl} \quad r^2 = 0.998 \end{aligned}$$

**WASTE DISPOSAL**  
This product is made to be used in professional laboratories.  
P501: Dispose of contents according to national/international regulations.

CRÉATININE-E FL	
CE 2H250	4 x 50 + 2 x 25 ml
CE 6U280	4 x 56 + 4 x 14 ml

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