**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ
HITACHI 911/912**

TEST:	CHOL	
APP. CODE:	263	
WAVELENGTH (Sec/Pri):	700 - 505	
ASSAY:	1-POINT	TIME: 10 POINT: 20
SAMPLE VOL:	NORMAL: 3 DECREASE: 2 INCREASE: 5	
	R1 VOLUME: 300 R2 VOLUME: 0 R3 VOLUME: 0 R4 VOLUME: 0	
ABS LIMIT:	32000 - INC	
PROZONE LIMIT:	0 - UPPER	
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)	
SD LIMIT:	0.250	
DUPLICATE LIMIT:	3%	
ST. 1 CONC:	0.0	
EXPECTED VALUE:	100 - 200	
UNIT:	mg/dl	
INSTR. FACTOR (y=ax+b):	a=1	b=0

**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ
OLYMPUS AU 400/480/600/640/680/2700 (Test code 865)**

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INTENDED USE
Reagent for quantitative in vitro determination of cholesterol in biological fluids.

PRINCIPLE OF THE METHOD
All cholesterol esters present in plasma are hydrolyzed quantitatively into free cholesterol and fatty acids by cholesterol esterase. In the presence of oxygen, free cholesterol is then oxidized by cholesterol oxidase to cholesten-4-ene-3-one and H₂O₂. The H₂O₂ reacts with p-chlorophenol and 4-aminoantipyrine in the presence of peroxidase to form a quinoneimine dye. The intensity of color formed is proportional to the cholesterol concentration and can be measured photometrically between 480 and 520 nm.

KIT COMPONENTS
For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label.
Keep away from direct light sources.

**CHOL R1 2H500 10 x 50 ml (liquid) white cap
6U448 8 x 56 ml (liquid) white cap**

Composition: Good's buffer pH 7.20, sodium cholate 8 mM, CHE ≥ 400 U/l, CHOD ≥ 200 U/l, POD ≥ 500 U/l, 4-AAP 0.6 mM, 4-chlorophenol 2 mM.

Store all components at 2-8°C.

REAGENT PREPARATION
Use separate reagent ready to use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-.
Caution: keep well refrigerated.

CHOL R1: Warning. Causes serious eye irritation (H319). Causes skin irritation (H315). Wear protective gloves. Eye protection (P280). IF ON SKIN: Wash with plenty of water (P302+P352). IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). If eye irritation persists: get medical advice (P337+P313).

N-acetylcysteine (NAC), metamilozolo and acetaminophen may cause interference in the Trinder reaction.^(1,2)
To avoid interference, the blood withdrawal should be performed before drug administration.

SPECIMEN
Serum, plasma EDTA.
Sample is stable 3 days at 2-8°C and 1 month at -20°C.

EXPECTED VALUES

desirable:	140 - 200 mg/dl
borderline/high risk:	200 - 240 mg/dl
high risk:	> 240 mg/dl

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:
QUANTINORM CHEMA
with normal or close to normal control values
QUANTIPATH CHEMA
with pathological control values.
If required, a multiparametric, human based calibrator is available:
AUTOCAL H

Please contact Customer Care for further information.

TEST PERFORMANCE

Linearity
the method is linear up to 700 mg/dl.
If the limit value is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)
the limit of detection is 1 mg/dl.

Interferences
no interference was observed by the presence of:
hemoglobin ≤ 500 mg/dl
bilirubin ≤ 15 mg/dl
lipids ≤ 850 mg/dl

Precision

intra-assay (n=10)	mean (mg/dl)	SD (mg/dl)	CV%
sample 1	101.50	1.84	1.80
sample 2	176.20	2.74	1.60

inter-assay (n=20)	mean (mg/dl)	SD (mg/dl)	CV%
sample 1	100.99	2.11	2.10
sample 2	176.51	2.23	1.30

Methods comparison
a comparison between Chema and a commercially available product gave the following results:

Colesterolo FL Chema = x
Cholesterol competitor = y
n = 100
 $y = 0.979x - 1.71 \text{ mg/dl} \quad r^2 = 0.995$

WASTE DISPOSAL

This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.