APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN HITACHI 911/912	
TEST:	CREA-E
APP. CODE:	340
WAVELENGTH (Sec/Pri):	700 - 546
ASSAY:	2 POINT END <i>TIME: 10</i> <i>POINT: 16 - 31</i> <i>DILUENT: Water</i>
SAMPLE VOL:	NORMAL: 5 DECREASE: 3 INCREASE: 8
	R1 VOLUME: 240 R2 VOLUME: 0 R3 VOLUME: 60 <i>DILUENT: 5</i> R4 VOLUME: 0
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)
SD LIMIT:	0.250
DUPLICATE LIMIT:	3%
ST. 1 CONC:	0.00
EXPECTED VALUE:	0.5 - 1.2
UNIT:	mg/dl
INSTR. FACTOR (y=ax+b):	a=1 b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN OLYMPUS AU 400/480/600/640/680/2700 (Test code 870)	
TEST NAME:	CREA-E
SAMPLE:	Volume 5 µl Dilution 0 µl
REAGENTS:	R1 Volume 240 µl Dilution 0 µl R2 Volume 60 µl Dilution 0 µl
WAVELENGTH:	Pri. 540 Sec. 700
METHOD:	END
REACTION SLOPE:	+
MEASURING POINT 1:	First 0 Last 27
MEASURING POINT 2:	First 0 Last 10
REAGENT OD LIMIT:	First L -0.1 First H 0.5 Last L -0.1 Last H 0.5
DYNAMIC RANGE:	L 0.05 H 50
CORRELATION FACTOR:	A 1 B 0
UNIT:	mg/dl
CALIBRATION TYPE:	AB
FORMULA:	Y = AX + B

 Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

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INTENDED USE
Reagent for quantitative in vitro determination of creatinine in biological fluids.

PRINCIPLE OF THE METHOD
Through a series of enzymatic reactions, creatinine is converted in glycine, whilst endogenous components such as creatine and sarcosine are eliminated in the first step of the sequence. The formed hydrogen peroxide reacts with TOPS in the presence of peroxidase, to give a quinoneimine dye. The intensity of color, measured at 546 nm, is proportional to creatinine concentration in the sample.

KIT COMPONENTS
For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label at 2-8°C.
Keep away from direct light sources.

CREA-E R1 2H250 4 x 50 ml (liquid) white cap
6U280 4 x 56 ml (liquid) white cap
CREA-E R2 2H250 2 x 25 ml (liquid) red cap
6U280 4 x 14 ml (liquid) red cap

Composition in the test: Creatinase ≥ 10 kU/l, Creatininase ≥ 10 kU/l, Sarcosine Oxidase ≥ 1 kU/l, Peroxidase ≥ 5 kU/l, TOPS ≥ 3 mM, 4-aminoantipyrine ≥ 20 mg/l.

Store all components at 2-8°C.

REAGENT PREPARATION
Use separate reagent ready to use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-.
Caution: keep well refrigerated.

PRECAUTIONS
Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

N-acetylcysteine (NAC), metamizole and acetaminophen may cause interference in the Trinder reaction.^(1,2)
To avoid interference, the blood withdrawal should be performed before drug administration.

SPECIMEN
Serum, plasma. Urine.
Creatinine is stable 24 hours at 2-8°C. Freeze samples for prolonged storage.
Dilute urine sample 1:100 with deionized water.

EXPECTED VALUES
Serum/plasma:
Men: 0.67 - 1.17 mg/dl (59 - 104 µmol/l)
Women: 0.51 - 0.95 mg/dl (45 - 84 µmol/l)
24h urine:
Men: 1000 - 2000 mg/24h (8.85 - 17.70 mmol/24h)
Women: 800 - 1800 mg/24h (7.08 - 15.93 mmol/24h)
Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION
It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:
QUANTINORM CHEMA
with normal or close to normal control values.
QUANTIPATH CHEMA
with pathological control values.

If required, a multiparametric, human based calibrator is available:
AUTOCAL H
Please contact Customer Care for further information.

TEST PERFORMANCE
Linearity
The method is linear up to 50 mg/dl.
If the value is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)
The limit of detection is 0.04 mg/dl.

Interferences
No interference was observed by the presence of:
hemoglobin ≤ 1000 mg/dl
bilirubin ≤ 28 mg/dl
lipids ≤ 1400 mg/dl
ascorbic acid ≤ 50 mg/dl

Precision
intra-assay (n=10) mean (mg/dl) SD (mg/dl) CV%
sample 1 0.98 0.01 1.10
sample 2 3.98 0.02 0.56
inter-assay (n=20) mean (mg/dl) SD (mg/dl) CV%
sample 1 0.97 0.03 3.28
sample 2 3.97 0.11 2.90

Methods comparison
A comparison between Chema and a commercially available product gave the following results with 99 samples:

$$\begin{aligned} \text{Creatinine competitor} &= x \\ \text{Creatinine Chema} &= y \\ y &= 1.004 x + 0.037 \text{ mg/dl} \quad r^2 = 0.998 \end{aligned}$$

WASTE DISPOSAL
This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.

CRÉATININE-E FL	
CE 2H250	4 x 50 + 2 x 25 ml
CE 6U280	4 x 56 + 4 x 14 ml

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