APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ <b>HITACHI 911/912</b>	
TEST:	<b>HDL</b>
APP. CODE:	<b>035</b>
WAVELENGTH (Sec/Pri):	<b>700 - 600</b>
ASSAY:	<b>2-POINT END</b> <i>TIME: 10</i> <i>POINT: 3 - 31</i>
SAMPLE VOL:	NORMAL: <b>3</b> DECREASE: <b>2</b> INCREASE: <b>5</b>
	R1 VOLUME: <b>270</b> R2 VOLUME: <b>90</b> R3 VOLUME: <b>0</b> R4 VOLUME: <b>0</b>
ABS LIMIT:	<b>32000 - INC</b>
PROZONE LIMIT:	<b>0 - UPPER</b>
CALIB METHOD:	<b>LINEAR</b> (POINT: 2 - SPAN: 2 - WEIGHT: 0)
SD LIMIT:	<b>0.250</b>
DUPLICATE LIMIT:	<b>3%</b>
ST. 1 CONC:	<b>0.0</b>
EXPECTED VALUE:	<b>30.0 - 90.0</b>
UNIT:	<b>mg/dl</b>
INSTR. FACTOR (y=ax+b):	a=1    b= 0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ <b>OLYMPUS AU 400/480/600/640/680/2700 (Test code 866)</b>	
TEST NAME:	<b>HDL</b>
SAMPLE:	Volume <b>3</b> µl      Dilution <b>0</b> µl
REAGENTS:	R1 Volume <b>270</b> µl      Dilution <b>0</b> µl R2 Volume <b>90</b> µl      Dilution <b>0</b> µl
WAVELENGTH:	Pri. <b>600</b> Sec. <b>700</b>
METHOD:	<b>END</b>
REACTION SLOPE:	<b>+</b>
MEASURING POINT 1:	First <b>0</b> Last <b>27</b>
MEASURING POINT 2:	First <b>0</b> Last <b>10</b>
REAGENT OD LIMIT:	First L <b>-0.1</b> First H <b>0.5</b> Last L <b>-0.1</b> Last H <b>0.5</b>
DYNAMIC RANGE:	L <b>1.0</b> H <b>220</b>
CORRELATION FACTOR:	A <b>1</b> B <b>0</b>
UNIT:	<b>mg/dl</b>
CALIBRATION TYPE:	<b>AB</b>
FORMULA:	<b>Y = AX + B</b>

 Chema Diagnostica  
Via Campania 2/4  
60030 Monsano (AN) - ITALY - EU  
phone +39 0731 605064  
fax +39 0731 605672  
e-mail: mail@chema.com  
website: http://www.chema.com

ITALIANO rev. 09/05/2019

HDL-direct FL	
HD 2H160	6 x 20 + 2 x 20 ml
HD 6U240	4 x 45 + 4 x 15 ml

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**INTENDED USE**  
Reagent for quantitative in vitro determination of HDL-cholesterol in biological fluids.

**PRINCIPLE OF THE METHOD**  
Anti human β-lipoprotein antibody in reagent R1 binds to lipoproteins (LDL, VLDL, and chylomicrons) other than HDL. The antigen-antibody complexes formed block enzyme reactions when reagent R2 is added. Cholesterol esterase and cholesterol oxidase in reagent R2 react only with HDL-C. Hydrogen peroxide produced by the enzyme reactions with HDL-C yields a blue color complex upon oxidative condensation of F-DAOS [N-ethyl-N-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxy-4-fluoroaniline, sodium salt] and 4-aminoantipyrine in the presence of peroxidase. By measuring the absorbance of the blue color complex produced, at the optimum wavelength of 593 nm, the HDL-C concentration in the sample can be calculated when compared with the absorbance of the HDL-C Calibrator.

**KIT COMPONENTS**  
**For in vitro diagnostic use only.**  
The components of the kit are stable until expiration date on the label.  
Keep away from direct light sources.

**HDL-C R1**    **2H160: 6 x 20 ml (liquid) white cap**  
**6U240: 4 x 45 ml (liquid) white cap**

Composizione: Good's buffer pH 7.0 30 mmol/l, 4-aminoantipirina 0.9 mmol/l, POD 2400 U/l, ascorbate ossidasi 2700 U/l, and anti human-lipoprotein antibody, blend of 5-chloro-2-methyl-2-H-isothiazol-3-one [EC No 247-500-7] and 2-methyl-2-H-isothiazol-3-one [EC No 220-239-6] (3:1) in concentration 0.0015-0.06%.

**HDL-C R2**    **2H160: 2 x 20 ml (liquid) red cap**  
**6U240: 4 x 15 ml (liquid) red cap**

Composizione: Good's buffer pH 7.0 30 mmol/l, colesterolo esterasi 4000 U/l, colesterolo ossidasi 20000 U/l, and F-DAOS 0.8 mmol/l.

Store all components at 2-8°C.

**REAGENT PREPARATION**  
Use separate reagent ready to use.  
Stability: up to expiration date on labels at 2-8°C.  
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-.  
Caution: keep well refrigerated.

**PRECAUTIONS**  
**HDL-C R1: Warning.** May cause an allergic skin reaction (H317). Avoid breathing vapours (P261). IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water (P303+P361+P353). If skin irritation or rash occurs: Get medical advice (P333+P313).

**HDL-C R2:** It is not classified as hazardous.

N-acetylcysteine (NAC), metamizole and acetaminophen may cause interference in the Trinder reaction.<sup>(1,2)</sup> To avoid interference, the blood withdrawal should be performed before drug administration.

**SPECIMEN**  
Use serum as a specimen. It is recommended to measure HDL-C immediately after collection. Ascorbic acid, bilirubin, and hemoglobin do not have a significant effect on the measurement.

**EXPECTED VALUES**  
Adult male: 35.3 - 79.5 mg/dl  
Adult female: 42.0 - 88.0 mg/dl

Each laboratory should establish appropriate reference intervals related to its population.

**QUALITY CONTROL AND CALIBRATION**  
It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:  
**QUANTINORM CHEMA**  
with normal or close to normal control values  
**QUANTIPATH CHEMA**  
with pathological control values.  
If required, a multiparametric, human based calibrator is available:  
**AUTOCAL H**

Please contact Customer Care for further information.

**TEST PERFORMANCE**  
**Linearity**  
the method is linear up to 220 mg/dl.  
If the limit is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

**Sensitivity/limit of detection (LOD)**  
the limit of detection is 1 mg/dl.

**Interferences**  
no interference was observed by the presence of:  
hemoglobin ≤ 500 mg/dl  
bilirubin (free) ≤ 50 mg/dl  
bilirubin (conjugated) ≤ 40 mg/dl  
ascorbic acid ≤ 50 mg/dl

**Precision**  
intra assay (n=10)    mean (mg/dl)    SD (mg/dl)    CV%  
sample 1    32.1    0.18    0.55  
sample 2    88.9    0.61    0.68

**Comparison**  
a comparison between HDL-direct FL and other methods has shown the following results:

HDL-C Chema = x  
HDL-C competitor = y  
n = 50  
  
y = 0.96x + 2.5 mg/dl    r<sup>2</sup> = 0.998

**WASTE DISPOSAL**  
This product is made to be used in professional laboratories.  
P501: Dispose of contents according to national/international regulations.

HDL-direct FL	
HD 2H160	6 x 20 + 2 x 20 ml
HD 6U240	4 x 45 + 4 x 15 ml

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#### BIBLIOGRAFIA / BIBLIOGRAPHY / BIBLIOGRAPHIE / BIBLIOGRAFÍA / БИБЛИОГРАФИЯ

1) N-acetylcysteine interference of Trinder-based assays. Genzen JR, Hunsaker JJ, Nelson LS, Faine BA, Krasowski MD. Clin Biochem. 2016 Jan;49(1-2):100-4
2) Drug interference in Trinder reaction. Wiewiorka O, Čermáková Z, Dastych M. Euromedlab 2017. ISSN 1437-4431
3) Rifai, N., Warnick, G.R. Ed. Laboratory Measurement of Lipids, Lipoproteins and Apolipoproteins AACCC Press. Washington, DC, USA, 1994
4) Burtis, C. A and Ashwood, E. R., Ed. Tietz Textbook of Clinical Chemistry, 2nd Ed., Saunders, Philadelphia, 1994.
5) Gordon, T., Castelli, W.P., Hjortland, M.C., et al., Am. J. Med 62,707 - 714, (1977)

<b>IVD</b>	dispositivo medico-diagnostico <i>in vitro</i> <i>in vitro</i> diagnostic medical device dispositif médical de diagnostic <i>in vitro</i> producto sanitario para diagnóstico <i>in vitro</i> <i>in vitro</i> диагностические медицинские устройства
<b>LOT</b>	numero di lotto batch code numéro de lot número de lote лот выпуска
<b>REF</b>	numero di catalogo catalogue number réfêrence catalogue número de catálogo номер по каталогу
	limite di temperatura temperature limit limite de température limite de temperatura диапазон температуры при хранении
	usare entro la data use-by date utiliser avant la date utilizar por fecha срок годности
	attenzione caution attention atención внимание
	consultare le istruzioni d'uso consult instructions for use consulter les instructions d'utilisation consultar las instrucciones de uso смотреть рабочие инструкции