MICROALBUMIN CONTROL SET

2 x 2 ml

REF MA S004 CH

LOT -----

 -----

| componente component composant componente | metodo method méthode método | valore value valeur valor | intervallo range intervalle intervalo | unità unit unidad |
|---|---|------------------------------------|--|-------------------------|
| LEV 1 Microalbumina Microalbumin Microalbumine Microalbúmina | metodo immunoturbidimetrico immunoturbidimetric method méthode immunoturbidimétrique método inmunoturbidimétrico | ---- | ---- | mg/l |
| LEV 2 Microalbumina Microalbumin Microalbumine Microalbúmina | metodo immunoturbidimetrico immunoturbidimetric method méthode immunoturbidimétrique método inmunoturbidimétrico | ---- | ---- | mg/l |

Italiano

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Usò previsto

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CONTROL SET

Intended use

MICROALBUMIN CONTROL SET is composed of a sera set of quality control of methods for determination of Microalbumin in urine, through immunoturbidimetric technique. The solutions are used to monitor accuracy and/or precision both in automated analyzers and in manual use.

Package components

Controls (liquid) 2 x 2 ml

Composition

MICROALBUMIN CONTROL SET is a liquid. It is a dilution of defibrinated human plasma with phosphate buffered saline. It also contains stabilizers and sodium azide (< 0,1%) as a biocide.

Assigned values

The value for microalbumin is based on RPPHS/ERM-DA470 from IFCC.

Storage and stability

Store at 2-8°C.

The controls are ready for use; after vial opening, sera are stable for 6 weeks if stored tightly closed at 2-8°C. Do not freeze.

Test procedure

Use supplied materials as provided.

Additional materials required: automatic micropipette, general laboratory equipment.

Precautions and warnings

For in vitro diagnostic use only.

Exercise the normal precautions required for handling all laboratory reagents. This product has been tested and found free from HbsAg and antibodies to HCV and HIV. However, because no test method can offer complete assurance that products derived from human sources will not transmit infectious agents, it is recommended that this product is handled with the same precautions used for patient specimens.

Disposal

This product is made to be used in professional laboratories. Please consult local regulations for a correct waste disposal.

P501: Dispose of contents according to national/international regulations.

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