APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN HITACHI 911/912	
TEST:	IGG
APP. CODE:	354
WAVELENGTH (Sec/Pri):	600
ASSAY:	2 POINT END <i>TIME: 10</i> <i>POINT: 16 - 31</i> <i>DILUENT: water</i>
SAMPLE VOL:	NORMAL: 3 DECREASE: 2 INCREASE: 5
	R1 VOLUME: 240 R2 VOLUME: 0 R3 VOLUME: 60 <i>DILUENT: 5</i> R4 VOLUME: 0
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LOGIT-LOG (5P) <i>(POINT: 6 - SPAN: 6 - WEIGHT: 0)</i>
SD LIMIT:	0.50
DUPLICATE LIMIT:	10%
ST. 1 CONC:	0.0
EXPECTED VALUE:	7.0 - 16.0
UNIT:	g/l
INSTR. FACTOR (y=ax+b):	a=1 b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN OLYMPUS AU 400/480/600/640/680/2700 (Test code 887)	
TEST NAME:	IGG
SAMPLE:	Volume 3 µl Dilution 0 µl
REAGENTS:	R1 Volume 240 µl Dilution 0 µl R2 Volume 60 µl Dilution 0 µl
WAVELENGTH:	Pri. 600 Sec.
METHOD:	END
REACTION SLOPE:	+
MEASURING POINT 1:	First 0 Last 27
MEASURING POINT 2:	First 0 Last 10
REAGENT OD LIMIT:	First L -0.1 First H 0.5 Last L -0.1 Last H 0.5
DYNAMIC RANGE:	L 0.01 H 30.0
CORRELATION FACTOR:	A 1 B 0
UNIT:	g/l
CALIBRATION TYPE:	5AB
FORMULA:	POLYGONAL

 Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

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INTENDED USE
Reagent for quantitative in vitro determination of IgG in biological fluids.

PRINCIPLE OF THE METHOD
Immunoglobulins G (IgG) selectively react with an anti-IgG antibody and form an immunocomplex. The produced turbidity is proportional to the concentration of IgG in the sample, and can be measured at the wavelength of 600 nm.

KIT COMPONENTS
For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label.
Keep away from direct light sources at 2-8°C..

IGG R1 **2H100 2 x 40 ml (liquid) white cap**
6U140 2 x 56 ml (liquid) white cap

Composizione: Buffer pH 7.50, PEG ≥ 2%, stabilizers and preservatives.

IGG R2 **2H100 2 x 10 ml (liquid) red cap**
6U140 2 x 14 ml (liquid) red cap

Composizione: Anti-human IgG antibody ≥ 2%, stabilizers and preservatives.

REAGENT PREPARATION
Use separate reagent ready to use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-.
Caution: keep well refrigerated.

PRECAUTIONS
Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

SPECIMEN
Serum, plasma.
Keep specimens away from direct light sources.
Samples are stable 7 days when stored at 2-8°C and 1 month at -20°C.

EXPECTED VALUES
Newborns 7.0-14.8 g/l
Adults 7.0-16.0 g/l

Each laboratory should establish appropriate reference intervals related to its population

QUALITY CONTROL AND CALIBRATION
It is suggested to perform an internal quality control. For this purpose the following human based control serum is available:
QUANTINORM CHEMA
with normal or close to normal control values.
If required, a multiparametric, human based calibrator is available:
REFERENCE P MULTICALIBRATOR
Please contact Customer Care for further information.

TEST PERFORMANCE
Measure interval
Measure interval depends on the concentration of the highest standard used for calibration.
If such a limit value is exceeded, it is suggested to dilute sample 1+4 with distilled water and to repeat the test, multiplying the result by 5.

Hook effect
No Hook effect is observed with concentrations lower than 180 g/l.

Sensitivity/limit of detection
The limit of detection is 0.02 g/l.

Interferences
No interference was observed by the presence of:
hemoglobin ≤ 1000 mg/dl
bilirubin ≤ 40 mg/dl
lipids ≤ 950 mg/dl
rheumatoid factor ≤ 630 IU/ml

intra-assay (n=10) mean (g/l) SD (g/l) CV%
sample 1 12.45 0.20 1.58
sample 2 16.01 0.19 1.16

inter-assay (n=20) mean (g/l) SD (g/l) CV%
sample 1 12.39 0.28 2.28
sample 2 15.98 0.28 1.77

Methods comparison
A comparison between CHEMA and a commercially available product gave the following results:

$$\begin{aligned} \text{IgG competitor} &= x \\ \text{IgG FL CHEMA} &= y \\ n &= 72 \\ y &= 0.985x - 0.48 \text{ g/l} \quad r^2 = 0.99 \end{aligned}$$

WASTE DISPOSAL
This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.