APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN HITACHI 911/912	
TEST:	IGM
APP. CODE:	355
WAVELENGTH (Sec/Pri):	340
ASSAY:	2 POINT END <i>TIME: 10</i> <i>POINT: 16 - 31</i> <i>DILUENT: water</i>
SAMPLE VOL:	NORMAL: 3 DECREASE: 2 INCREASE: 5
	R1 VOLUME: 240 R2 VOLUME: 0 R3 VOLUME: 60 <i>DILUENT: 5</i> R4 VOLUME: 0
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LOGIT-LOG (5P) <i>(POINT: 6 - SPAN: 6 - WEIGHT: 0)</i>
SD LIMIT:	0.50
DUPLICATE LIMIT:	10%
ST. 1 CONC:	0.0
EXPECTED VALUE:	0.0 - 2.3
UNIT:	g/l
INSTR. FACTOR (y=ax+b):	a=1 b= 0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN OLYMPUS AU 400/480/600/640/680/2700 (Test code 888)	
TEST NAME:	IGM
SAMPLE:	Volume 3 µl Dilution 0 µl
REAGENTS:	R1 Volume 240 µl Dilution 0 µl R2 Volume 60 µl Dilution 0 µl
WAVELENGTH:	Pri. 340 Sec.
METHOD:	END
REACTION SLOPE:	+
MEASURING POINT 1:	First 0 Last 27
MEASURING POINT 2:	First 0 Last 10
REAGENT OD LIMIT:	First L -0.1 First H 0.5 Last L -0.1 Last H 0.5
DYNAMIC RANGE:	L 0.02 H 4.0
CORRELATION FACTOR:	A 1 B 0
UNIT:	g/l
CALIBRATION TYPE:	5AB
FORMULA:	POLYGONAL

 Chema Diagnostica
 Via Campania 2/4
 60030 Monsano (AN) - ITALY - EU
 phone +39 0731 605064
 fax +39 0731 605672
 e-mail: mail@chema.com
 website: http://www.chema.com

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INTENDED USE
Reagent for quantitative in vitro determination of IgM in biological fluids.

PRINCIPLE OF THE METHOD
Immunoglobulins M (IgM) selectively react with an anti-IgM antibody and form an immunocomplex. The produced turbidity is proportional to the concentration of IgM in the sample and can be measured at the wavelength of 340 nm.

KIT COMPONENTS
For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label.
Keep away from direct light sources at 2-8°C.

IGM R1 **2H100 2 x 40 ml (liquid) white cap**
6U140 2 x 56 ml (liquid) white cap

Composizione: Buffer pH 7.50, PEG ≥ 2%, stabilizers and preservatives.

IGM R2 **2H100 2 x 10 ml (liquid) red cap**
6U140 2 x 14 ml (liquid) red cap

Composizione: Anti-human IgM antibody ≥ 2%, stabilizers and preservatives.

REAGENT PREPARATION
Use separate reagent ready to use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-.
Caution: keep well refrigerated.

PRECAUTIONS
IGM R1: Danger. Causes serious eye damage (H318).
Wear protective gloves. Eye protection (P280).
 **IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). Immediately call a doctor (P310).

IGM R2: It is not classified as hazardous.

SPECIMEN
Serum, plasma.
Keep specimens away from direct light sources.
Samples are stable 7 days when stored at 2-8°C and 1 month at -20°C.

EXPECTED VALUES	
Newborns	0.05-0.3 g/l
Adults	0.4-2.3 g/l

Each laboratory should establish appropriate reference intervals related to its population

QUALITY CONTROL AND CALIBRATION
It is suggested to perform an internal quality control. For this purpose the following human based control serum is available:

QUANTINORM CHEMA
with normal or close to normal control values.
If required, a multiparametric, human based calibrator is available:
REFERENCE P MULTICALIBRATOR

Please contact Customer Care for further information.

TEST PERFORMANCE
Measure interval
Measure interval depends on the concentration of the highest standard used for calibration.
If such a limit value is exceeded, it is suggested to dilute sample 1+4 with distilled water and to repeat the test, multiplying the result by 5.

Hook effect
No Hook effect is observed with concentrations lower than 30 g/l.

Sensitivity/limit of detection
The limit of detection is 0.016 g/l.

Interferences
No interference was observed by the presence of:
hemoglobin ≤ 1000 mg/dl
bilirubin ≤ 45 mg/dl
lipids ≤ 770 mg/dl
rheumatoid factor ≤ 630 IU/ml

Precision	mean (g/l)	SD (g/l)	CV%
intra-assay (n=10)			
sample 1	0.70	0.005	0.69
sample 2	1.40	0.009	0.66
inter-assay (n=20)			
sample 1	0.70	0.024	3.38
sample 2	1.40	0.061	4.32

Methods comparison
A comparison between CHEMA and a commercially available product gave the following results:

$$y = 1.186x - 0.058 \text{ g/l} \quad r^2 = 0.99$$

WASTE DISPOSAL
This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.

IgM FL	
GM 2H100	2 x 40 + 2 x 10 ml
GM 6U140	2 x 56 + 2 x 14 ml

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