APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN HITACHI 911/912	
TEST:	<b>TRF</b>
APP. CODE:	<b>356</b>
WAVELENGTH (Sec/Pri):	<b>505</b>
ASSAY:	<b>2 POINT END</b> <i>TIME: 10</i> <i>POINT: 16 - 31</i> <i>DILUENT: water</i>
SAMPLE VOL:	NORMAL: <b>2</b> DECREASE: <b>2</b> INCREASE: <b>5</b>
	R1 VOLUME: <b>240</b> R2 VOLUME: <b>0</b> R3 VOLUME: <b>60</b> <i>DILUENT: 5</i> R4 VOLUME: <b>0</b>
ABS LIMIT:	<b>32000 - INC</b>
PROZONE LIMIT:	<b>0 - UPPER</b>
CALIB METHOD:	<b>LOGIT-LOG (5P)</b> <i>(POINT: 6 - SPAN: 6 - WEIGHT: 0)</i>
SD LIMIT:	<b>0.50</b>
DUPLICATE LIMIT:	<b>10%</b>
ST. 1 CONC:	<b>0.0</b>
EXPECTED VALUE:	<b>1.0 - 4.0</b>
UNIT:	<b>g/l</b>
INSTR. FACTOR (y=ax+b):	a=1    b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN OLYMPUS AU 400/480/600/640/680/2700 (Test code 889)	
TEST NAME:	<b>TRF</b>
SAMPLE:	Volume <b>2 µl</b> Dilution <b>0 µl</b>
REAGENTS:	R1 Volume <b>240 µl</b> Dilution <b>0 µl</b> R2 Volume <b>60 µl</b> Dilution <b>0 µl</b>
WAVELENGTH:	Pri. <b>520</b> Sec.
METHOD:	<b>END</b>
REACTION SLOPE:	<b>+</b>
MEASURING POINT 1:	First <b>0</b> Last <b>27</b>
MEASURING POINT 2:	First <b>0</b> Last <b>10</b>
REAGENT OD LIMIT:	First L <b>-0.1</b> First H <b>0.5</b> Last L <b>-0.1</b> Last H <b>0.5</b>
DYNAMIC RANGE:	L <b>0.01</b> H <b>9.0</b>
CORRELATION FACTOR:	A <b>1</b> B <b>0</b>
UNIT:	<b>g/l</b>
CALIBRATION TYPE:	<b>5AB</b>
FORMULA:	<b>POLYGONAL</b>

 Chema Diagnostica  
Via Campania 2/4  
60030 Monsano (AN) - ITALY - EU  
phone +39 0731 605064  
fax +39 0731 605672  
e-mail: mail@chema.com  
website: http://www.chema.com

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**INTENDED USE**  
Reagent for quantitative in vitro determination of transferrin in biological fluids.

**PRINCIPLE OF THE METHOD**  
Transferrin selectively reacts with an anti-transferrin antibody and forms an immunocomplex. The produced turbidity is proportional to the concentration of transferrin in the sample, and can be measured at the wavelength of 505 nm.

**KIT COMPONENTS**  
**For in vitro diagnostic use only.**  
The components of the kit are stable until expiration date on the label.  
Keep away from direct light sources.

**TRF R1**    **2H100: 2 x 40 ml (liquid) white cap**  
**6U140: 2 x 56 ml (liquid) white cap**

Composizione: Buffer pH 7.00, PEG ≥ 2%, stabilizers and preservatives.

**TRF R2**    **2H100: 2 x 10 ml (liquid) red cap**  
**6U140: 2 x 14 ml (liquid) red cap**

Composizione: Anti-human transferrin antibody ≥ 2%, stabilizers and preservatives.

Store all components at 2-8°C.

**REAGENT PREPARATION**  
Use separate reagent ready to use.  
Stability: up to expiration date on labels at 2-8°C.  
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-.  
Caution: keep well refrigerated.

**PRECAUTIONS**  
Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.  
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

**SPECIMEN**  
Serum, plasma.  
Keep specimens away from direct light sources.  
Samples are stable 7 days when stored at 2-8°C and 1 month at -20°C.

EXPECTED VALUES	
Newborns	1.2-2.5 g/l
Adults	2.0-3.6 g/l

Each laboratory should establish appropriate reference intervals related to its population

**QUALITY CONTROL AND CALIBRATION**

It is suggested to perform an internal quality control. For this purpose the following human based control serum is available:

**QUANTINORM CHEMA**  
with normal or close to normal control values.  
If required, a multiparametric, human based calibrator is available:

**REFERENCE P MULTICALIBRATOR**  
Please contact Customer Care for further information.

**TEST PERFORMANCE**

**Measure interval**  
Measure interval depends on the concentration of the highest standard used for calibration.  
If such a limit value is exceeded, it is suggested to dilute sample 1+4 with distilled water and to repeat the test, multiplying the result by 5.

**Hook effect**  
No Hook effect is observed with concentrations lower than 40 g/l.

**Sensitivity/limit of detection**  
The limit of detection is 0.01 g/l.

**Interferences**  
No interference was observed by the presence of:  
hemoglobin ≤ 1000 mg/dl  
bilirubin ≤ 32 mg/dl  
lipids ≤ 1200 mg/dl  
rheumatoid factor ≤ 630 IU/ml

Precision	mean (g/l)	SD (g/l)	CV%
intra-assay (n=10)			
sample 1	2.20	0.02	1.15
sample 2	4.13	0.04	0.99

inter-assay (n=10)	mean (g/l)	SD (g/l)	CV%
sample 1	2.19	0.05	2.36
sample 2	4.12	0.08	2.05

**Methods comparison**  
A comparison between CHEMA and a commercially available product gave the following results:

$$\begin{aligned} \text{Transferrin competitor} &= x \\ \text{Transferrin FL CHEMA} &= y \\ n &= 82 \\ y &= 1.004x + 0.018 \text{ g/l} \quad r^2 = 0.991 \end{aligned}$$

**WASTE DISPOSAL**

This product is made to be used in professional laboratories.  
P501: Dispose of contents according to national/international regulations.

<b>TRANSFERRINE FL</b>	
TF 2H100	2 x 40 + 2 x 10 ml
TF 6U140	2 x 56 + 2 x 14 ml

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