APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN HITACHI 911/912	
TEST:	CRP
APP. CODE:	352
WAVELENGTH (Sec/Pri):	340
ASSAY:	2 POINT END <i>TIME: 10</i> <i>POINT: 16 - 31</i> <i>DILUENT: water</i>
SAMPLE VOL:	NORMAL: 11 DECREASE: 8 INCREASE: 15
	R1 VOLUME: 200 <i>DILUENT: 5</i> R2 VOLUME: 0 R3 VOLUME: 25 <i>DILUENT: 5</i> R4 VOLUME: 0
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LOGIT-LOG (5P) <i>(POINT: 6 - SPAN: 6 - WEIGHT: 0)</i>
SD LIMIT:	0.50
DUPLICATE LIMIT:	10%
ST. 1 CONC:	0.0
EXPECTED VALUE:	0.0 - 5.0
UNIT:	mg/l
INSTR. FACTOR (y=ax+b):	a=1 b= 0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN OLYMPUS AU 400/480/600/640/680/2700 (Test code 885)	
TEST NAME:	CRP
SAMPLE:	Volume 11 µl Dilution 0 µl
REAGENTS:	R1 Volume 200 µl Dilution 0 µl R2 Volume 25 µl Dilution 0 µl
WAVELENGTH:	Pri. 340 Sec.
METHOD:	END
REACTION SLOPE:	+
MEASURING POINT 1:	First 0 Last 27
MEASURING POINT 2:	First 0 Last 10
REAGENT OD LIMIT:	First L -0.1 First H 0.5 Last L -0.1 Last H 0.5
DYNAMIC RANGE:	L 1.0 H 200
CORRELATION FACTOR:	A 1 B 0
UNIT:	mg/l
CALIBRATION TYPE:	5AB
FORMULA:	POLYGONAL

 Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

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INTENDED USE
Reagent for quantitative in vitro determination of C-reactive protein in biological fluids.

PRINCIPLE OF THE METHOD
C-Reactive Protein (CRP) selectively reacts with an anti-CRP antibody and forms an immunocomplex. The produced turbidity is proportional to the concentration of CRP in the sample, and can be measured at the wavelength of 340 nm.

KIT COMPONENTS
For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label.
Keep away from direct light sources.

CRP R1 **2H135: 3 x 40 ml (liquid) white cap**
6U126: 2 x 56 ml (liquid) white cap

Composizione: Buffer pH 7.50, PEG ≥ 2%, stabilizers and preservatives.

CRP R2 **2H135: 1 x 15 ml (liquid) red cap**
6U126: 1 x 14 ml (liquid) red cap

Composizione: Anti-human CRP antibody ≥ 2%, stabilizers and preservatives.

Store all components at 2-8°C.

REAGENT PREPARATION
Use separate reagent ready to use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-.
Caution: keep well refrigerated.

PRECAUTIONS
CRP R1: Danger. Causes serious eye damage (H318).
Wear protective gloves. Eye protection (P280).
 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). Immediately call a doctor (P310).

CRP R2: It is not classified as hazardous.

SPECIMEN
Serum, plasma.
Keep specimens away from direct light sources.
Samples are stable 7 days when stored at 2-8°C and 1 month at -20°C.

EXPECTED VALUES
Adulti < 5 mg/l

Each laboratory should establish appropriate reference intervals related to its population

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. If required, a human based calibrator is available:
CRP CALIBRATOR

Please contact Customer Care for further information.

TEST PERFORMANCE
Measure interval
Measure interval depends on the concentration of the highest standard used for calibration.
If such a limit value is exceeded, it is suggested to dilute sample 1+4 with distilled water and to repeat the test, multiplying the result by 5.

Hook effect
No Hook effect is observed with concentrations lower than 2600 mg/l.

Sensitivity/limit of detection
The limit of detection is 1.0 mg/l.

Interferences
No interference was observed by the presence of:
hemoglobin ≤ 1000 mg/dl
bilirubin ≤ 30 mg/dl
lipids ≤ 1300 mg/dl
rheumatoid factor ≤ 415 IU/ml

Precision

intra-assay (n=10)	mean (mg/l)	SD (mg/l)	CV%
sample 1	28.6	0.44	1.54
sample 2	69.0	0.71	1.03

inter-assay (n=20)	mean (mg/l)	SD (mg/l)	CV%
sample 1	28.0	0.88	3.13
sample 2	69.2	1.87	2.70

Methods comparison
A comparison between CHEMA and a commercially available product gave the following results:

$$\begin{aligned} \text{CRP competitor} &= x \\ \text{CRP FL CHEMA} &= y \\ n &= 105 \\ y &= 0.96x + 0.07 \text{ mg/l} \quad r^2 = 0.998 \end{aligned}$$

WASTE DISPOSAL
This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.

CRP FL	
RP 2H135	3 x 40 + 1 x 15 ml
RP 6U126	2 x 56 + 1 x 14 ml

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