APPLICAZIONE / APPLICATION / APLICACIÓN HITACHI 911/912				
TEST:		TRF		
APP. CODE:		356		
WAVELENGTH (Sec/Pri):		505		
ASSAY:		2 POINT END	TIME: 10 POINT: 16 - 31 DILUENT: water	
SAMPLE VOL:		NORMAL: 2 DECREASE: 2 INCREASE: 5	R1 VOLUME: 240 R2 VOLUME: 0 R3 VOLUME: 60 R4 VOLUME: 0 DILUENT: 5	
ABS LIMIT:		32000 - INC		
PROZONE LIMIT:		0 - UPPER		
CALIB METHOD:		LOGIT-LOG (5P) (POINT: 6 - SPAN: 6 - WEIGHT: 0)		
SD LIMIT:		0.50		
DUPLICATE LIMIT:		10%		
ST. 1 CONC:		0.0		
EXPECTED VALUE:		1.0 - 4.0		
UNIT:		g/l		
INSTR. FACTOR (y=ax+b):		a=1	b= 0	

APPLICAZIONE / APPLICATION / APLICACIÓN OLYMPUS AU 400/480/600/640/680/2700 (Test code 889)			
TEST NAME:		TRF	
SAMPLE:	Volume 2 µl	Dilution 0 µl	
REAGENTS:	R1 Volume 240 µl	Dilution 0 µl	
	R2 Volume 60 µl	Dilution 0 µl	
WAVELENGTH:		Pri. 520 Sec.	
METHOD:		END	
REACTION SLOPE:		+	
MEASURING POINT 1:	First 0	Last 27	
MEASURING POINT 2:	First 0	Last 10	
REAGENT OD LIMIT:	First L -0.1	First H 0.5	
	Last L -0.1	Last H 0.5	
DYNAMIC RANGE:		L 0.01	H 9.0
CORRELATION FACTOR:	A 1	B 0	
UNIT:		g/l	
CALIBRATION TYPE:		5AB	
FORMULA:		POLYGONAL	

Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

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INTENDED USE
Reagent for quantitative in vitro determination of transferrin in biological fluids.

PRINCIPLE OF THE METHOD
Transferrin selectively reacts with an anti-transferrin antibody and forms an immunocomplex. The produced turbidity is proportional to the concentration of transferrin in the sample, and can be measured at the wavelenght of 505 nm.

KIT COMPONENTS	
For in vitro diagnostic use only. The components of the kit are stable until expiration date on the label. Keep away from direct light sources.	
TRF R1	2H100: 2 x 40 ml (liquid) white cap 6U140: 2 x 56 ml (liquid) white cap

Composition:	Buffer pH 7.00, PEG ≥ 2%, stabilizers and preservatives.
--------------	--

TRF R2	2H100: 2 x 10 ml (liquid) red cap 6U140: 2 x 14 ml (liquid) red cap
--------	--

Composition:	Anti-human transferrin antibody ≥ 2%, stabilizers and preservatives.
--------------	--

Store all components at 2-8°C.

REAGENT PREPARATION
Use separate reagent ready to use. Stability: up to expiration date on labels at 2-8°C. Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-. Caution: keep well refrigerated.

PRECAUTIONS
Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow. Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

SPECIMEN
Serum, plasma. Keep specimens away form direct light sources. Samples are stable 7 days when stored at 2-8°C and 1 month at -20°C.

EXPECTED VALUES	
Newborns	1.2-2.5 g/l
Adults	2.0-3.6 g/l

Each laboratory should establish appropriate reference intervals related to its population

QUALITY CONTROL AND CALIBRATION
It is suggested to perform an internal quality control. For this purpose the following human based control serum is available: QUANTINORM CHEMA with normal or close to normal control values. If required, a multiparametric, human based calibrator is available:
REFERENCE P MULTICALIBRATOR

Please contact Customer Care for further information.

TEST PERFORMANCE
Measure interval Measure interval depends on the concentration of the highest standard used for calibration. If such a limit value is exceeded, it is suggested to dilute sample 1+4 with distilled water and to repeat the test, multiplying the result by 5.

Hook effect No Hook effect is observed with concentrations lower than 40 g/l.
--

Sensitivity/limit of detection The limit of detection is 0.01 g/l.

Interferences No interference was observed by the presence of: hemoglobin ≤ 1000 mg/dl bilirubin ≤ 32 mg/dl lipids ≤ 1200 mg/dl rheumatoid factor ≤ 630 IU/ml
--

Precision intra-assay (n=10) sample 1 sample 2	mean (g/l) 2.20 4.13	SD (g/l) 0.02 0.04	CV% 1.15 0.99
---	----------------------------	--------------------------	---------------------

inter-assay (n=10) sample 1 sample 2	mean (g/l) 2.19 4.12	SD (g/l) 0.05 0.08	CV% 2.36 2.05
--	----------------------------	--------------------------	---------------------

Methods comparison A comparison between CHEMA and a commercially available product gave the following results:

Transferrin competitor = x Transferrin FL CHEMA = y n = 82	
y = 1.004x + 0.018 g/l	r ² = 0.991

WASTE DISPOSAL
This product is made to be used in professional laboratories. P501: Dispose of contents according to national/international regulations.

TRANSFERRINE FL	
TF 2H100	2 x 40 + 2 x 10 ml
TF 6U140	2 x 56 + 2 x 14 ml

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