APPLICAZIONE / APPLICATION / APLICACIÓN HITACHI 911/912			
TEST: CRP			
APP. CODE: 352			
WAVELENGTH (Sec/Pri): 340			
ASSAY: 2 POINT END		TIME: 10 POINT: 16 - 31 DILUENT: water	
SAMPLE VOL:	NORMAL: 11		
	DECREASE: 8		
	INCREASE: 15		
	R1 VOLUME: 200	DILUENT: 5	
	R2 VOLUME: 0		
	R3 VOLUME: 25	DILUENT: 5	
	R4 VOLUME: 0		
ABS LIMIT: 32000 - INC			
PROZONE LIMIT: 0 - UPPER			
CALIB METHOD: LOGIT-LOG (5P) (POINT: 6 - SPAN: 6 - WEIGHT: 0)			
SD LIMIT: 0.50			
DUPLICATE LIMIT: 10%			
ST. 1 CONC: 0.0			
EXPECTED VALUE: 0.0 - 5.0			
UNIT: mg/l			
INSTR. FACTOR (y=ax+b): a= 1 b= 0			

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN OLYMPUS AU 400/480/600/640/680/2700 (Test code 885)			
TEST NAME: CRP			
SAMPLE:	Volume 11 µl		Dilution 0 µl
REAGENTS:	R1 Volume 200 µl		Dilution 0 µl
	R2 Volume 25 µl		Dilution 0 µl
WAVELENGHT: Pri. 340 Sec.			
METHOD: END			
REACTION SLOPE: +			
MEASURING POINT 1:	First 0		Last 27
MEASURING POINT 2:	First 0		Last 10
REAGENT OD LIMIT:	First L -0.1	First H 0.5	
	Last L -0.1	Last H 0.5	
DYNAMIC RANGE: L 1.0 H 200			
CORRELATION FACTOR: A 1 B 0			
UNIT: mg/l			
CALIBRATION TYPE: 5AB			
FORMULA: POLYGONAL			

Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

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INTENDED USE

Reagent for quantitative in vitro determination of C-reactive protein in biological fluids.

PRINCIPLE OF THE METHOD

C-Reactive Protein (CRP) selectively reacts with an anti-CRP antibody and forms an immunocomplex. The produced turbidity is proportional to the concentration of CRP in the sample, and can be measured at the wavelength of 340 nm.

KIT COMPONENTS

For in vitro diagnostic use only.

The components of the kit are stable until expiration date on the label.

Keep away from direct light sources.

CRP R1 **2H135: 3 x 40 ml (liquid) white cap**
6U126: 2 x 56 ml (liquid) white cap

Composition: Buffer pH 7.50, PEG ≥ 2%, stabilizers and preservatives.

CRP R2 **2H135: 1 x 15 ml (liquid) red cap**
6U126: 1 x 14 ml (liquid) red cap

Composition: Anti-human CRP antibody ≥ 2%, stabilizers and preservatives.

Store all components at 2-8°C.

REAGENT PREPARATION

Use separate reagent ready to use.

Stability: up to expiration date on labels at 2-8°C.

Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-.

Caution: keep well refrigerated.

PRECAUTIONS

CRP R1: Danger. Causes serious eye damage (H318).
 Wear protective gloves. Eye protection (P280). IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). Immediately call a doctor (P310).

CRP R2: It is not classified as hazardous.

SPECIMEN

Serum, plasma.

Keep specimens away form direct light sources.

Samples are stable 7 days when stored at 2-8°C and 1 month at -20°C.

EXPECTED VALUES

Adults < 5 mg/l

Each laboratory should establish appropriate reference intervals related to its population

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control.
If required, a human based calibrator is available:

CRP CALIBRATOR

Please contact Customer Care for further information.

TEST PERFORMANCE

Measure interval

Measure interval depends on the concentration of the highest standard used for calibration.

If such a limit value is exceeded, it is suggested to dilute sample 1+4 with distilled water and to repeat the test, multiplying the result by 5.

Hook effect

No Hook effect is observed with concentrations lower than 2600 mg/l.

Sensitivity/limit of detection

The limit of detection is 1.0 mg/l.

Interferences

No interference was observed by the presence of:

hemoglobin	≤ 1000 mg/dl
bilirubin	≤ 30 mg/dl
lipids	≤ 1300 mg/dl
rheumatoid factor	≤ 415 IU/ml

Precision

intra-assay (n=10)	mean (mg/l)	SD (mg/l)	CV%
sample 1	28.6	0.44	1.54
sample 2	69.0	0.71	1.03

inter-assay (n=20)	mean (mg/l)	SD (mg/l)	CV%
sample 1	28.0	0.88	3.13
sample 2	69.2	1.87	2.70

Methods comparison

A comparison between CHEMA and a commercially available product gave the following results:

CRP competitor = x
CRP FL CHEMA = y
n = 105

y = 0.96 x + 0.07 mg/l r² = 0.998

WASTE DISPOSAL

This product is made to be used in professional laboratories.

P501: Dispose of contents according to national/international regulations.

CRP FL	
RP 2H135	3 x 40 + 1 x 15 ml
RP 6U126	2 x 56 + 1 x 14 ml

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