APPLICAZIONE / APPLICATION / APLICACIÓN HITACHI 911/912			
TEST: IGA			
APP. CODE: 353			
WAVELENGTH (Sec/Pri): 600			
ASSAY: 2 POINT END		TIME: 10 POINT: 16 - 31 DILUENT: water	
SAMPLE VOL:	NORMAL: 2		
	DECREASE: 2		
	INCREASE: 5		
	R1 VOLUME: 240		
	R2 VOLUME: 0		
	R3 VOLUME: 60		
	R4 VOLUME: 0		
ABS LIMIT: 32000 - INC		DILUENT: 5	
PROZONE LIMIT: 0 - UPPER			
CALIB METHOD:		LOGIT-LOG (5P) (POINT: 6 - SPAN: 6 - WEIGHT: 0)	
SD LIMIT: 0.50			
DUPLICATE LIMIT: 10%			
ST. 1 CONC: 0.0			
EXPECTED VALUE: 0.0 - 4.0			
UNIT: g/l			
INSTR. FACTOR (y=ax+b): a=1 b= 0			

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN OLYMPUS AU 400/480/600/640/680/2700 (Test code 886)			
TEST NAME: IGA			
SAMPLE:		Volume 2 µl	Dilution 0 µl
REAGENTS:		R1 Volume 240 µl R2 Volume 60 µl	Dilution 0 µl Dilution 0 µl
WAVELENGTH: Pri. 600 Sec.			
METHOD: END			
REACTION SLOPE: +			
MEASURING POINT 1:		First 0	Last 27
MEASURING POINT 2:		First 0	Last 10
REAGENT OD LIMIT:		First L -0.1 Last L -0.1	First H 0.5 Last H 0.5
DYNAMIC RANGE:		L 0.01	H 6.0
CORRELATION FACTOR:		A 1	B 0
UNIT: g/l			
CALIBRATION TYPE: 5AB			
FORMULA: POLYGONAL			

Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

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INTENDED USE
Reagent for quantitative in vitro determination of IgA in biological fluids.

PRINCIPLE OF THE METHOD
Immunoglobulins A (IgA) selectively react with an anti-IgA antibody and form an immunocomplex. The produced turbidity is proportional to the concentration of IgA in the sample, and can be measured at the wavelenght of 600 nm.

KIT COMPONENTS
For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label.
Keep away from direct light sources.

IGA R1 2H100 2 x 40 ml (liquid) white cap
6U140 2 x 56 ml (liquid) white cap

Composition: Buffer pH 7.50, PEG ≥ 2%, stabilizers and preservatives.

IGA R2 2H100 2 x 10 ml (liquid) red cap
6U140 2 x 14 ml

(liquid) red cap

Composition: Anti-human IgA antibody ≥ 2%, stabilizers and preservatives.

REAGENT PREPARATION
Use separate reagent ready to use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-.
Caution: keep well refrigerated.

PRECAUTIONS
Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.
Perform the test according to the general “Good Laboratory Practice” (GLP) guidelines.

SPECIMEN
Serum, plasma.
Keep specimens away form direct light sources.
Samples are stable 7 days when stored at 2-8°C and 1 month at -20°C.

EXPECTED VALUES	
Newborns	0-0.02 g/l
Adults	0.7-4.0 g/l

Each laboratory should establish appropriate reference intervals related to its population

QUALITY CONTROL AND CALIBRATION
It is suggested to perform an internal quality control. For this purpose the following human based control serum is available:

QUANTINORM CHEMA
with normal or close to normal control values.
If required, a multiparametric, human based calibrator is available:
REFERENCE P MULTICALBRATOR

Please contact Customer Care for further information.

TEST PERFORMANCE
Measure interval
Measure interval depends on the concentration of the highest standard used for calibration.
If such a limit value is exceeded, it is suggested to dilute sample 1+4 with distilled water and to repeat the test, multiplying the result by 5.

Hook effect
No Hook effect is observed with concentrations lower than 38 g/l.

Sensitivity/limit of detection
The limit of detection is 0.01 g/l.

Interferences
No interference was observed by the presence of:
hemoglobin ≤ 1000 mg/dl
bilirubin ≤ 27 mg/dl
lipids ≤ 1290 mg/dl
rheumatoid factor ≤ 630 IU/ml

Precision			
intra-assay (n=10)	mean (g/l)	SD (g/l)	CV%
sample 1	1.46	0.02	1.04
sample 2	2.38	0.02	0.94

inter-assay (n=20)	mean (g/l)	SD (g/l)	CV%
sample 1	1.46	0.07	4.99
sample 2	2.40	0.10	4.16

Methods comparison
A comparison between CHEMA and a commercially available product gave the following results:

IgA competitor = x IgA FL CHEMA = y n = 21	
y = 0.97x + 0.04 g/l	r² = 0.99

WASTE DISPOSAL
This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.

IgA FL		
GA 2H100	2 x 40 + 2 x 10 ml	
GA 6U140	2 x 56 + 2 x 14 ml	

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