FOSFATASI ALCALINA FL IFCC

AF 2H400	8 x 40 + 4 x 20 ml
AF 6U420	6 x 56 + 6 x 14 ml

USO

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emoglobina	≤ 400 mg/dl
bilirubina	≤ 40 mg/dl
lipidi	≤ 900 mg/dl

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INTENDED USE

Reagent for quantitative in vitro determination of alkaline phosphatase in biological fluids.

PRINCIPLE OF THE METHOD

The enzyme alkaline phosphatase (EC 3.1.3.1, ortho-phosphoric-monoester phosphohydrolase) hydrolyzes the 4-NPP to release 4-nitrophenol, under alkaline conditions. The 4-nitrophenol formed is detected spectrophotometrically at 405 nm to give a measurement of alkaline phosphatase activity in the sample.

The present method has been made according to IFCC.

KIT COMPONENTS**For in vitro diagnostic use only.**

The components of the kit are stable until expiration date on the label.

Keep away from direct light sources.

ALP IFCC R1 2H400: 8 x 40 ml (liquid) white cap
6U420: 6 x 56 ml (liquid) white cap

ALP IFCC R2 2H400: 4 x 20 ml (liquid) red cap
6U420: 6 x 14 ml (liquid) red cap

Composition in the test: 2-amino-2-methyl-1-propanol buffer 0.35 M pH 10.40 (30°C), magnesium acetate 2 mM, zinc sulfate 1 mM, HEDTA 2 mM, 4-NPP 16 mM.

Store all components at 2-8°C.

REAGENT PREPARATION

Use separate reagent ready to use.

Stability: up to expiration date on labels at 2-8°C.

Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-.

Caution: keep well refrigerated.

PRECAUTIONS

ALP IFCC R1: Danger. Causes serious eye damage (H318). Causes skin irritation (H315). IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. (P305+P351+P338). Wear protective gloves / eye protection / face protection (P280). Immediately call a POISON CENTER / doctor (P310). Wash with water thoroughly after handling (P264).

ALP IFCC R2: It is not classified as hazardous.

SPECIMEN

Serum, plasma (heparinate only).

Sera kept at room temperatures usually show a slight but real increase in activity, which varies from 1% over a 6-h period to 3 to 6% over a 1 to 4 days period. Even in sera stored at refrigerator temperature, activity increases slowly. In frozen sera, activity decreases but slowly recovers after thawing the serum.

A similar enhancement of activity, but of greater magnitude, occurs with reconstituted lyophilized preparations, such as those available as control sera or calibrators. In reconstituted material the increases with storage at 4 and 20°C are about 10 and 30%, respectively. Enhancement of activity continues for several days, but at a decreasing rate. The cause of this phenomenon is not known but may be attributed to renaturation of partially denatured enzyme or to dissociation, on warming, of a phosphate-lipoprotein complex or a multimer of the enzyme that was formed in the freeze-drying process.

EXPECTED VALUES

Men:	40 - 129 U/l	(0.67 - 2.15 µkat/l)
Women:	35 - 104 U/l	(0.58 - 1.74 µkat/l)

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

QUANTINORM CHEMA - MULTINORM CHEMA

with normal or close to normal control values

QUANTIPATH CHEMA - MULTIPATH CHEMA

with pathological control values.

If required, a multiparametric, human based calibrator is available:

AUTOCAL H

Please contact Customer Care for further informations.

TEST PERFORMANCE**Linearity**

the method is linear up to 3000 U/l.

If a ΔA/min of 0.500 is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)

the limit of detection is 5.2 U/l.

Interferences

no interference was observed by the presence of:

hemoglobin	≤ 400 mg/dl
bilirubin	≤ 40 mg/dl
lipids	≤ 900 mg/dl

Precision

intra-assay (n=10)	mean (U/l)	SD (U/l)	CV%
sample 1	84.40	2.41	2.86
sample 2	222.40	5.74	2.58

inter-assay (n=20)	mean (U/l)	SD (U/l)	CV%
sample 1	86.66	2.66	3.07
sample 2	210.39	6.08	2.89

Methods comparison

a comparison between Chema and a commercially available product gave the following results:

$$y = 1.03x - 2.57 \quad r^2 = 0.998$$

WASTE DISPOSAL

This product is made to be used in professional laboratories.

P501: Dispose of contents according to national/international regulations.

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ
HITACHI 911/912

TEST: **ALP IFCC**

APP. CODE: **343**

WAVELENGTH (Sec/Pri): **700 - 415**

ASSAY: **RATE-A** *TIME: 10*
POINT: 20 - 30
DILUENT: water

SAMPLE VOL: NORMAL: **7**
DECREASE: **5**
INCREASE: **10**

R1 VOLUME: **280** *DILUENT: 5*
R2 VOLUME: **0**
R3 VOLUME: **70** *DILUENT: 5*
R4 VOLUME: **0**

ABS LIMIT: **32000 - INC**

PROZONE LIMIT: **0 - UPPER**

CALIB METHOD: **LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)**

SD LIMIT: **0.250**

DUPLICATE LIMIT: **3%**

ST. 1 CONC: **0.0**

EXPECTED VALUE: **30 - 130**

UNIT: **U/l**

INSTR. FACTOR (y=ax+b): a=1 b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ
OLYMPUS AU 400/480/600/640/680/2700 (Test code 873)

TEST NAME: **ALP IFCC**

SAMPLE: Volume **7 µl** Dilution **0 µl**

REAGENTS: R1 Volume **280 µl** Dilution **0 µl**
R2 Volume **70 µl** Dilution **0 µl**

WAVELENGTH: Pri. **410** Sec. **700**

METHOD: **RATE**

REACTION SLOPE: **+**

MEASURING POINT 1: First **14** Last **27**

MEASURING POINT 2: First Last

REAGENT OD LIMIT: First L **-0.1** First H **1.0**
Last L **-0.1** Last H **1.0**

DYNAMIC RANGE: L **5** H **3000**

CORRELATION FACTOR: A **1** B **0**

LINEARITY LIMIT: **15%**

UNIT: **U/l**

CALIBRATION TYPE: **AB**

FORMULA: **Y = AX + B**

Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com