APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ <b>HITACHI 911/912</b>	
TEST:	<b>LDL</b>
APP. CODE:	<b>059</b>
WAVELENGTH (Sec/Pri):	<b>700 - 600</b>
ASSAY:	<b>2-POINT END</b> <i>TIME: 10</i> <i>POINT: 16 - 31</i>
SAMPLE VOL:	NORMAL: <b>3</b> DECREASE: <b>2</b> INCREASE: <b>5</b>
	R1 VOLUME: <b>270</b> R2 VOLUME: <b>0</b> R3 VOLUME: <b>90</b> R4 VOLUME: <b>0</b>
ABS LIMIT:	<b>32000 - INC</b>
PROZONE LIMIT:	<b>0 - UPPER</b>
CALIB METHOD:	<b>LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)</b>
SD LIMIT:	<b>0.250</b>
DUPLICATE LIMIT:	<b>3%</b>
ST. 1 CONC:	<b>0.0</b>
EXPECTED VALUE:	<b>70.0 - 220.0</b>
UNIT:	<b>mg/dl</b>
INSTR. FACTOR (y=ax+b):	a=1    b= 0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ <b>OLYMPUS AU 400/480/600/640/680/2700 (Test code 867)</b>	
TEST NAME:	<b>LDL</b>
SAMPLE:	Volume <b>3 µl</b> Dilution <b>0 µl</b>
REAGENTS:	R1 Volume <b>270 µl</b> Dilution <b>0 µl</b> R2 Volume <b>90 µl</b> Dilution <b>0 µl</b>
WAVELENGTH:	Pri. <b>600</b> Sec. <b>700</b>
METHOD:	<b>END</b>
REACTION SLOPE:	<b>+</b>
MEASURING POINT 1:	First <b>0</b> Last <b>27</b>
MEASURING POINT 2:	First <b>0</b> Last <b>10</b>
REAGENT OD LIMIT:	First L <b>-0.1</b> First H <b>0.5</b> Last L <b>-0.1</b> Last H <b>0.5</b>
DYNAMIC RANGE:	L <b>1.0</b> H <b>500</b>
CORRELATION FACTOR:	A <b>1</b> B <b>0</b>
UNIT:	<b>mg/dl</b>
CALIBRATION TYPE:	<b>AB</b>
FORMULA:	<b>Y = AX + B</b>

 Chema Diagnostica  
Via Campania 2/4  
60030 Monsano (AN) - ITALY - EU  
phone +39 0731 605064  
fax +39 0731 605672  
e-mail: mail@chema.com  
website: http://www.chema.com

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**INTENDED USE**  
Reagent for quantitative in vitro determination of LDL-cholesterol in biological fluids.

**PRINCIPLE OF THE METHOD**  
When a sample is mixed with Reagent R1, the protecting reagent binds to LDL and protects LDL from enzyme reactions. Cholesterol esterase (CHE) and cholesterol oxidase (CO) react with non-LDL lipoproteins [chylomicrons (CM), very low density lipoproteins (VLDL) and HDL]. Hydrogen peroxide produced by the enzyme reactions with non-LDL cholesterol is decomposed by catalase in Reagent R1. When Reagent R2 is added, the protecting reagent is removed from LDL and catalase is inactivated. In this second process, CHE and CO react only with LDL-C. Hydrogen peroxide produced by the enzyme reactions with LDL-C yields a color complex upon oxidase condensation with N-(3-sulfopropyl)-3-methoxy-5-methylaniline (HMMPs) and 4-aminoantipyrine (4AA) in the presence of peroxidase (POD). By measuring the absorbance of the blue color complex produced at approximately 600 nm, the LDL-C concentration in the sample can be calculated when compared vs. the absorbance of the LDL-C Calibrator.

**KIT COMPONENTS**  
**For in vitro diagnostic use only.**  
The components of the kit are stable until expiration date on the label. Keep away from direct light sources.

**LDL-C R1**    **2H080: 3 x 20 ml (liquid) white cap**  
**6U120: 2 x 45 ml (liquid) white cap**

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Store at 2-8°C. Do not freeze.

REAGENT PREPARATION

Use separate reagents ready to use.  
Stability: up to expiration date on labels at 2-8°C.  
Stability since first opening of vials: 30 days at 2-8°C.

**PRECAUTIONS**  
**LDL-C R1: Warning.** May cause an allergic skin reaction (H317). Wear protective gloves/protective clothing/eye protection/face protection (P280). IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower] (P303+P361+P353). If skin irritation or rash occurs: Get medical advice (P333+P313).

**LDL-C R2:** It is not classified as hazardous.

**SPECIMEN**  
Serum, plasma heparinate.  
Anticoagulants such as heparin, citrate, oxalate and EDTA do not have significant influences on the assay when they are used in their usual amounts.  
Samples with triglyceride concentrations exceeding 1000 mg/dl should be diluted and reanalyzed.

Use fresh specimens. Do not use specimens that repeated freeze-thaw, because lipoproteins may be denatured.

**EXPECTED VALUES**  
normal values: 76 - 218 mg/dl  
NCEP ATP's Decision cut-off points for LDL-C:  
desirable: < 130 mg/dl  
borderline high risk for CHD 130 - 159 mg/dl  
high risk for CHD ≥ 160 mg/dl

**QUALITY CONTROL AND CALIBRATION**  
It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

**QUANTINORM CHEMA**  
with normal or close to normal control values  
**QUANTIPATH CHEMA**  
with pathological control values.  
If required, a multiparametric, human based calibrator is available:  
**AUTOCAL H**  
Please contact Customer Care for further information.

**TEST PERFORMANCE**  
**Linearity**  
the method is linear up to 500 mg/dl.  
If the limit is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

**Sensitivity/limit of detection (LOD)**  
the limit of detection is 1 mg/dl.

**Interferences**  
no interference was observed by the presence of:  
hemoglobin ≤ 500 mg/dl  
bilirubin (free) ≤ 50 mg/dl  
bilirubin (conjugated) ≤ 40 mg/dl  
ascorbic acid ≤ 50 mg/dl  
lipid ≤ 2000 mg/dl

**Precision**  
intra assay (n=21)    mean (mg/dl)    SD (mg/dl)    CV%  
sample 1            85.6            1.46            1.71  
sample 2            129.6           2.28            1.76  
  
inter-assay (n=9)    mean (mg/dl)    SD (mg/dl)    CV%  
sample 3            87.8            1.69            1.92  
sample 4            129.6           1.99            1.53

**Methods comparison**  
a comparison between LDL-direct FL and a CDC Reference method (beta-quantification) gave the following results:

LDL-direct FL Chema = x  
Reference method = y  
n = 25

$y = 1.0015x - 0.715 \text{ mg/dl} \quad r^2 = 0.986$

**WASTE DISPOSAL**  
This product is made to be used in professional laboratories.  
P501: Dispose of contents according to national/international regulations.